DONJOY®

ULTRASLING IV ER

ENGLISH	2
ESPAÑOL	4
DEUTSCH	6
ITALIANO.....	8
FRANÇAIS.....	10
NEDERLANDS	12
PORTUGUÊS	14
ČEŠTINA	16
SLOVENČINA.....	18
DANSK.....	20
SVENSKA	22
SUOMI	24
NORSK.....	26
简体中文	28
日本語	30

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






DJO, LLC
5919 Sea Otter Place,
Suite 200
Carlsbad, CA 92010
USA

ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE:

The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INTENDED USE / INDICATIONS: For protocols involving immobilization in external rotation • Glenohumeral dislocation/subluxation • Capsular Shifts • Use after posterior shoulder stabilizations • Post-operative Bankart procedures • Post-operative positioning after release of severe anterior capsule contracture • Soft tissue strains/repairs

PERFORMANCE CHARACTERISTICS: Providing immobilization or controlled movement of the limb or body segment.

CONTRAINDICATIONS: None.

WARNINGS AND PRECAUTIONS:

- If you have skin disorders/injuries in the relevant part of the body, the product should be fitted and worn only after consultation with a medical professional.
- Do not use over open wounds.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
- Do not use this device if it was damaged and/or packaging has been opened.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INFORMATION: To ensure proper fit, product should be sized on patient pre-operatively. Product is shipped with strapping and sling configured to address post-op or rehab for a right shoulder/arm. The symbol for Right ® on cushion should be facing up. To convert for left shoulder/arm, detach sling from abductor cushion; flip cushion over so waist buckle is to front and the Left (L) symbol is facing up; attach sling to outside of cushion with Velcro strips.

1. Release waist belt on cushion, Velcro strips at wrist position on front of sling and forearm strap (Fig. A).
2. Align cushion on the injured side at waist/hip level with elbow fixed at 90°. Place cushion so line on top of cushion is parallel to the front of your body (if you stand at a table, the line should be parallel to the table).
3. Bring waist strap from cushion around back of waist and secure buckle at front of pillow (Fig. B). Adjust waist strap and trim as necessary for proper fit.
4. Place forearm into sling with arm as far back as possible in the sling (Fig. C). Secure Velcro straps at wrist position by first anchoring the strap closest to cushion, on cushion. Feed outside Velcro strap through inside D-ring and secure strap on outside of sling (Fig. D). Place small mesh thumb strap (i) between thumb and forefinger and secure on sling. Secure forearm strap (ii) to sling and cushion (Fig. E).
5. Using your free arm reach behind your body and slip arm through shoulder strap (Fig. F). This motion is similar to a putting on a backpack. One strap should be over the unaffected shoulder (clavicular area) and the other under the arm (straps B & C). Secure end of clip on Strap A to mate at front of cushion (Fig. G). Adjust and trim all straps for proper fit by removing Y-tab on end of straps, trim to proper length, replace Y-tab.

FOR O.R. APPLICATION: Complete application steps 1 – 4 followed by: Unfasten one end of axilla (underarm) strap. This strap is strap C for right arm or strap B for left arm. Bring shoulder strap (strap A) over and around neck area and bring axilla (underarm) strap under arm and feed through flat O-ring. Secure strap with Velcro tab.

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6. Adjust ALL straps for adequate stabilization of injured shoulder. Remove Y-tabs on strap ends and trim as needed. Replace Y-tabs and secure straps.
7. For exercise and rehab detach the outside Velcro strap at wrist area to increase the range of motion.

NOTE: Large, flat O-rings on front (Pectoral/Clavicle) and back (Scapula) should be positioned for maximum patient comfort.

STEP DOWN OPTION FOR REHAB: To use the UltraSling IV ER without the Abductor Cushion, remove waist strap from Cushion, thread waist strap through small D-rings located on front and back of sling. Secure strap around waist with Velcro tabs. Trim as needed.

CLEANING INSTRUCTIONS: Hand wash in cold water 30°C/85°F. Hang dry.

MATERIAL CONTENTS: Nylon, Polyester, Spandex, Polyurethane.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

RX ONLY.

FOR SINGLE PATIENT USE ONLY

**NOT MADE WITH NATURAL RUBBER
LATEX.**

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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