

enovis[™]

Creating Better Together[™]



CATALOGUE

2024/25

enovis[™]

Creating Better Together[™]

CATALOGUE

2024/25

WORKING WITH ENOVIS

Enovis® is a world leading orthopaedic company focused on helping people maintain and regain the joy of natural mobility. Offering a diverse portfolio of products and solutions to prevent injury, rehabilitate muscles, bones and joints and treat pain, we strive to help people live their lives to the fullest. We believe that our powerful brands, comprehensive range of products, extensive distribution network, unwavering focus on quality, innovation and customer service, and our strong relationships with orthopaedic and physical therapy professionals have all contributed to our leading market position. We have award winning manufacturing facilities across the globe with industry leading product design teams to continue the development and production of new products.

Our customer service team in Guildford is recognised in the industry as providing excellence in the support of our customers. Enovis Therapies and Appliances Managers throughout the UK and Ireland provide detailed product and sales support.

ENOVIS CUSTOMER CARE

In addition to fast and accurate delivery, you can expect a friendly and knowledgeable Customer Service Department to take your order and answer any questions. Each of our Therapies and Appliances Managers serves and supports an individual geographic territory, so you can count on expertise specific to your area. Should your area Account Manager be unavailable, any of our other representatives will be equally helpful and able to serve your needs.

You can be assured your order will ship the same day, providing the item(s) requested are available and order is placed before 2:30 PM. We look forward to serving you between the hours of 8.00^{AM}-4.00^{PM}.

HOW TO ORDER

To order from UK:

UK Freephone: 0800 587 0857

Landline: +44 (0)1483 459659

Fax: +44 (0)1483 459470

E-Mail: ukorders@enovis.com

QUESTIONS WE WILL ASK

Account Address: We will need to take the full delivery and invoice address.

Product: Please state the name, part number, description and size.

Purchase/sales order number: Please have this number if you call after receiving your order

SPECIAL REQUIREMENTS

The majority of orders will ship within 24 hours. Custom products and special orders (logos, replacement parts, special neoprene colours, etc.) are shipped within 6 weeks.

WARRANTY

All products have a limited warranty from the date of purchase against defects in materials and workmanship. Warranty varies with some products. Please contact Customer Services with specific questions.

RETURN POLICY

Before returning or exchanging any product, please call the Customer Service Department for a Return Authorisation Number. They will assist you with product specific return policies.

Please see page 134 for full terms and conditions.

OUR BRANDS

AIRCAST®

 chattanooga®

DONJOY®

PROCARE®



POPULAR PRODUCTS

[AIRSELECT STANDARD]page 53



ORDER NOW

[HEEL WEDGES]page 55



ORDER NOW

[XCELTRAX TALL]page 57



ORDER NOW

[ACTYFOOT]page 63



ORDER NOW

[PODALUX 2]page 45



ORDER NOW

[REDDIE BRACE]page 24



ORDER NOW

[CRYOCUFF]page 94



ORDER NOW

[ROAM OA]page 29



ORDER NOW

[CLINICAL SHOULDER IMMOBILISER]page 81



ORDER NOW

[X-ROM]page 39



ORDER NOW

ENOVIS ECO EFFORTS

We are highly conscious of the environmental impact of medical devices, and are delighted to share our Carbon Reduction Plan, which outlines Enovis UK's commitment to Net Zero emissions by 2050.

You can scan the QR code below to view. Here are several key points.

CURRENT CARBON REDUCTION PROJECTS:

- New Emissions Monitoring System to track Greenhouse Gas performance
- EV Charging points at UK Office
- Increased hybrid meetings to reduce emissions from travel
- Increased availability of 'bulk' packed products to reduce shipping emissions

TARGETS:

- Replace gas heating in our UK warehouse to electric alternatives by 2027
- All UK company vehicles to be hybrid by 2027 and EV from 2034
- Reduce electricity emissions at our Guildford HQ by 12% per year
- Reduce emissions from flights by 5% per year

HOW YOU CAN HELP:

There are also several ways in which hospitals can help us reduce our impact on the environment. See below our recommendations:

- Consolidate orders (E.g. order 20 Boots in one go rather than 2 or 3 a time)
- Order bulk packs of products to reduce multiple shipments, emissions and packaging
- Only provide patients with specific products when totally necessary
- Identify charities and organisations who may be able to help re-use equipment





ENOVIS EDUCATION & RESOURCES

At Enovis we are passionate about Creating Better Together.

That is better patient outcomes, better innovations for healthcare providers, and better resources to help clinicians offer the best service possible.

See below for a range of resources:

CUSTOMER SERVICES

- We have an experienced customer services team who are here to assist with product queries and orders. Should you require any assistance please contact via the below:
- Phone: **01483 459 659**
- Email: ukorders@enovis.com (Click the email to contact us)

LOCAL SALES REPRESENTATIVES

- Your hospital or clinic will have a dedicated local representative on hand to provide product training, facilitate trials, demonstrations, and ensure a smooth transition when introducing a new product. Please contact us for your local representative details.

YOUTUBE

- We have a dedicated YouTube channel with product fitting videos for both clinicians and patients which you can find [BELOW](#)

ENOVIS BLOG

- Filled with useful insight backed by clinical evidence, our blog is for clinicians looking to widen their knowledge of medical devices and therapies. This includes research on knee and ankle bracing, and modalities such as muscle stimulation and shockwave therapy. Click [BELOW](#) to visit

ENOVIS INSTITUTE

- A fantastic free resource for clinicians, complete with E-Learning, Research Library, Webinars and more. [SIGN UP HERE](#)

FITTING GUIDES AND BROCHURES

- For our key products, we have PDF patient fitting guides, as well as printed leaflets for patients (for a select number of products). Please contact your local representative for more information or email ukmarketing@enovis.com


YOUTUBE

[Click Here](#)



ENOVIS BLOG

[Click Here](#)






ENOVIS INSTITUTE

[Click Here](#)



THE BEST FIT

We at Enovis have categorised our diverse portfolio into three areas, **SMART**, **STANDARD** and **VALUE**, to help you to identify the product best suited to your requirements.



THE CATEGORIES:

SMART Effective products that accelerate recovery with high compliance and simplified ease-of-use

STANDARD Everyday, essential products at a competitive price

VALUE Simple solutions for optimising the cost-benefit ratio under tight economic constraints

	SMART	STANDARD	VALUE
ACCELERATED RECOVERY	○	○	○
INNOVATIVE SOLUTIONS	○	○	○
PATIENT-CENTRIC DESIGN	○	○	
PATIENTS WITH ACTIVE LIFESTYLES	○	○	
CLINICALLY PROVEN	○	○	
RANGE OF INDICATIONS	○	○	○
TIME-PROVEN PRODUCTS	○	○	○
BULK PRICING	○	○	○
ECONOMY		○	○

ONLINE PURCHASE

Private clinics and patients can purchase many Enovis products from the UK and Irish online shop, Motioncare. Visit www.motioncare.co.uk or www.motioncareireland.com. Clinics registering for an online account may also be entitled to discounts off the listed prices.



CONTENTS

[THE BEST FIT]

THE BEST FIT.....	9
ONLINE PURCHASE.....	9

[KNEE & LOWER LIMB]

4-POINTS-OF-LEVERAGE™ BRACING SYSTEM.....	16
HOW DONJOY 4-POINTS-OF- LEVERAGE PROTECTS ACL.....	17
FOURCEPOINT HINGE TECHNOLOGY.....	17
KNEE BRACE LADDER.....	18
LIGAMENT DEFICIENCIES.....	19
MEASURING GUIDE.....	19
DEFIANCE® PRO ●.....	20
DEFIANCE COLOUR OPTIONS.....	21
DEFIANCE HINGE OPTIONS.....	21
RENEGADE™ ●.....	22
4-TITUDE® ●.....	22
FULLFORCE™ ●.....	23
ARMOR™ FOURCEPOINT ●.....	23
ARMOR™ ●.....	24
DELUXE HINGED KNEE ●.....	24
REDDIE® BRACE ●.....	24
PLAYMAKER® II ●.....	25
PLAYMAKER XPERT ●.....	25
RIGID KNEE LIGAMENT BRACES ACCESSORIES.....	26
OSTEOARTHRITIS (OA).....	27
ROAM OA ●.....	29
OA DEFIANCE ●.....	30
OA ADJUSTER™ 3 ●.....	30
OA NANO™ ●.....	31
CLIMA-FLEX™ OA ●.....	32
OA REACTION™ WEB ●.....	32
STRAPILAX KNEE ●.....	33
ROTULAX ●.....	33
FORTILAX KNEE ●.....	34
INFRAPATELLAR BAND ●.....	34
REACTION™ WEB ●.....	35
PATELAX ●.....	35
TRU-PULL® LITE ●.....	36
STABILAX ●.....	36
GENULAX ●.....	37
X-ROM™ ●.....	39
X-ACT™ ROM KNEE ●.....	40
IMMO KNEE 3V (3-PANEL KNEE IMMOBILISER) ●.....	40
IMMO AT4 3V KNEE IMMOBILISER ●.....	41
LEG ELEVATOR ●.....	41
VERSAROM HIP™ ●.....	41

[FOOT]

ACTYTOE ●	44
REMEDY PRO™ ●	45
SQUARED TOE POST-OP SHOE ●	45
PODALUX II ●	45
ROCKER CAST BOOT ●	46
MED/SURG SHOE ●	46
POST-OP SHOE ●	47
OFF-LOADING DIABETIC SHOE ●	47
AIRLIFT™ PTTD BRACE ●	48
PODOUS BOOT ●	49
PLANTAR FASCIITIS NIGHT SPLINT ●	49
PROWEDGE™ NIGHT SPLINT ●	49
PODALIB® ●	49

[WALKER BOOTS]

AIRSELECT™ ELITE ●	52
AIRSELECT™ STANDARD ●	53
AIRSELECT™ SHORT ●	53
XP DIABETIC WALKER™ ●	54
XP DIABETIC WALKER REPLACEMENT PARTS ●	54
TUBE STRETCH SOCK ●	55
HYGIENE COVER	55
WEATHER COVER ●	55
10° WEDGE KIT ●	55
MAXTRAX XCEL ROM ●	56
NEXTEP™ CONTOUR WALKER ●	56
NEXTEP™ CONTOUR SHORTIE WALKER ●	56
XCELTRAX™ TALL ●	57
XCELTRAX™ ANKLE ●	57
MINITRAX™ ●	58
SHOELIFT™ ●	58

[WALKER BOOTS]

ACTYFOOT™ ●	63
A60™ ANKLE SUPPORT ●	64
AIRHEEL™ ●	64
AIRSPORT™ ●	64
AIRSPORT+® ●	64
AIR-STIRRUP® UNIVERSE™ ●	65
AIR-STIRRUP® ●	65
FUSIOGEL ●	65
FUSIOLIGHT II ●	66
STRAPILAX ANKLE ●	66
MALOLAX™ ●	67
FORTILAX™ ANKLE ●	67

[SPINE]

XTEND 174 CERVICAL COLLAR ●	70
TRANSITIONAL 172 CERVICAL COLLAR ●	70
ADAMS COLLAR ●	71
PREMIUM PANEL ELASTIC BINDER ●	71
ANATOMICAL RANGE OF SPINAL SUPPORT	72
GENERAL INDICATIONS & RANGE OF SPINAL SUPPORT	73
ISOFORM™ SIO ●	74
ISOFORM™ LO ●	75
ISOFORM™ LO+ ●	75
ISOFORM™ LSO ●	75
ISOFORM™ LSO+ ●	75
ISOFORM TLSO ●	76
ISOFORM POSTURAL EXTENSION ●	76
ACTISTRAP ●	77
SPINOSTRAP ●	77
IMMOSTRAP™ ●	77
MATERNITY BELT ●	77

BEST FIT KEY

See page 9 for BEST FIT guide

- SMART
- STANDARD
- VALUE

CONTENTS

[SHOULDER & ARM]

SULLY® SHOULDER STABILISER ●	80
HUMERAL FRACTURE BRACE/ OVER THE SHOULDER ●	80
HUMERAL STABILISING SYSTEM (HSS ●)	80
SHOULDER IMMOBILISER WITH FOAM STRAPS ●	81
CLINICAL SHOULDER IMMOBILISER ●	81
ULTRASLING QUADRANT™ ●	82
ULTRASLING PRO™ ●	82
CONDILAX™ ●	83
X-ACT ROM ELBOW ●	84
MAYO CLINIC ELBOW BRACE ●	84
CHIEFTAIN ARM SLING ●	85
BLUE VOGUE ARM SLING ●	85
HEALTHCARE BEAR ARM SLING ●	85
PNEUMATIC ARMBAND ●	85

[WRIST & HAND]

MANULAX ●	88
ERGOFORM™ ●	88
COMFORTFORM™ WRIST ●	89
COMFORTFORM™ WRIST WITH ABDUCTED THUMB ●	89
UNIVERSAL WRIST-O-PRENE ●	90
UNIVERSAL THUMB-O-PRENE™ ●	90
ELASTIC WRIST BRACE ●	90
WRIST WRAPS ●	91
CURVED ALUMINIUM FINGER SPLINTS ●	91
STAXX PLASTIC FINGER SPLINTS ●	91

[HOT/COLD THERAPY]

DURA*SOFT™ SHOULDER WRAP ●	93
DURA*SOFT™ KNEE WRAP ●	93
CRYO/CUFF™ IC COOLER ●	94
CRYO/CUFF™ GRAVITY COOLER ●	94
CRYO/CUFF™ WRAPS & CUFFS ●	95
ICEMAN® CLEAR3™ ●	96
ICEMAN® MCGUIRE KNEE WRAP ●	96
ICEMAN® FOREFOOT WRAP ●	96
ICEMAN® HIP WRAP-ON PAD ●	97
ICEMAN® WRAP-ON COLD PAD ●	97
ICEMAN® ANKLE CUFF ●	97
ICEMAN® SHOULDER WRAP-ON PAD ●	97
HYDROCOLLATOR ●	98
HOTPAC TERRY COVERS	99
DIAL THERMOMETER	99
HYDROCOLLATOR HOTPACS	100
COLPAC®	101
FLEXIPAC	101
COLPAC® CHILLING UNIT	101

[CMF]

CMF SPINALOGIC® ●	104
CMF OL1000™ DUAL COIL ●	104
CMF OL1000™ SINGLE COIL ●	105
CMF OL1000™ SINGLE COIL ●	105

[VASCULAR]

VENAFLOW® ELITE	107
VENAPURE®	108

[CLINICAL PRODUCTS]

INTELECT® FOCUS SHOCKWAVE 112

INTELECT® RPW 2 114

MOBILE 2 RPW SHOCKWAVE 116

RPW SHOCKWAVE ACCESSORIES 116

INTELECT® RPW LITE 117

RPW TRANSMITTERS 118

RUBBER ELECTRODES 119

ELECTRODE SPONGES 119

NYLATEX INDIVIDUAL 119

INTELECT® MOBILE 2 120

LIGHTFORCE® THERAPY LASERS 122

XLI/XPI/FXI 123

LIGHTFORCE® TREATMENT HEADS 124

VITALSTIM® PLUS 126

TENS/NMES 127

WIRELESS PRO 4 CHANNEL (FULL) 128

CEFAR TENS 130

REHAB X2 130

DURA-STICK 130

DURA-STICK® PLUS 131

DURA-STICK® PREMIUM 132

LUMBAR HOME TRACTION DEVICE 132

CERVICAL HOME TRACTION DEVICE 132

ARTROMOT®-K1 133

[MISCELLANEOUS]

PRODUCT WARRANTY 133

TERMS & CONDITIONS 134

BEST FIT KEY

See **page 9** for BEST FIT guide

- SMART
- STANDARD
- VALUE



DEFIANCE PRO

KNEE & LOWER LIMB

KNEE BRACING AT THE HEART OF ENOVIS HISTORY

While knee braces are a common sight, that hasn't always been the case. The modern knee brace, used to prevent injury, has only been around since 1967.

The late 60s were important for the evolution of knee braces. In 1967, a professor of physiology, Dr. Robert F. McDavid, invented the first lateral knee brace designed to prevent injury or reinjury. His brace provided lateral protection of the knee. Then during the 1967 season and Super Bowl III, Hall of Fame New York Jets' Quarterback Joe Namath played with a now-famous knee brace. Jack Castiglia of the Lenox Hill Brace Shop with along with noted sports physician Dr. James Nicholas designed the brace. This allowed Namath to keep playing despite being plagued by knee problems. Namath's knee brace was so iconic that it is now on display at the Pro Football Hall of Fame.

During the 1970s, more prophylactic knee braces began being tested and tried by more NFL players in an effort to reduce career ending knee injuries. In 1978, Philadelphia Eagles' offensive line captain Mark Nordquist and lawyer Ken Reed founded a small company in a Carlsbad, California garage. They named the new company DonJoy® after their wives, Donna and Joy.

The first Enovis products were simple sleeves made of sewn-together neoprene that were pulled over the knee as well as the ankle and elbow joints for support. Fast forward to today where knee braces come in a variety of shapes and sizes, designed for all manner of knee injuries and prevention.

COLD THERAPY KNEE OPTIONS

For more details see pages 92-101



DURA*SOFT™ KNEE WRAP



CRYO/CUFF™ KNEE



CRYO/CUFF™



ICEMAN™ MCGUIRE KNEE WRAP

4-POINTS-OF-LEVERAGE™ BRACING SYSTEM

The 4-Points-of-Leverage is a dynamic bracing system made up of a rigid cuff and strap configuration, whereby a posterior force is applied to the tibia, preventing anterior movement. This reduces the strain on the anterior cruciate ligament (ACL).

THE 4-POINTS-OF-LEVERAGE BRACE DESIGN:

- Significantly reduces ACL strain*
- Provides a constant dynamic load



POINT OF LEVERAGE #1

- The tibia is anchored by a cuff placed at the bottom of the posterior calf

POINT OF LEVERAGE #2

- The femur is anchored by a cuff placed at the top of the anterior thigh

POINT OF LEVERAGE #3


- A strap across the back of the lower thigh pushes the femur anteriorly

POINT OF LEVERAGE #4

- A strap placed on the anterior tibia applies an active constant load to prevent anterior tibia translation

*Fleming et al. (2000), The influence of functional Knee Bracing on the ACL Strain Biomechanics in Weight bearing and Nonweightbearing Knees, The American Journal of Sports Medicine, Vol. 6, p. 815-824

HOW DONJOY 4-POINTS-OF LEVERAGE PROTECTS ACL



Patients with an ACL deficient knee have an anterolateral instability that causes the knee to give way.

It is very difficult to control instability in an ACL deficient knee in a direct way using a knee brace.

Indirect stabilisation of the knee is achieved with the biomechanical 4 point design of the DonJoy brace.



The 4-point system is clinically proven to be the most effective way to prevent anterior translation of the tibia.

A DonJoy brace will provide a high level of stability for the ACL deficient knee by effectively reducing or preventing the anterior tibial translation through the use of the 4-points-of-leverage.

FOURCEPOINT HINGE TECHNOLOGY

Aggressive eccentric contractions of the quadriceps produce significant anterior tibial translation and can injure the ACL or permanently stretch out an ACL graft (DeMorat et al., AJSM, 2004). Additional studies show quadricep-intensive activities should be avoided at joint angles, particularly the last 30 degrees of extension, when strains are at a maximum. Enovis has developed a hinge to address maximising ACL and ACL graft protection in these "at risk" positions.

Small knee flexion angles at initial contact during landing tasks may contribute to ACL injuries. The FourcePoint Hinge can be used to control joint range of motion by increasing the knee flexion angle during landing, resulting in a decrease of anterior shear force on the ACL by 9-13% (Garret and Yu, AJSM, 2004).



The DonJoy FourcePoint utilises a spring mechanism in the hinge to apply gradually increasing resistance during knee extension. The resistance engages in the last 25 degrees of extension (relative to the extension stop).

The hinge has three levels of resistance with the option to turn off resistance if appropriate.

THE RESISTANCE SERVES THREE CRITICAL ROLES:

- First, it reduces the time spent near full extension or in the "at risk" position
- Second, it increases the posterior load applied by the brace on the tibia to prevent anterior tibial translation
- Third, it eliminates the extension shock felt when a patient extends into a 10 degree standard rigid stop


FourcePoint Hinge Technology works to enhance DonJoy's 4-Points-of-Leverage design by damping knee joint extension, which improves the mechanical performance of the brace and reduces shear forces at the knee. Addressing the "at-risk" position, the hinge resistance engages in the last 25 degrees of knee extension.

The FourcePoint hinge combined with the 4-Points-of-Leverage cuff and strapping design provides a more comfortable brace that reduces anterior shear forces at the knee. This is beneficial for patients as it reduces strain on the healing ACL graft, providing stability to an unstable, ACL deficient knee or for prophylactic use.


KNEE BRACE LADDER

Understand the product choices you have available to your knee injury patients from DonJoy.
 These braces can be purchased directly from Enovis via www.motioncare.co.uk*

ACL MCL LCL BRACING LADDER




CI BRACING LADDER




*Playmaker Xpert for MCL and LCL instability only

PCL BRACING LADDER



OA LADDER



LIGAMENT DEFICIENCIES

BRACE	LIGAMENT DEFICIENCIES SUPPORTED					BRACE OPTIONS/CHARACTERISTICS			
	ACL	MCL	PCL	LCL	COMBINED INSTABILITY (ACL & PCL)	HINGE	TYPE	STRAPPING OPTIONS	CALF CUFF OPTIONS
DEFIANCE PRO	○	○	○	○	○	Polycentric FourcePoint Drop-Lock OA	Rigid, customised	ACL, PCL or CI	Standard/Short
FULLFORCE	○					FourcePoint	Rigid	ACL	Standard/Short
RENEGADE	○	○		○	○	Polycentric	Rigid	ACL or CI	Standard/Short
ARMOR	○	○	○	○	○	Polycentric	Rigid	ACL, PCL or CI	Short
ARMOR FOURCEPOINT	○	○		○		FourcePoint	Rigid	ACL	Short
4-TITUDE			○			Polycentric	Rigid	PCL	Standard/Short
PLAYMAKER II	○	○		○		Polycentric	Soft, wraparound or sleeve	N/A	N/A
PLAYMAKER XPERT		○		○		Polycentric	Soft, wraparound	N/A	N/A
DELUXE HINGED KNEE		○		○		Polycentric	Soft, wraparound	N/A	N/A
REDDIE BRACE		○		○		Polycentric	Soft, wraparound	N/A	N/A

MEASURING GUIDE



CLICK HERE

WATCH THE MEASURING VIDEO FOR ASSISTANCE

Measure at three points on the leg and use the chart to determine the size required.

Thigh 15cm above mid-patella (kneecap) · Knee, mid-patella (kneecap) · Calf 15cm below mid-patella (kneecap)

If one of the three measurements is a size up, select that size

For example, Thigh Size 2, Knee Size 3, Calf Size 2. Order Size 3

Note that for OA braces, the off-loading hinge is located on the opposite side of the affected compartment.

SIZE	THIGH CIRC. 15cm ABOVE MID-PATELLA	KNEE CENTRE CIRC. MID-PATELLA	CALF CIRC. 15cm BELOW MID-PATELLA
1 (XS)	33 – 39 cm	30.5 – 33cm	25.5 – 30.5cm
2 (S)	39 – 47 cm	33 – 35.5cm	30.5 – 35.5cm
3 (M)	47 – 53 cm	35.5cm – 38cm	35.5cm – 40.5cm
4 (L)	53 – 60 cm	38 – 43cm	40.5 – 47cm
5 (XL)	60 – 67 cm	43 – 48.25cm	47 – 50.75cm
6 (XXL)	67 – 75 cm	48.25cm – 53.25cm	50.75cm – 56cm
7 (XXXL)	75 – 81 cm	53.25 – 58.5cm	56 – 61cm



DONJOY

DEFIANCE® PRO

DonJoy's Defiance® has been the world's leading custom knee brace for over 20 years. Over 2 million Defiance braces have been prescribed to elite professional athletes, middle-school soccer players, and athletes of all shapes, sizes, ages and genders.

Now, Defiance PRO pushes our dominance even farther, boasting the greatest strength-to-weight ratio of any knee brace on the market.

Professional athletes choose Defiance PRO for the ultimate confidence and freedom to move.

INDICATIONS

Moderate to severe ACL, PCL, CI, MCL and LCL instabilities. The OA Defiance is suitable for the active patient with osteoarthritis.



CHARACTERISTICS

- New low-profile design: The trusted high-performance carbon fiber frame now lighter and 5mm lower profile than Defiance Classic
- New C-6 liner material: Soft to the touch, moisture wicking and anti-microbial
- Less than 600g: Custom carbon fiber frame is built tough
- New silicone condyle pads: Comfortably help keep the brace in place
- Memory foam condyle pads: Contour to the knee (not shown)
- Swiveling internally mounted straps: Intimate fit around the anatomy for a lighter, more comfortable feel and improved anti-migration
- New soft touch strap ends: Flexible and strong hold keeps straps secure

ORDERING INFORMATION

SIZE	THIGH CIRC.	CALF CIRC.
Minimum	31.8 cm (12 1/2")	26 cm (10 1/4")
Maximum	97 cm (38 1/8")	66.5 cm (26 1/4")

Each Defiance is custom made. Measurements taken with CCMI Mark III

DEFIANCE PRO COLOUR OPTIONS



DEFIANCE HINGE OPTIONS



FOURCEPOINT ACL:

This patented hinge design helps protect the ACL by incorporating a dampening mechanism in the last 25° of extension to help keep the knee out of the at risk position and increase the posterior tibial load.

FOURCEPOINT PLUS:

The Plus hinge provides ACL protection benefits of the FourcePoint hinge with a reinforced steel outer plate ideally suited for the rigors of high-impact or collision sports

OA ADJUSTABLE HINGE:

The custom OA frame offers 1 to 7 degrees of off-load. The telescoping hinge helps provide additional off-loading up to 3 degrees. The memory foam condyle pads are included for increased patient comfort and compliance.

STANDARD LOW-PROFILE HINGE:

The polycentric hinge stabilizes the knee while matching natural knee kinematics. Polycentric hinges also maximize the anterior/posterior and varus/valgus stabilizing forces of the brace. This can be used for ACL, PCL, combined instability, MCL, or LCL indications.

DONJOY

RENEGADE™



INDICATIONS

Recommended for moderate to severe ACL, CI/MCL & LCL instabilities, offering comfortable support for low to mid-range sporting activities.

CHARACTERISTICS

- Ultra soft suede liners, moulded strap tabs and tapered tibia padding promotes improved patient comfort and fit
- Available in standard and short calf lengths
- Available in 7 sizes, XS-XXXL
- Lightweight, low profile, aircraft quality aluminium frame with ergonomic swooping thigh and calf cuff
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

PART NUMBER		DESCRIPTION
LEFT	RIGHT	
11-0890-x	11-0891-x	ACL, Standard Calf
11-0906-x	11-0905-x	CI, Standard Calf

Sizing information is on page 19

DONJOY

4-TITUDE®



INDICATIONS

PCL Instabilities.

For active, sedentary and adolescent patients. Ideally suited for non-contact sports and moderate activity, e.g. walking, golf.

CHARACTERISTICS

- 4 Points of Leverage system
- Available in 7 sizes, XS-XXXL
- Lightweight (500 gram) frame with mouldable cuffs for custom fit
- Low profile hinges also offer flexion-extension control
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

PART NUMBER		DESCRIPTION
LEFT	RIGHT	
11-0833-x-06060	11-0832-x-06060	PCL, Standard Calf

Sizing information is on page 19

DONJOY

FULLFORCE™



INDICATIONS

Moderate to severe ACL, MCL & LCL instabilities.
For active lifestyles involving mid-range sports e.g. tennis, netball, cricket.

CHARACTERISTICS

- Enhanced 4-Points-of-Leverage with ACL strapping
- Available in 7 sizes, XS-XXXL
- Low profile FourcePoint hinge
- Strong, mouldable aircraft quality aluminium frame with carbon fibre look finish
- Internally mounted swivelling strap tabs to help accommodate leg movement
- Swooping thigh cuff for inner thigh clearance & bilateral use
- Optional short or standard calf length
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

PART NUMBER		DESCRIPTION
LEFT	RIGHT	
11-3221-x	11-3220-x	ACL, Short Calf
11-0259-x	11-0258-x	ACL, Standard Calf

13-8925 FourcePoint Allen Key

Sizing information is on page 19

DONJOY

ARMOR™ FOURCEPOINT



INDICATIONS

Moderate to severe ACL, MCL and LCL instabilities.
For active lifestyles involving extreme sports e.g. skiing, motocross, contact sports.

CHARACTERISTICS

- Enhanced 4-Points-of-Leverage with ACL strapping
- Available in 7 sizes, XS-XXXL
- Low profile FourcePoint hinge with steel outer plate with optional flexion and extension control
- Mouldable, medial swooping frame for custom fit
- Short calf for boot clearance
- Strong aircraft quality aluminium frame with internally mounted strap tabs
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

PART NUMBER		THIGH CIRCUMFERENCE	DESCRIPTION
LEFT	RIGHT		
11-1441-1	11-1440-1	33cm-39cm	XS
11-1441-2	11-1440-2	39cm-47cm	S
11-1441-3	11-1440-3	47cm-53cm	M
11-1441-4	11-1440-4	53cm-60cm	L
11-1441-5	11-1440-5	60cm-67cm	XL
11-1441-6	11-1440-6	67cm-75cm	XXL
11-1441-7	11-1440-7	75cm-81cm	XXXL

13-8925 FourcePoint Allen Key

Sizing information is on page 19



INDICATIONS

Moderate to severe ACL, PCL, CI, MCL and LCL instabilities. For active lifestyles involving extreme sports e.g. skiing.

CHARACTERISTICS

- 4-Points-of-Leverage system
- Strapping configurations available for ACL, PCL or CI instabilities
- Available in 7 sizes, XS-XXXL
- Mouldable, medial swooping frame for custom fit and short calf for boot clearance
- Low profile hinges also offer flexion-extension control
- Strong aircraft quality aluminium frame to reduce brace weight and improve comfort
- Internally mounted strap tabs
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

PART NUMBER		DESCRIPTION
LEFT	RIGHT	
11-1005-x-26000	11-1004-x-26000	ACL, Short Calf
11-1026-x-26000	11-1025-x-26000	PCL, Short Calf
11-1030-x-26000	11-1029-x-26000	CI, Short Calf

Sizing information is on page 19



DONJOY
DELUXE
HINGED KNEE



INDICATIONS

Mild to moderate medial lateral control. Ideal for patients who have allergies or sensitivity to neoprene.

CHARACTERISTICS

- Lightweight neoprene alternative fabric that allows for improved airflow
- Polycentric hinges allow for extension control 10°, 20°, 30° and 40° (10° installed) and flexion control 45°, 60°, 75° and 90°
- Fits left or right knee
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	THIGH CIRC.	SIZE
81-07582	33cm-39cm	XS
81-07583	39cm-47cm	S
81-07585	47cm-53cm	M
81-07587	53cm-60cm	L
81-07588	60cm-67cm	XL
81-07589	67cm-75cm	XXL
81-07589-10	75cm-81cm	XXXL

Measurement taken 15cm above mid-patella



PROCARE
REDDIE® BRACE



INDICATIONS

Ideal for mild MCL and/or LCL sprains, mild arthritis, post arthroscopic meniscectomy, or mild global instabilities.

CHARACTERISTICS

- 3mm neoprene wraparound configuration suitable for the older patient
- Removable dual axis polycentric hinges provide medial/lateral support
- Popliteal cut-out for extra comfort in flexion
- Fits left or right knee
- Available in 8 sizes XS-XXXXL
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	THIGH CIRC.	SIZE
79-82392	34cm-39cm	XS
79-82393	39cm-45cm	S
79-82395	45cm-52cm	M
79-82397	52cm-58cm	L
79-82398	58cm-64cm	XL
79-82399	64cm-71cm	XXL
79-82399-10	71cm-77cm	XXXL
79-82399-11	77cm-83cm	XXXXL

Measurement taken 15cm above mid-patella

PLAYMAKER® II



Playmaker II
Pull-On Sleeve

INDICATIONS

Suitable for mild ACL, PCL, MCL and LCL instabilities. Suitable for moderate activity such as walking and golf. For use in warmer environments.

CHARACTERISTICS

- Advanced high tenacity nylon/spandex material promotes a cooling effect while delivering a lighter and slimmer profile
- New moulded strap ends are designed for easy application and durability
- Enhanced low profile hinge allows for a range of motion control
- Popliteal cut-out for extra comfort in flexion
- Fits left or right knee
- Available in 7 sizes XS-XXXL
- 6 month warranty

ORDERING INFORMATION

PART NUMBER		THIGH CIRCUMFERENCE	SIZE
PULL-ON SLEEVE	WRAPAROUND		
11-3495-1	11-3497-1	33cm-39cm	XS
11-3495-2	11-3497-2	39cm-47cm	S
11-3495-3	11-3497-3	47cm-53cm	M
11-3495-4	11-3497-4	53cm-60cm	L
11-3495-5	11-3497-5	60cm-67cm	XL
11-3495-6	11-3497-6	67cm-75cm	XXL
11-3495-7	11-3497-7	75cm-81cm	XXXL

Measurement taken 15cm above mid-patella

PLAYMAKER XPERT



INDICATIONS

For mild medial/lateral support of the knee. For use in warmer environments.

CHARACTERISTICS

- Advanced high tenacity nylon/spandex material promotes a cooling effect while delivering a lighter and slimmer profile
- Removable polycentric hinges and aluminium uprights
- Popliteal cut-out for extra comfort in flexion
- Fits left or right knee
- 6 month warranty
- Latex-free

ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
82-0023-1	33cm-39cm	XS
82-0023-2	39cm-47cm	S
82-0023-3	47cm-53cm	M
82-0023-4	53cm-60cm	L
82-0023-5	60cm-67cm	XL
82-0023-6	67cm-75cm	XXL
82-0023-7	75cm-81cm	XXXL

Measurement taken 15cm above mid-patella

RIGID KNEE LIGAMENT BRACES ACCESSORIES




FULL EXTENSION STOP

PART NO.	SIZE
11-0398-9-00000	U



LYCRA UNDERGARMENT

PART NO.	SIZE
11-0075-x-00000	S-XXL



SPORTS COVER

PART NO.		SIZE
BRACE LENGTH - STANDARD	BRACE LENGTH - SHORT	
11-0016-x-06000	11-0103-x-06000	S-XXL



SUSPENSION STRAP KIT

PART NO.	SIZE
11-0543-1-06000	XS
11-0543-2-06000	S-M
11-0543-3-06000	L-XXL



NEOPRENE UNDERGARMENT

PART NO.	SIZE
11-0122-x	XS-XXL



AIR CONDYLE KIT

PART NO.	SIZE
25-1923-0-06000	U



IMPACT GUARD

PART NO.	SIZE
11-1006-x-06000	S-XXL



CONDYLE KIT

PART NO.	SIZE
25-0271-0-00000	U



FLOAM CONDYLE KIT

PART NO.	SIZE
11-1107-0-00000	U



PINCH GUARD

PART NO.		SIZE
ARMOR BRACE	DEFIANCE BRACE	
11-1095-9-06000	11-1367-9-06000	U

OSTEOARTHRITIS (OA)



We are at a unique stage in our healthcare system where patients are being diagnosed with knee OA at a younger age, and there is increasing demand for conservative treatment options. OA is the fourth leading disability in the world with 8.75 million people in the UK seeking treatment for OA*. This combined with a growing active aging population sees the healthcare system not being able to continue to offer total knee replacements as a primary intervention.

OA knee bracing is a conservative management alternative.

*Reference-Arthritis Research UK



DONJOY[®]

THE NEXT GENERATION OF
OFFLOADER BRACES FROM DONJOY

CONTACT US TO ARRANGE A
DEMONSTRATION

[ARRANGE A DEMO HERE](#)

AVAILABLE SOON



DONJOY
OA GO



AVAILABLE JANUARY



DONJOY
ROAM OA



MORE INFO HERE

NEW



DONJOY
ROAM OA™

The newest innovation in unloader braces from DonJoy

INDICATIONS

For patients with symptoms, physical findings or radiographic evidence of moderate to severe unicompartmental knee osteoarthritis, post-operative meniscal tears only where unicompartmental unloading is indicated, or varusvalgus knee conditions that may benefit from unicompartmental unloading.



CHARACTERISTICS

- Easy mag-align clips.
- Set and forget technology
- Boa®* dial for just the right fit
- Slim profile hinge
- Anti-migration technology
- Patent-pending comfort straps
- Hand Washable
- Backed by DonJoy®
- Easy to wear
- Only 6 off-the-shelf SKUs
- Covered by most insurance plans

ORDERING INFORMATION

PART NUMBER		THIGH CIRC.	KNEE CENTER	CALF CIRC.	SIZE
RIGHT MEDIAL/ LEFT LATERAL	LEFT MEDIAL/ RIGHT LATERAL				
11-1726-1-06000	11-1725-1-06000	33 - 47 cm	30.5 - 35.5 cm	25.5 - 35.5 cm	XS/S
11-1726-3-06000	11-1725-3-06000	47 - 59.5 cm	35.5 - 35.5 cm	35.5 - 45.5 cm	M/L
11-1726-5-06000	11-1725-5-06000	59.5 - 81 cm	43.2 - 58.4 cm	45.5 - 61 cm	XL/XXL/XXXL

DONJOY

OA DEFIANCE



INDICATIONS

Unicompartmental OA with or without combined ACL, PCL, Combined, MCL and LCL instabilities and activities of daily living.

CHARACTERISTICS

- Customised design, with a wide range of colour choices
- Reinforced carbon composite frame
- Durable support for moderate to severe levels of osteoarthritis.
- This polyaxial hinge comes with a telescoping condyle pad
- Allows patient-controlled additional unloading on the base load of 1 to 7° built into the frame.
- Condyle pad is included for increased patient comfort and compliance
- Suited for patients with active lifestyles or athletes.

ORDERING INFORMATION

SIZE	THIGH CIRCUMFERENCE	CALF CIRCUMFERENCE
Minimum	31.8 cm (12 1/2")	26 cm (10 1/4")
Maximum	97 cm (38 1/8")	66.5 cm (26 1/4")

Each Defiance is custom made. Measurements taken with CCMI Mark III

DONJOY

OA ADJUSTER™ 3



INDICATIONS

For patients with moderate to severe unicompartmental OA, needing significant off-load and ligament stability in a tough durable frame.

CHARACTERISTICS

- Nano MAG™ technology on the off-loading side, combined with the durability of the aluminum frame, provides lightweight function with optimal strength
- Clinically Proven-DonJoy's Adjuster technology provides adjustable off-loading with the highest level of condylar separation
- 4-Points-of-Leverage system-provides clinically proven ligament stability combined with low-profile hinges and double upright design allow a contoured, personalised fit for all levels of activity
- Available in 7 sizes, XS-XXXL
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

PART NUMBER		AFFECTED COMPARTMENT
LEFT	RIGHT	
11-1591-x	11-1590-x	Medial Compartment
11-1593-x	11-1592-x	Lateral Compartment

Note: Adjustment hinges are on opposite side of affected compartment. Sizing information is on page 19

PART NUMBER	DESCRIPTION
13-3120-0-06000	Replacement Key
25-1376-0-00000	Stop kit, extension 20°, 30°, 45° flexion 60°, 75° and full lockout



DONJOY
OA NANO™

Combining Nano MAG™, a new hi-tech material, with DonJoy's proven adjuster technology, OA Nano is the world's lightest functional OA knee brace.

INDICATIONS

For mild to moderate unicompartmental OA, the OA Nano is ideal for patients who want to stay active and move freely.



CHARACTERISTICS

- The flexibility of the magnesium frame providing an intimate fit and increased comfort with Nano MAG (on the loading side) provides strength for off-loading
- DonJoy's clinically proven three point off-loading technology shifts the compressive loads from the degenerative area to the healthy compartment

- 4 Points of Leverage and contoured frame design enhances brace suspension for all levels of activity
- Available in 7 sizes, XS-XXXL
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

PART NUMBER		AFFECTED COMPARTMENT
LEFT	RIGHT	
11-1215-x	11-1214-x	Medial Compartment
11-1217-x	11-1216-x	Lateral Compartment

Note: Adjustment hinges are on opposite side of affected compartment. Sizing information on page 15

PART NUMBER	DESCRIPTION
13-3120-0-06000	Replacement Key
25-1376-0-00000	Stop kit, extension 20°, 30°, 45° flexion 60°, 75° and full lockout

DONJOY

CLIMA-FLEX™ OA



AVAILABLE ON NHS
PRESCRIPTION



INDICATIONS

For patients with mild to moderate OA knee pain. Ideal for light activities such as walking, golf, or hiking.

CHARACTERISTICS

- Clima-Flex OA offers a new level of comfortable, cooling relief through compression and off-loading
- New DonJoy C-6 Fabric technology delivers breathable, moisture-wicking comfort that's soft to the touch, antimicrobial and incorporates temperature control
- Keyless off-loading enables on-the-fly hinge adjustment to empower patients with convenient, comforting relief from OA knee pain
- Available in 7 sizes, XS-XXXL
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	AFFECTED COMPARTMENT
11-8810-x	Right Medial/Left Lateral
11-8811-x	Left Medial/Right Lateral

Note: Off-loading hinge is on opposite side of affected compartment. Sizing information is on page 19

DONJOY

OA REACTION™ WEB



AVAILABLE ON NHS
PRESCRIPTION



INDICATIONS

Suitable for early to mild unicompartmental osteoarthritis of the knee & patella. The ideal "first-time" choice, allowing patients to get accustomed to wearing a brace and making it easier to transition into a more rigid OA brace when the time comes. OA Reaction Web is also a good option for patients wishing to use it as a lighter alternative to their usual rigid brace.

CHARACTERISTICS

- Four way stretch wrap-around design
- Silicone web technology
- Low profile, lightweight, discreet
- Single aluminium upright
- Silicone condyle pad
- Fitted undersleeve
- Available in 7 sizes, XS-XXXL
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	AFFECTED COMPARTMENT
82-7426-x	Right Medial/Left Lateral
82-7427-x	Left Medial/Right Lateral

Note: Off-loading hinge is on opposite side of affected compartment. Sizing information is on page 19

STRAPILAX KNEE



INDICATIONS

Support and compression following acute and chronic knee injury and mild knee osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Mild MCL/LCL sprains. Mild instability. Prevention of injury

CHARACTERISTICS

- Adjustable knee strapping works like functional tape to help provide added compression, stability, and a custom fit
- Anatomically contoured knee pad to aid stability of the patella, massaging effect, and to help improve proprioception
- Anatomical 3D knit design to help provide effective and targeted compression (20-36mmHg, medical standard) and support
- Stretch zones over the patella and the back of the knee to aid breathability, comfort, and freedom of movement
- Low profile with non-slip silicone bands
- Not made with natural rubber latex

ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
82-0711-0	30cm-34cm	XXS/Paediatric
82-0711-1	35cm-38cm	XS
82-0711-2	39cm-42cm	S
82-0711-3	43cm-46cm	M
82-0711-4	47cm-49cm	L
82-0711-5	50cm-53cm	XL
82-0711-6	54cm-57cm	XXL
82-0711-7	58cm-61cm	XXXL

Measurement taken 16cm above mid-patella

ROTULAX



INDICATIONS

Support and compression following acute and chronic knee injury and mild knee osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Patellofemoral pain. Mild Instability. Prevention of injury

CHARACTERISTICS

- Anatomically contoured knee pad to aid stability of the patella, massaging effect, and to help improve proprioception
- Anatomical 3D knit design to help provide effective and targeted compression (20-36mmHg, medical standard) and support
- Lightweight, breathable, elastic, soft knit for enhanced comfort
- Stretch zones over the patella and the back of the knee to aid breathability, comfort, and freedom of movement
- Low profile with non-slip silicone bands
- Not made with natural rubber latex
- Easy to clean, machine washable

ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
82-0710-1	35cm-38cm	XS
82-0710-2	39cm-42cm	S
82-0710-3	43cm-46cm	M
82-0710-4	47cm-49cm	L
82-0710-5	50cm-53cm	XL
82-0710-6	54cm-57cm	XXL
82-0710-7	58cm-61cm	XXXL

Measurement taken 15cm above mid-patella

DONJOY

FORTILAX KNEE



INDICATIONS

Support and compression following acute and chronic knee injury and mild knee osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Prevention of injury

CHARACTERISTICS

- Anatomical 3D knit design to help provide effective and targeted compression (20-36mmHg, medical standard) and support
- Lightweight, breathable, elastic, soft knit for enhanced comfort
- Stretch zones over the patella and the back of the knee to aid breathability, comfort, and freedom of movement
- Low profile
- Not made with natural rubber latex
- Easy to clean, machine washable

ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
82-0709-0	30cm-34cm	XXS/Paediatric
82-0709-1	35cm-38cm	XS
82-0709-2	39cm-42cm	S
82-0709-3	43cm-46cm	M
82-0709-4	47cm-49cm	L
82-0709-5	50cm-53cm	XL
82-0709-6	54cm-57cm	XXL
82-0709-7	58cm-61cm	XXXL

Measurement taken 15cm above mid-patella

AIRCAST

INFRAPATELLAR BAND



INDICATIONS

Osgood Schlatter disease; patella tendonitis.

CHARACTERISTICS

- An aircell focuses compression on the patella tendon, helping to decrease stress at the tibial tubercle
- Fits left or right knee
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	CIRCUMFERENCE	SIZE
08A-B	25.5cm-43cm	U

Measurement taken just below the kneecap

REACTION™ WEB



INDICATIONS

Chondromalacia patella, Osgood-Schlatter disease, patella OA, mild bicompartmental OA, quadriceps and patella tendonitis/tendinosis.

CHARACTERISTICS

- Flexible web design absorbs shocks and peak loads by dispersing energy across the knee
- The web design dynamically stabilises the patella on all sides with dual axis hinges for additional support
- Fitted under-sleeve for comfort
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
82-0215-2	33cm-46cm	XS/S
82-0215-3	47cm-59cm	M/L
82-0215-4	60cm-75cm	XL/XXL
82-0215-5	>75cm	XXXL

Measurement taken 15cm above mid-patella

ACCESSORIES	DESCRIPTION	SIZE
11-0541-x	Replacement Undersleeve	S to XXXL

Sizing information is on page 19

PATELAX



INDICATIONS

Support and compression following injury to the patellar tendon. Proprioceptive and neuromuscular stimulation through compression. Overuse injuries such as Osgood-Schlatter disease, patellofemoral pain syndrome, patellar tendinitis and patellar maltracking.

CHARACTERISTICS

- Anatomical 3D knit design helps provide effective and targeted compression (20-36mmHg, medical standard) and support
- Breathable, elastic, soft knit for enhanced wearing comfort
- Stretch zones over the patella and the back of the knee for added breathability, comfort, and freedom of movement
- Compression-reduced edges help divert pressure at the ends of the support to provide a comfortable fit
- Anatomically contoured knee pad aid stability of the patella, massaging effect, and to help improve proprioception
- Patellar strap with pad helps relieve tension on the patellar tendon

ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
82-0731-1	35cm-38cm	XS
82-0731-2	39cm-42cm	S
82-0731-3	43cm-46cm	M
82-0731-4	47cm-49cm	L
82-0731-5	50cm-53cm	XL
82-0731-6	54cm-57cm	XXL
82-0731-7	58cm-61cm	XXXL

Measurement taken 15cm above mid-patella

DONJOY

TRU-PULL® LITE



INDICATIONS

Recommended for acute mild to moderate anterior knee pain, support for patella instability, lateral subluxation and maltracking.

CHARACTERISTICS

- Patented dual durometer (silicone) buttress system provides consistent corrective force on the patella
- Bifurcated strap dynamically pulls the patella into proper alignment
- Semi-flexible removable nylon hinges included
- 21cm long support giving maximum power to pull patella into alignment
- Fabric is hypoallergenic, breathable, and lightweight
- 6 month warranty

ORDERING INFORMATION

PART NUMBER		THIGH CIRCUMFERENCE	SIZE
LEFT	RIGHT		
11-0261-1	11-0260-1	33cm-39cm	XS
11-0261-2	11-0260-2	39cm-47cm	S
11-0261-3	11-0260-3	47cm-53cm	M
11-0261-4	11-0260-4	53cm-60cm	L
11-0261-5	11-0260-5	60cm-67cm	XL
11-0261-6	11-0260-6	67cm-75cm	XXL
11-0261-7	11-0260-7	75cm-81cm	XXXL

Measurement taken 15cm above mid-patella

DONJOY

STABILAX



INDICATIONS

Support and compression following acute and chronic knee injury and mild knee osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Mild MCL/LCL sprains. Mild instability and hypermobility. Prevention of injury

CHARACTERISTICS

- Removable semi-rigid bilateral polycentric hinges
- Adjustable knee strapping works like functional tape to help provide added compression, stability, and a custom fit
- Anatomically contoured knee pad to aid stability of the patella, massaging effect, and to help improve proprioception
- Anatomical 3D knit design to help provide effective and targeted compression (20-36mmHg, medical standard) and support
- Stretch zones over the patella and the back of the knee to aid breathability, comfort, and freedom of movement
- Low profile with non-slip silicone bands
- Not made with natural rubber latex

ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
82-0712-1	35cm-38cm	XS
82-0712-2	39cm-42cm	S
82-0712-3	43cm-46cm	M
82-0712-4	47cm-49cm	L
82-0712-5	50cm-53cm	XL
82-0712-6	54cm-57cm	XXL
82-0712-7	58cm-61cm	XXXL

Measurement taken 15cm above mid-patella




DONJOY GENULAX

XXX

INDICATIONS

Support and compression following acute and chronic knee injury and mild knee osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Mild MCL/LCL sprains. Mild instability and hypermobility. Mild meniscus injuries. Prevention of injury



CHARACTERISTICS

- Removable bilateral polycentric hinges with rigid arms help offer enhanced medial/lateral stability
- Adjustable knee strapping works like functional tape to help provide added compression, stability, and a custom fit
- Anatomically contoured knee pad to aid stability of the patella, massaging effect, and to help improve proprioception

- Stretch zones over the patella and the back of the knee to aid breathability, comfort, and freedom of movement
- Low profile with non-slip silicone bands/Not made with natural rubber latex

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
82-0713-1	35cm-38cm	XS
82-0713-2	39cm-42cm	S
82-0713-3	43cm-46cm	M
82-0713-4	47cm-49cm	L
82-0713-5	50cm-53cm	XL
82-0713-6	54cm-57cm	XXL
82-0713-7	58cm-61cm	XXXL
82-0709-7	58cm-61cm	XXXL

Measurement taken 15cm above mid-patella

POST OPERATIVE
KNEE & HIP





DONJOY
X-ROM™



INDICATIONS

To aid immobilization and protected range of motion associated with ACL, PCL, LCL, and MCL surgeries and meniscal repairs

**CLICK TO VIEW
KEY FEATURES**



CHARACTERISTICS

- Quickly fit patients right out of the box
- Precisely control ROM with settings in 10° increments from -10° to 90° extension and -10° to 120° flexion
- Quick-lock design: Easily lock at 10° increments from -10° to 30°
- Strap, tab and buckle design make post-op application and adjustment fast and simple
- Easily accommodate patients from 4'11" to 6'5"

- Four independently telescoping sliders with easy push button control settings help enhance brace placement, stability and comfort while avoiding the incision and peroneal nerve
- Bendable hinge bars help enhance varus/valgus contouring
- Secure contact points
- Strong, lightweight design helps reduce brace migration
- Anti-microbial pads

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-2181-9	X-ROM™ Post-Op Knee Brace	U*

*Length adjustability: 48 to 68.3cm

**CLICK HERE FOR
FITTING ADVICE VIDEO**

DONJOY

X-ACT™ ROM KNEE



INDICATIONS

Immobilisation and protected range of motion associated with ACL, PCL, LCL and MCL surgeries and meniscal repairs.

CHARACTERISTICS

- Quick and easy to apply, offers a comfortable, contoured fit and provides precise protection during healing
- The streamlined, durable aluminium hinges provide precise protection and range of motion control in 10° increments
- Extension settings from -10° to 90°
- Flexion settings from -10° to 120°
- Easy quick-lock at 0°, 15° and 30°
- The swivelling quick-clip buckles attach easily and provide a close-fitting contour around the leg
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
-------------	-------------	------

11-2151-9	X-Act ROM Knee	U*
-----------	----------------	----

Length adjustability: 48cm to 63.5cm *Maximum thigh circumference: 71cm

DONJOY

PLAYMAKER II X-ROM



INDICATIONS

Paediatric Patients: for immobilization and protected range of motion associated with ACL, PCL, LCL, and MCL surgeries and meniscal repairs.

CHARACTERISTICS

- Precise ROM Protection. Control ROM with settings in 10-degree increments from -10 degrees to 90 degrees extension and -10 degrees to 120 degrees flexion.
- Quick-lock Design. Easily lock at 10-degree increments from -10 to 30 degrees.
- Comfortable Fit. High-tenacity (Hi10™) nylon/spandex fabric offers unmatched comfort and fit.
- Breathable Spacer Mesh Fabric. Lightweight spacer fabric on front and back panels promotes a “cool” environment support.
- Fast, Easy Application. Molded strap tabs allow quick and easy application and strap adjustments to deliver a secure and comfortable fit.
- Anti-microbial Coating. Neutralizes odor while increasing comfort and durability.
- Customizable Design. Trimmable, adjustable strap and buckle design allows off-the-shelf and customizable applications.

ORDERING INFORMATION

PART NUMBER	SIZE
-------------	------

11-0021-1	XS
-----------	----

11-0021-2	S
-----------	---



DONJOY
IMMO KNEE 3V
(3-PANEL KNEE
IMMOBILISER)



INDICATIONS

Provides post-operative knee support and immobilisation.

CHARACTERISTICS

- Aluminium bilateral stays and one posterior aluminium stay for increased stability
- Suitable for a wide variety of leg sizes
- Fits left or right leg
- Contains Latex

ORDERING INFORMATION

PART NUMBER	LENGTH	QUANTITY	SIZE
82-96105	40cm (16")	1	U
79-96105-10	40cm (16")	10	U
82-96106	50cm (20")	1	U
79-96106-10	50cm (20")	10	U
82-96107	60cm (24")	1	U
79-96107-10	60cm (24")	10	U
82-96108	65cm (26")	1	U



DONJOY
IMMO AT4 3V KNEE
IMMOBILISER



INDICATIONS

Pre or post-operative total immobilisation of the knee joint.

CHARACTERISTICS

- Lightweight (less than 600 grams) universal 3 panel immobiliser
- Malleable & removable bilateral aluminium stays and a posterior stay
- Perforated foam for enhanced comfort
- Internal anti-migration strap
- Fits left or right leg
- Contains Latex

ORDERING INFORMATION

PART NUMBER	LENGTH	SIZE
82-43-30	30cm (12")	U/Paediatric
82-43-01	40cm (16")	U
82-43-02	50cm (20")	U
82-43-03	60cm (24")	U
82-43-04	65cm (26")	U

Measure inseam of leg minus 20cm (8")



PROCARE
LEG ELEVATOR



INDICATIONS

For patients with leg casts, or those requiring stable elevation and internal/external rotation.

CHARACTERISTICS

- Full foam construction to elevate and maintain proper knee and lower leg position at 45° angle
- Supplied vacuum sealed for easy storage
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER	LENGTH x HEIGHT	SIZE
79-90191	80cm x 25.5cm	U



DONJOY
VERSAROM HIP™



INDICATIONS

DonJoy® VersaROM™ Hip is indicated for use with post-op patients who have undergone a surgical procedure and need an orthosis until the component is integrated and soft tissue heals.

CHARACTERISTICS

- Integrated belt compression with cinching straps
- Malleable hip and thigh frames
- Precise ROM Adjustments
- Range of Motion allows controlled hip flexion and extension adjustments in 15° increments
- Telescoping Height Adjustment
- Laterally placed abduction paddle with three-point pressure
- Customized Fit
- Lightweight, slim-profile design

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
11-3335	DonJoy® VersaROM™ Hip

FOOT



ACTYTOE

FOOT BRACING



FOOT

THE BEST FIT GUIDE

POST-OP SHOES		
SMART	STANDARD	VALUE
-	Remedy Pro	Post-Op Shoe
-	Squared Toe Post-Op Shoe	Med/Surg Shoe
-	Off-Loading Diabetic Shoe	Rocker Cast Boot
-	Podalux II	-

SPLINTS		
SMART	STANDARD	VALUE
AirLift PTTD Brace	AirHeel	Podous Boot
-	PodaLib	Plantar Fasciitis Night Splint
-	-	ProWedge Night Splint
-	-	Superlite A.F.O.



AIRCAST
ACTYTOE



INDICATIONS

- Pain relief in patients with symptomatic Hallux Valgus (bunion)
- Mild or moderate mal-positioning of the big toe
- Maintaining position after Hallux Valgus (bunion) surgery



CHARACTERISTICS

- Low-profile design
- Locking clip for splint immobilization (0° and 30° dorsi- and plantar flexion)
- Memory foam cushion
- Adjustable correction angle (0° and 10°)
- Adjustable, thin and soft straps
- Anatomically contoured splint
- Breathable and lightweight materials

- For day and night wear
- Suitable for walking
- Available in 3 sizes (S, M, L)
- For left and right foot
- Can be worn with wide footwear
- Not made with natural rubber latex

ORDERING INFORMATION

PART NUMBER	SHOE SIZE			SIZE
	EU	UK MEN	UK WOMEN	
82-04S	< 36	< 4	< 3.5	S
82-04M	35 - 42	3.5 - 8	3 - 7.5	M
82-04L	41 - 46	7.5 - 11.5	7 - 11	L



INDICATIONS

Reducing forefoot pressure, forefoot trauma, post-surgical healing and diabetic ulcers or wounds.

CHARACTERISTICS

- 15° wedge design reduces forefoot pressure by shifting weight to midfoot and heel
- Easy forefoot closure
- Durable outer sole, padded liner and adjustable ankle strap
- Seamless, padded heel
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE MENS		SIZE
	MENS	WOMENS	
79-81722	<4.5	<4	XS
79-81723	4.5 - 6.5	4 - 6	S
79-81725	6.5 - 8.5	6 - 8	M
79-81727	8.5 - 10.5	8 - 10	L
79-81728	>10.5	>10	XL



PROCARE
SQUARED TOE
POST-OP SHOE



INDICATIONS

Post-operative treatment of hallux valgus, toe deformities and protocols where maximum immobilisation and protection of the forefoot are required.

CHARACTERISTICS

- Squared design allows additional room for bulky dressings
- High ankle strapping for ease of ambulation and pressure reduction
- Rigid rocker sole with skid resistant sole
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE MENS		SIZE
	MENS	WOMENS	
79-81232	3 - 4.5	2 - 3.5	XS
79-81233	5 - 6.5	4 - 5.5	S
79-81235	7 - 8.5	6 - 7.5	M
79-81237	9 - 11.5	8 - 10.5	L
79-81238	>11.5	>10.5	XL



DONJOY
PODALUX II



INDICATIONS

Post-op, post-trauma foot and toe applications and for patients requiring extra stability, more specifically after hallux valgus and ensuring the foot, toes and forefoot are held in a proper position during rehabilitation.

CHARACTERISTICS

- Wide, solid skid resistant sole for extra stability
- Removable rigid foot plate for additional firmness and arch support
- Lined with breathable micro-mesh textile for patient comfort
- Two contact closure straps for easy application and adjustment
- Toe piece included to cover up the toes if necessary
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	SHOE MENS	SHOE WOMENS	SIZE
82-0042-1	3 - 4	2.5 - 3.5	S
82-0042-2	4.5 - 5.5	4 - 5	M
82-0042-3	6 - 7	6 - 6.5	L
82-0042-4	7.5 - 8	7 - 7.5	XL
82-0042-5	8.5 - 10.5	8 - 10.5	XXL
82-0042-6	11 - 13	11 - 13	XXXL

PROCARE

ROCKER CAST BOOT



INDICATIONS

For use with lower leg casted fractures.

CHARACTERISTICS

- E.V.A. (Ethyl Vinyl Acetate) rocker sole, in a sandal style with canvas upper
- Easy loop/lock fastening
- Paediatric sizes XXS and XXXS have a child friendly print on the canvas upper material
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	LENGTH	TOE WIDTH	SIZE
79-81110	11cm	6cm	XXXS
79-81111	13cm	7.5cm	XXS
79-81112	17cm	9cm	XS
79-81113	22cm	10.5cm	S
79-81114	25cm	11cm	S/M
79-81115	27cm	12cm	M
79-81116	29cm	12.5cm	M/L
79-81117	31cm	13cm	L
79-81118	33cm	15cm	XL

PROCARE

MED/SURG SHOE



INDICATIONS

Post-op, post-trauma applications where some mobility of the foot is allowed.

CHARACTERISTICS

- Loop/lock contact closure straps
- Comfortable contoured heel with foam mesh upper and semi-flexible skid resistant sole
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE MENS	SIZE
79-81132	<4	XS
79-81133	4 - 8	S
79-81135	8 - 11	M
79-81137	11 - 14	L
79-81138	>14	XL

PART NUMBER	UK SHOE SIZE WOMENS	SIZE
79-81142	<1	XS
79-81143	1 - 3.5	S
79-81145	3.5 - 5.5	M
79-81147	5.5 - 7.5	L

POST-OP SHOE



INDICATIONS

Post foot surgery, metatarsal and soft tissue injuries.

CHARACTERISTICS

- Rigid E.V.A. (Ethyl Vinyl Acetate) rocker sole and contoured heel with foam/nylon mesh upper
- Protective tongue provides maximum comfort
- Easy contact closure
- Available in female (narrower) and male (wider) sizes
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE MENS	SIZE
79-90183	6-8	S
79-90185	8-10	M
79-90187	10-12	L
79-90188	>12	XL

PART NUMBER	UK SHOE SIZE WOMENS	SIZE
79-90192	<1	XS
79-90193	1.5 - 2	S
79-90195	2.5 - 5	M
79-90197	5.5 - 8	L

OFF-LOADING DIABETIC SHOE



INDICATIONS

Pressure relief during the treatment of diabetic ulcers.

CHARACTERISTICS

- Tri-laminate IMPAX grid insole designed to redistribute load away from ulcers with easy to remove IMPAX foam blocks
- IMPAX grid allows for a continuous plantar surface
- Lightweight breathable material enhances patient's comfort and promotes compliance
- Squared toe box and wide foot bed to leave ample room for wound dressings
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE MENS	UK SHOE SIZE WOMENS	SIZE
79-81513	3.5 - 5.5	4 - 6	S
79-81515	6 - 8	6.5 - 8.5	M
79-81517	8.5 - 10.5	9 - 10	L
79-81518	11 - 13.5	>10.5	XL

ORDERING INFORMATION - REPLACEMENT INSOLE

PART NUMBER	UK SHOE SIZE MENS	UK SHOE SIZE WOMENS	SIZE
79-81523	3.5 - 5.5	4 - 6	S
79-81525	6 - 8	6.5 - 8.5	M
79-81527	8.5 - 10.5	9 - 10	L
79-81528	11 - 13.5	>10.5	XL



AIRCAST

AIRLIFT™ PTTD BRACE



INDICATIONS

Phase I and initial phase II of posterior tibial tendon dysfunction, first symptoms of flat foot.



CHARACTERISTICS

- Adjustable arch aircell for individual support and comfort
- Anatomically designed shells for secure ankle stability
- Integrated inflation pump

- Step-in design and simple two strap application
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER		UK SHOE SIZE	SIZE
LEFT	RIGHT		
02PSL-1	02PSR-1	<6	S
02PML-1	02PMR-1	6 - 8	M
02PLL-1	02PLR-1	>8	L



PROCARE PODOUS BOOT



INDICATIONS

Prevent lower extremity disorders associated with trauma or immobility including: pressure necrosis, ankle contractures and post-operative hip abduction or rotation.

CHARACTERISTICS

- Lightweight, quality plastic shell with breathable poly/pile liner
- Adjustable toe extension helps protect toes from bedding
- Anti-rotation bar may be repositioned to assist with internal or external rotation
- Latex free
- Replacement liners available as a pack of 3

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MENS	WOMENS	
79-90550	<9.5	<9	U
79-90551	>9.5	>9	XL
REPLACEMENT LINERS		PACK SIZE	SIZE
79-90560		Pack of 3	U
79-90561		Pack of 3	XL



PROCARE PLANTAR FASCIITIS NIGHT SPLINT



INDICATIONS

Controlled application of a prolonged stretch on the plantar fascia and Achilles tendon.

CHARACTERISTICS

- Lightweight, sturdy one-piece flexible shell with a terry liner
- Adjustable bilateral dorsiflexion straps
- Padded calf and ankle straps for additional comfort
- Includes optional toe wedge for added stretch of the plantar fascia
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MENS	WOMENS	
79-97753	<5	<6	S
79-97755	5 - 8	6 - 8.5	M
79-97757	8 - 10.5	8.5 - 11	L
79-97758	>10.5	11 - 13.5	XL



PROCARE PROWEDGE™ NIGHT SPLINT



INDICATIONS

Plantar fasciitis and Achilles tendonitis.

CHARACTERISTICS

- Soft, padded contoured posterior shell which provides passive dorsiflexion
- Removable liner may be washed for patient hygiene
- Includes optional toe wedge for added stretch of the plantar fascia
- Night time immobilisation and stretching
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MENS	WOMENS	
79-81403	Paediatric	Paediatric	S
79-81405	2 - 6	3 - 7	M
79-81407	6.5 - 9	7 - 10	L
79-81408	>9	>10	XL



DONJOY PODALIB®



INDICATIONS

Drop foot, neurological weakness or paralysis of the foot-lifting muscles, helps prevent lateral imbalance without restricting movement, improving ambulation.

CHARACTERISTICS

- Quick closure buckles helps to adjust and secure the tension of the elastic strap and foot flexion angle
- Adjustable elastic strap with metal hook adapts to all shoe types with laces, buckles or eyelets
- Optional foot cuff for barefoot use or swimming
- Anatomical ankle cuff with two aluminium stays for improved support
- Posterior foam pad for comfort on the Achilles Tendon area
- Contains latex

ORDERING INFORMATION

PART NUMBER	ANKLE CIRC.	ANKLE HEIGHT	SIZE
82-0013-1	13cm - 19cm	10cm	XS
82-0013-2	13cm - 20cm	13cm	S
82-0013-3	20cm - 27cm	13cm	M
82-0013-4	27cm - 34cm	13cm	L



WALKER BOOTS

AIRSELECT ELITE

WALKER BOOTS

THE BEST FIT GUIDE

	WALKERS		
	SMART	STANDARD	VALUE
WALKERS	AirSelect Elite	AirSelect Standard/Short	Nextep Contour Air Walker
	-	-	Nextep Contour Walker/ Shortie Walker
	-	-	XcelTrax Tall/Ankle

COLD THERAPY ANKLE OPTIONS

For more details see pages 92-101





AIRCAST

AIRSELECT™ ELITE



AirSelect puts comfort, control and healing in the hands of the patient with the most advanced pneumatic walking boot available. Beneath its sleek appearance are patented technological advances that work together to enhance treatment and improve outcomes. The integrated inflation system lets patients inflate each aircell individually, to maximise comfort and minimise oedema.

INDICATIONS

Specifically designed for stable fractures of the lower leg, foot, and ankle; severe ankle sprains; and post-operative use. Clinically proven effective in healing fractures and managing oedema.



CHARACTERISTICS

- The most advanced walking boot, engineered to provide the ultimate in protection, comfort, and oedema control
- Overlapping posterior and distal duplex aircells provide intermittent compression, reducing oedema 3 times faster
- Pre-inflated front panel aircell and custom inflation system provide a full-contact fit
- Enclosed toe cover provides protection from hazards and the elements
- Latex free
- Fits left or right leg
- 6 month warranty

PART NUMBER	SHOE SIZE		SIZE
	EU	UK	
01EP-XS	<35.5	<3	XS
01EP-S	35.5 - 39.5	3 - 6	S
01EP-M	39.5 - 43	6 - 9	M
01EP-L	43 - 47	9 - 12	L
01EP-XL	>47	>12	XL

LEARN MORE ABOUT AIRSELECT

CLICK FOR LINK TO FITTING VIDEO

AIRCAST**AIRSELECT™
STANDARD****INDICATIONS**

Stable fracture of foot/ankle/lower leg, severe ankle sprain, post-operative immobilisation.

CHARACTERISTICS

- Semi-rigid shell that supports the limb while providing protection
- Two adjustable distal aircells that provide compression and support
- Integrated inflation system for ease of use
- Lab-tested rocker sole encourages natural gait
- Latex free
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER	SHOE SIZE		SIZE
	EU	UK	
01EF-XS	<35.5	<3	XS
01EF-S	35.5 - 39.5	3 - 6	S
01EF-M	39.5 - 43	6 - 9	M
01EF-L	43 - 47	9 - 12	L
01EF-XL	>47	>12	XL

AIRCAST**AIRSELECT™ SHORT****INDICATIONS**

Metatarsal fractures, stress and stable fracture of the foot, bunionectomy, hallux valgus, soft tissue injury, severe ankle sprain.

CHARACTERISTICS

- The Short walking boot has a durable, semi-rigid shell that helps support the limb while providing protection
- The pre-inflated front panel and two customisable distal aircells provide compression and support
- Integrated inflation system for ease of use
- SoftStrike technology absorbs and dissipates shock, while the lightweight and lab-tested rocker sole encourages natural gait
- Latex free
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER	SHOE SIZE		SIZE
	EU	UK	
01ES-XS	<35.5	<3	XS
01ES-S	35.5 - 39.5	3 - 6	S
01ES-M	39.5 - 43	6 - 9	M
01ES-L	43 - 47	9 - 12	L
01ES-XL	>47	>12	XL

AIRCAST

XP DIABETIC WALKER™



INDICATIONS

Stable foot and/or ankle fracture; severe ankle sprain; post-operative; Charcot foot; ulcer in malleolus, hind-foot and heel.

CHARACTERISTICS

- Specifically designed to meet the needs of the higher risk diabetic patient with a high rocker sole for maximum off-loading of the diabetic foot
- Dual density insole to help eliminate pressure points
- Graduated compression to help reduce oedema
- Regulated aircell inflation with the specialised hand bulb with pressure gauge
- Latex free
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	SHOE SIZE		SIZE
	EU	UK	
01PD-S	35 - 38	4 - 5.5	S
01PD-M	39 - 42	6 - 8	M
01PD-L	>42	>8	L
01PD-XL	>45	>10.5	XL

AIRCAST

XP DIABETIC WALKER REPLACEMENT PARTS



INDICATIONS

Replacement parts for use with XP Diabetic Walker System.

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
01PG	Replacement Pump with Pressure Gauge	U
4003	Replacement Diabetic Sock	U
01R-S	Diabetic Replacement Kit (includes 1 sock and 1 insole)	S
01R-M	Diabetic Replacement Kit (includes 1 sock and 1 insole)	M
01R-L	Diabetic Replacement Kit (includes 1 sock and 1 insole)	L
01R-XL	Diabetic Replacement Kit (includes 1 sock and 1 insole)	XL
01PI-S	IMPAX Grid Insole Replacement Kit (includes 1 sock and 1 insole)	S
01PI-M	IMPAX Grid Insole Replacement Kit (includes 1 sock and 1 insole)	M
01PI-L	IMPAX Grid Insole Replacement Kit (includes 1 sock and 1 insole)	L
01PI-XL	IMPAX Grid Insole Replacement Kit (includes 1 sock and 1 insole)	XL

NEW ORDER CODE



AIRCAST
TUBE STRETCH
SOCK



INDICATIONS

To be used with Aircast walker boots to create a barrier between the boot and the skin, to absorb any sweat and decrease friction.

CHARACTERISTICS

- Cotton sock for use with AirSelect Elite, Short and Standard
- Fits left or right foot
- Latex free

ORDERING INFORMATION

PART NUMBER	SIZE
0129A-XL-TN	U



AIRCAST
HYGIENE COVER

INDICATIONS

Hygienic cover, intended for use on Aircast walking braces. This cover is intended for indoor use. It prevents surfaces from being soiled by the soles of the walker (carpets, furniture, beds). Not suitable for outdoor use.

CHARACTERISTICS

- Made of breathable nylon material
- Easy to use
- Anti slip sole

ORDERING INFORMATION

PART NUMBER	WALKING BRACE SIZE
0130A-PED/S	XS/S
0130A-M/L	M/L
0130A-XL	XL



AIRCAST
WEATHER COVER



INDICATIONS

Designed to be worn with the Aircast walking braces to protect against bad weather.

CHARACTERISTICS

- Made of durable water repellent material
- Available in 2 heights

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	WALKING BRACE SIZE	SIZE
01WXT-S		Paediatric/S	S
01WXT-L	Tall	M/L	L
01WXT-XL		XL	XL
01WXS-S		Paediatric/S	S
01WXS-L	Short	M/L	L
01WXS-XL		XL	XL



AIRCAST
10° WEDGE KIT



INDICATIONS

The Aircast 10° wedges are designed: to comfortably position the ankle joint in plantar flexion during Achilles tendon rehabilitation for use with the Aircast AirSelect Standard and Elite walker boots

CHARACTERISTICS

- Customized height adjustments - the kit contains 3 wedges:
 - 3 wedges = angle of 30° of plantar flexion
 - 2 wedges = angle of 20° of plantar flexion
 - 1 wedge = angle of 10° of plantar flexion
- Lightweight but strong hollow design to keep additional weight to a minimum
- Double-sided adhesive tape included with each wedge for easy application and secure fixation
- Bilateral: fits either foot
- Latex free
- Sustainable design made from recycled material

ORDERING INFORMATION

PART NUMBER	AIRSELECT WALKER BOOT SIZE
82-0110	XS-S-M
82-0210	L-XL

MAXTRAX XCEL ROM



INDICATIONS

- Post-operative or post-traumatic immobilization of the ankle or foot
- Severe ankle and/or foot sprain
- Conservative treatment of ankle and/or foot fractures
- Severe ankle tendinopathy
- Conservative treatment of Achilles tendon rupture

CHARACTERISTICS

- Adjustable ROM hinge: Allows ROM settings in 7.5° increments between 45° plantar-flexion and 30° dorsiflexion with locking option
- Unique pre-shaped double upright design: Modifies pressure on the malleoli and lower leg and ensures a better fit
- Edgeless rocker sole: Helps promote a natural, stable gait
- Shock-absorbing sole: Reduces the impact of heel-strike during ambulation
- Padded removable liners: For cool comfort and easy cleaning and replacement
- Adjustable hook-and-loop straps: Help reduce pressure on the injury site
- XS-XL sizes, tall and ankle versions: Accommodate a range of leg sizes and heights

ORDERING INFORMATION

PART NUMBER		SHOE SIZE	SIZE
ANKLE	TALL		
82-102-XS	82-104-XS	<3.5	XS
82-102-S	82-104-S	4 - 6.5	S
82-102-M	82-104-M	7 - 10	M
82-102-L	82-104-L	10 - 12	L
82-102-XL	82-104-XL	>12	XL



PROCARE NEXTEP™ CONTOUR WALKER



INDICATIONS

Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg or forefoot.

CHARACTERISTICS

- Secure wraparound shell housing provides immobilisation and support for ligament and tendon repair
- Eliminates need for serial casting
- Brushed nylon/foam liner with detachable toe cover
- Rocker bottom with non-skid sole
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MENS	WOMENS	
79-95063	<3 - 6	2 - 5	S
79-95065	6.5 - 9	5.5 - 8	M
79-95067	9.5 - 12.5	8.5 - 13.5	L



PROCARE NEXTEP™ CONTOUR SHORTIE WALKER



INDICATIONS

Provides immobilisation for acute ankle injuries, stable fractures or post-op forefoot surgeries of the foot.

CHARACTERISTICS

- Lightweight, durable design
- Shorter one-piece wraparound posterior shell and anterior tibial plate
- Brushed nylon/foam liner for cool comfort during rehabilitation
- Cushioned, non-skid rocker sole
- Fits left and right leg

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MENS	WOMENS	
79-95083	<3 - 6	2 - 5	S
79-95085	6.5 - 9	5.5 - 8	M
79-95087	9.5 - 12.5	8.5 - 13.5	L

PROCARE

XCELTRAX™ TALL



INDICATIONS

Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

CHARACTERISTICS

- Lightweight design
- Unique pre-shaped ergonomic frame helps to modify the level of pressure on the injury at both the malleoli and lower leg
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER		UK SHOE SIZE	SIZE
PACK OF 1	PACK OF 16		
79-95492-1		<3.5	XS
79-95493-1	79-95493-16	4 - 6.5	S
79-95495-1	79-95495-16	7 - 10	M
79-95497-1	79-95497-16	10 - 12	L
79-95498-1		>12	XL

CLICK FOR LINK TO FITTING VIDEO

PROCARE

XCELTRAX™ ANKLE



INDICATIONS

Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

CHARACTERISTICS

- Short version
- Lightweight design
- Unique pre-shaped ergonomic frame helps to modify the level of pressure on the injury at both the malleoli and lower leg
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER		UK SHOE SIZE	SIZE
PACK OF 1	PACK OF 16		
79-95502-1		<3.5	XS
79-95503-1	79-95503-16	4 - 6.5	S
79-95505-1	79-95505-16	7 - 10	M
79-95507-1	79-95507-16	10 - 12	L
79-95508-1		>12	XL

PROCARE
MINITRAX™



INDICATIONS

Ideal support and comfort for your smallest patients following acute ankle sprains, metatarsal fractures, stress and stable fractures of the lower leg.

CHARACTERISTICS

- Child friendly design elements to make wearing a boot less intimidating
- Lightweight design features anatomically shaped uprights to facilitate fit
- Low profile rocker sole features a durable, non-slip tread
- Easy to apply and adjust with a simple strapping and liner system
- Fits left or right leg

ORDERING INFORMATION

PART	AGE	FOOT	LEG	SHOE	SIZE
79-95510	1 - 2 ½	<15cm	<23cm	<7.5	S
79-95520	2 ½ - 5	15-17.5cm	<28cm	7.5 - 10	M
79-95530	5 - 8	17.5-20cm	<36cm	10.5 - 1 adult	L

PROCARE
SHOELIFT™



INDICATIONS

Simple, fast and easy to apply, the Shoelift™ uses an innovative design that helps to equalize a patient's limb length and reduce body strain while walking in a cast or walking boot.

CHARACTERISTICS

- Easily removable
- Recommended for securely fastened closed toe shoe or sneaker
- Available in three different sizes
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MENS	WOMENS	
79-95603	6.5 - 8	4 - 7	S
79-95605	8.5 - 11	7 - 9.5	M
79-95607	>12	>10	L





ANKLE

ACTYFOOT

ANKLE BRACING

For more than 30 years, medical professionals worldwide have trusted Aircast® for technological breakthroughs in patient care. Rooted in scientific methods, each Aircast product is developed using the concept of “functional management.” Aircast set a new standard of care for sprains and other injuries with patented technology and the use of graduated pneumatic compression.

THE BEST FIT GUIDE

PREVENTIVE (CHRONIC ANKLE INSTABILITY)

SMART	STANDARD	VALUE
-	A60	-
-	AirSport	-
-	ActyFoot	-
-	FortiLax Ankle	-
-	MaloLax	-

ANKLE SPRAIN **GRADE I**

SMART	STANDARD	VALUE
-	AirSport	FusioGel
-	StrapiLax Ankle	FusioLight II

ANKLE SPRAIN **GRADE II**

SMART	STANDARD	VALUE
AirSport+	ActyFoot	Air-Stirrup Universe
Air-Stirrup	Actylight	-

ANKLE SPRAIN **GRADE III**

SMART	STANDARD	VALUE
AirSelect Elite	AirSelect Standard	-

COLD THERAPY ANKLE OPTIONS

For more details see pages 92-101



CRYO/CUFF® ANKLE



ICEMAN® ANKLE WRAP





AIRCAST
ACTYFOOT™



The ActyFoot™ ankle support's modular design allows the user to control the level of ankle stability as required.

Its removable lateral stay offers increased immobilization and protection to aid recovery following injury or for added confidence during activity. Without the stay, the unique strap system provides sufficient support for sports and daily activities while allowing greater freedom of movement.

INDICATIONS

Designed to provide support, stabilization and partial immobilization of the ankle. Acute lateral ligament injuries. Acute ankle sprains and ligament ruptures. Prevention of further injuries in pre-operative chronic ankle instability. Protection of ankle joint



CHARACTERISTICS

- Three integrated straps provide proprioceptive feedback and work together with semi-rigid elements to stabilize the ankle and restrict excessive range of movement.
- Detachable lateral stay allows brace to be modified for use in both acute and chronic injury phases
- Low profile. Fits in most shoes or trainers and can be worn over a sock for additional comfort
- Can be fully opened to allow simple, secure application
- Calf strap provides additional stability to the ankle joint to help improve overall balance and control
- Cross strap further helps minimize the risk of ankle inversion
- Forefoot strap allows wearer to adjust compression as required to help manage swelling

ORDERING INFORMATION

PART NUMBER		HEEL CIRCUMFERENCE	SIZE
LEFT	RIGHT		
82-02A-SL	82-02A-SR	<30cm	S
82-02A-ML	82-02A-MR	30cm-34cm	M
82-02A-LL	82-02A-LR	>34cm	L

Measurement taken across heel (calcaneus) to front of foot (talus)

**CLICK TO SEE
ACTYFOOT IN ACTION**



AIRCAST
A60™ ANKLE
SUPPORT



INDICATIONS

Chronic instability, prophylaxis.

CHARACTERISTICS

- A 60° stabiliser guards against rollover and helps prevent sprains
- Breath-O-Prene® fabric keeps wearer cool and dry
- Sleek, anatomic design for a comfortable bulk-free fit
- Simple single strap application and adjustment
- Latex free
- Now also available in white
- 6 month warranty

ORDERING INFORMATION

PART NUMBER				SHOE SIZE	SIZE
BLACK		WHITE			
LEFT	RIGHT	LEFT	RIGHT		
02TSL	02TSR	03TSL	03TSR	<6	S
02TML	02TMR	03TML	03TMR	6 - 8	M
02TLL	02TLR	03TLL	03TLR	>8	L



AIRCAST
AIRHEEL™



INDICATIONS

Chronic and acute Achilles tendonitis, calcaneal pain, plantar fasciitis.

CHARACTERISTICS

- Pulsating pneumatic compression minimises swelling and enhances circulation
- Lightweight materials and low profile design increase comfort
- Simplified design speeds application
- Fits left or right foot
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE	SIZE
09AS	<6	S
09AM	6 - 8	M
09AL	>8	L



AIRCAST
AIRSPORT™



INDICATIONS

Acute ankle sprain, post-operative (tendon rupture and ankle fracture), chronic instability, prophylaxis.

CHARACTERISTICS

- Foam-filled aircell cushions which support the ankle
- Semi-rigid encased shell offers proven Air-Stirrup support
- Integral forefoot and shin wraps combined with cross-strap for ankle stabilisation
- Heel width adjustment for easy application
- Latex free

ORDERING INFORMATION

PART NUMBER		UK SHOE SIZE	SIZE
LEFT	RIGHT		
02MXSL	02MXSR	<3	XS
02MSL	02MSR	3 - 5	S
02MML	02MMR	6 - 8	M
02MLL	02MLR	9 - 11	L
02MXLL	02MXLR	>11	XL



AIRCAST
AIRSPORT+®



INDICATIONS

Acute ankle sprain, post-operative (tendon rupture and ankle fracture), chronic instability, prophylaxis.

CHARACTERISTICS

- Boa® lacing system provides adjustable, uniform compression without creating pressure points
- Foam-filled aircell cushions which support the ankle
- Anterior cross strap and forefoot wrap provide additional compression and stability to the joint
- Contains latex

ORDERING INFORMATION

PART NUMBER		UK SHOE SIZE	SIZE
LEFT	RIGHT		
82-0266-1	82-0265-1	<3	XS
82-0266-2	82-0265-2	3 - 5	S
82-0266-3	82-0265-3	6 - 8	M
82-0266-4	82-0265-4	9 - 11	L
82-0266-5	82-0265-5	>11	XL



AIRCAST
AIR-STIRRUP®
UNIVERSE™



INDICATIONS

Ankle sprains and chronic ankle instability.

CHARACTERISTICS

- Patented Duplex aircell system to reduce swelling and pain by enhanced circulation
- Universal semi-rigid shell design fits left or right ankle
- Universal size
- Latex free
- 6 month warranty
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	SIZE
02E	U



AIRCAST
AIR-STIRRUP®



INDICATIONS

Acute ankle distortion grade I, II & III, post-operative application, chronic ankle instability.

CHARACTERISTICS

- Semi-rigid, anatomically designed shells for protection, comfort and the prevention of inversion/eversion
- Patented Duplex aircell system to enhance circulation and reduce swelling
- Pre-inflated aircells for easy application
- Streamlined to fit in shoes for early protected weight-bearing
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER		PATIENT HEIGHT	BRACE LENGTH	SIZE
LEFT	RIGHT			
02JL	02JR	Small children	15cm	XS (Paed.)
02CL	02CR	<157cm (<5'2")	22cm	S
02BL	02BR	157cm-167cm (5'2"-5'6")	23cm	M
02AL	02AR	>162cm (>5'4")	27cm	L

DONJOY
FUSIOGEL



NEW



INDICATIONS

Designed to support the ankle for mild to moderate ankle sprains and help with chronic ankle instabilities.

CHARACTERISTICS

- Reusable cold therapy rel liners help reduce swelling
- Comfortable foam and gel liners help reduce pressure on malleoli
- Semi-rigid shells help minimize ankle rotation
- Anatomical shell design contours the malleoli for comfortable compression
- Moisture-absorbing, antibacterial Thermo-Tex™ treatment helps prevent odours
- Lightweight, low profile design
- Fits left or right ankle

ORDERING INFORMATION

PART NUMBER	XXX	XXX
82-0053-1	17cm	S
82-0053-2	23cm	M
82-0053-3	25.5cm	L

DONJOY

FUSIOLIGHT II



INDICATIONS

Designed to support the ankle for mild to moderate ankle sprains and help with chronic ankle instabilities.

CHARACTERISTICS

- Comfortable foam and gel liners help reduce pressure on malleoli
- Semi-rigid shells help minimize ankle rotation
- Anatomical shell design contours the malleoli for comfortable compression
- Moisture-absorbing, antibacterial Thermy-Tex™ treatment helps prevent odours
- Lightweight, low profile design
- Fits left or right ankle

ORDERING INFORMATION

PART NUMBER	BRACE HEIGHT	SIZE
82-0052-1	17cm	S
82-0052-2	23cm	M
82-0052-3	25.5cm	L

NEW



AIRCAST

STRAPILAX ANKLE



INDICATIONS

Support and compression following acute and chronic ankle injury and mild ankle and foot osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Grade I ankle sprains. Prevention of injury.

CHARACTERISTICS

- Adjustable ankle strapping works like functional tape to help provide added compression, stability, and a custom fit
- Two anatomically contoured ankle pads help protect, stabilize, and relieve the malleoli
- Anatomical 3D knit design to help provide effective and targeted compression (20-36mmHg, medical standard) and support
- Lightweight, breathable, elastic, soft knit for enhanced comfort
- Low profile, fits easily in a shoe
- Not made with natural rubber latex
- Easy to clean, machine washable

ORDERING INFORMATION

PART NUMBER	ANKLE CIRCUMFERENCE	SIZE
82-0716-0	17cm-18cm	XXS/Paediatric
82-0716-1	19cm-20cm	XS
82-0716-2	21cm-22cm	S
82-0716-3	23cm-24cm	M
82-0716-4	25cm-26cm	L
82-0716-5	27cm-28cm	XL
82-0716-6	29cm-30cm	XXL

Measurement taken above malleoli



DONJOY
MALOLAX™



INDICATIONS

Support and compression following acute and chronic ankle injury and mild ankle and foot osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Prevention of injury.

CHARACTERISTICS

- Two anatomically contoured ankle pads help protect, stabilize, and relieve the malleoli
- Anatomical 3D knit design to help provide effective and targeted compression (20-36mmHg, medical standard) and support
- Lightweight, breathable, elastic, soft knit for enhanced comfort
- Low profile, fits easily in a shoe
- Not made with natural rubber latex
- Easy to clean, machine washable

ORDERING INFORMATION

PART NUMBER	ANKLE CIRCUMFERENCE	SIZE
82-0715-1	19cm-20cm	XS
82-0715-2	21cm-22cm	S
82-0715-3	23cm-24cm	M
82-0715-4	25cm-26cm	L
82-0715-5	27cm-28cm	XL
82-0715-6	29cm-30cm	XXL

Measurement taken above malleoli

DONJOY
FORTILAX™ ANKLE



INDICATIONS

Support and compression following acute and chronic ankle injury and mild ankle and foot osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Prevention of injury.

CHARACTERISTICS

- Anatomical 3D knit design to help provide effective and targeted compression (20-36mmHg, medical standard) and support
- Lightweight, breathable, elastic, soft knit for enhanced comfort
- Low profile, fits easily in a shoe
- Not made with natural rubber latex
- Easy to clean, machine washable

ORDERING INFORMATION

PART NUMBER	ANKLE CIRCUMFERENCE	SIZE
82-0714-0	17cm-18cm	XXS/Paediatric
82-0714-1	19cm-20cm	XS
82-0714-2	21cm-22cm	S
82-0714-3	23cm-24cm	M
82-0714-4	25cm-26cm	L
82-0714-5	27cm-28cm	XL
82-0714-6	29cm-30cm	XXL

Measurement taken above malleoli

SPINE



SPINOSTRAP

SPINE BRACING



THE BEST FIT GUIDE

SPINAL SUPPORT		
SMART	STANDARD	VALUE
IsoFORM Postural Extension	Actistrap	Premium Panel Elastic Binders
IsoFORM TLSO	Immostrap	Maternity Belt
Spinostrap	IsoFORM SIO	Premium Panel Elastic Binders
-	IsoFORM LO/LO+	-
-	IsoFORM LSO/LSO+	-

NECK		
SMART	STANDARD	VALUE
-	XTEND 174	Transitional 172
-	C1 Collar	Adams Collar

SPINE

COLD THERAPY BACK OPTIONS

For more details see pages 92-101



CRYO/CUFF™ BACK/HIP/RIB

PROCARE

XTEND 174 CERVICAL COLLAR



INDICATIONS

Cervical spine immobilisation.

CHARACTERISTICS

- Helps provide substantial motion restriction while helping to minimise painful pressure points
- Enhanced patient comfort using soft, anti-microbial, breathable foam padding
- Collar comes with a set of replacement pads

ORDERING INFORMATION

PART NUMBER	MODEL	COLLAR	NECK CIRC.	SIZE
79-83206	Paed 0-6 mths			Paed 0
79-83227	Paed 6mths-2 yrs	2.9cm	≤ 23cm	Paed 1
79-83228	Paed 2-6 yrs	3.2cm	20cm - 28cm	Paed 2
79-83222	Paed 6-12yrs	4.4cm	23cm - 33cm	Paed 3
79-83221	Adult	2.5cm	25cm - 50cm	Super Short
79-83232	Adult	3cm	20cm - 35cm	XSmall
79-83233	Adult	5cm	25cm - 50cm	Short
79-83235	Adult	6.3cm	25cm - 50cm	Regular
79-83237	Adult	7.5cm	25cm - 50cm	Tall
79-83238	Adult	4.5cm	35cm - 60cm	Stout

PROCARE

TRANSITIONAL 172 CERVICAL COLLAR



INDICATIONS

Restricts cervical spine flexion, extension and rotation to promote patient recovery.

CHARACTERISTICS

- Enhanced comfort using soft, breathable foam padding helps promote and encourage patient compliance
- Removable open cell foam pads wick moisture away from the skin to ensure patient hygiene and comfort
- X-ray, CT and MRI lucent
- Collar comes with a set of replacement pads

ORDERING INFORMATION

PART NO.	MODEL	COLLAR	NECK CIRC.	SIZE
79-83272	Ped 1 - 18 mths	1.9cm	10.5cm - 30.5cm	PED-1
79-83273	Ped 9 - 24 mths	3.2cm	28.0cm - 33.0cm	PED-2
79-83275	Ped 1 - 3 yrs	5.0cm	29.2cm - 35.6cm	PED-3
79-83277	Ped 2 - 5 yrs	5.7cm	30.5cm - 36.8cm	PED-4
79-83278	Ped 3 - 6 yrs	6.4cm	33.0cm - 40.6cm	PED-5
79-83363	Adult	5.5cm	33cm - 55cm	S
79-83365	Adult	7.5cm	33cm - 55cm	M
79-83367	Adult	9.5cm	33cm - 55cm	L
79-83368	Adult	11.5cm	33cm - 55cm	XL
79-83238	Adult	4.5cm	35cm - 60cm	Stout

PROCARE

ADAMS COLLAR



INDICATIONS

Soft tissue damage.

CHARACTERISTICS

- Polystyrene beads covered by 100% cotton, natural colour stockinette for comfortable and secure fit
- Loop contact closure for ease of use
- Polystyrene beads permit a quick wash and dry
- An optional rigid plastic stay is included with the collar to provide additional support
- Sold as a pack of 10

ORDERING INFORMATION

PART NUMBER	NECK CIRCUMFERENCE	SIZE
79-83263	27cm -34cm	S
79-83265	31cm-38cm	M
79-83267	35cm-42cm	L
79-83260	Mixed pack (contains 2 x small, 5 x medium and 3 x large collars)	

PROCARE

PREMIUM PANEL ELASTIC BINDER



INDICATIONS

Compression and support for post natal or post abdominal surgery, abdominal strains and weakness.


CHARACTERISTICS

- Constructed with panel elastic and soft cotton type lining for patient comfort
- Universal style support with contact closure
- Multi panel design helps prevent panel rolling
- Available in 22.5cm (9") 3 panel and 30cm (12") 4 panel heights

ORDERING INFORMATION

PART NUMBER		WAIST CIRCUMFERENCE
9" 3-PANEL BINDER	12" 4-PANEL BINDER	
79-89070	79-89090	76-114 cm (30-45")
79-89071	79-89091	114-157 cm (45-62")

ANATOMICAL RANGE OF SPINAL SUPPORT



Designed to aid superior lower back protection and stabilization, the new DonJoy® IsoFORM™ back braces help provide support and pain relief for a comprehensive range of chronic and acute indications.

ISOFORM FEATURE

Adjustable pull tabs help provide targeted compression and help reduce rotation of brace. The pull tab system independently tightens top and bottom allowing patients to localise compression.

ISOFORM FEATURE

Low profile design is easy to apply and remove. The semi-universal sizing and adjustable belt wings allow clinicians to easily customize the fit for a wide range of patients.

GENERAL INDICATIONS & RANGE OF SPINAL SUPPORT



DONJOY BACK BRACE PRODUCT LINE	S1-S5	L1-L5	L1-L5	T9-S1	T9-S1	T2/T3-S1	T2/T3-S1
INDICATIONS FOR USE	ISOFORM SIO	ISOFORM LO	ISOFORM LO+	ISOFORM LSO	ISOFORM LSO+	ISOFORM POSTURAL EXTENSION	ISOFORM TLSO
SACROILIAC PAIN	○						
LOW BACK PAIN	○	○	○	○	○	○	○
LOW BACK SPRAINS/STRAINS	○	○	○	○	○	○	○
LUMBAR DISC DISPLACEMENT		○	○	○	○	○	○
OSTEOPOROSIS		○	○	○	○	○	○
DISC HERNIATION & DEGENERATION		○	○	○	○	○	○
POST-OPERATIVE LAMINECTOMY		○	○	○	○	○	○
POST-OPERATIVE DISECTOMY		○	○	○	○	○	○
IDET PROCEDURE		○	○	○	○	○	○
SPONDYLOLISTHESIS		○	○	○	○	○	○
POST-OPERATIVE FUSION				○	○	○	○
SPINAL STENOSIS				○	○	○	○
COMPRESSION FRACTURES						○	○
KYPHOSIS						○	○
SCOLIOSIS							○



DONJOY

ISOFORM™ BACK BRACES

Designed to aid superior lower back protection and stabilization, DonJoy® IsoFORM™ back braces help provide support and pain relief for a comprehensive range of chronic and acute indications. Supplement support with hot/cold packs for additional relief.



ISOFORM™ SIO



INDICATIONS

Flexible, low-profile belt for SI stabilization

CHARACTERISTICS

- Single ergonomic pull tab compression system helps reduce SI joint motion
- Lightweight design for comfortable extended wear
- Comfortably fits waist sizes from 71-152 cm (28-60") and up to 168 cm (66")

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-1692-2	IsoFORM SIO	S/M (Hip 71-122cm)
11-1692-4	IsoFORM SIO	L/XL (Hip 122-152cm)



ISOFORM™ LO



INDICATIONS

Low-profile design to help with minor lower back pain

CHARACTERISTICS

- Targeted compression helps support fatigued muscles that spasm and cause pain
- Easy to apply and remove
- Helps provide sagittal control and support from L1 to L5
- Comfortably fits waist sizes from 71-122 cm (28-48") and up to 168 cm (66")

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-1693-2	IsoFORM LO	S/M (Hip 71-122cm)
11-1693-4	IsoFORM LO	L/XL (Hip 122-152cm)
11-1699	Hot/Cold Pack Accessory	-



ISOFORM™ LO+



INDICATIONS

Similar design as the LO option with an anterior panel to help with more trunk support

CHARACTERISTICS

- Acute and chronic low back pain, lower back sprains or strains, osteoporosis, spondylolisthesis

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-1694-2	IsoFORM LO+	S/M (Hip 71-122cm)
11-1694-4	IsoFORM LO+	L/XL (Hip 122-152cm)
11-1699	Hot/Cold Pack Accessory	-



ISOFORM™ LSO



INDICATIONS


Trunk stability with motion restriction and pain relief aid

CHARACTERISTICS

- Limited lumbar flexion and extension helps restrict movement
- Targeted compression helps support fatigued muscles that spasm and cause pain
- Low-profile brace for comfortable extended wear
- Easy to apply and remove
- Helps provide superior sagittal control and support from T9-S1
- Comfortably fits waist sizes from 71-122 cm (28-48") and up to 168 cm (66")

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-1695-2	IsoFORM LSO	S/M (Hip 71-122cm)
11-1695-4	IsoFORM LSO	L/XL (Hip 122-152cm)
11-1699	Hot/Cold Pack Accessory	-



ISOFORM™ LSO+



INDICATIONS

Similar design as the LSO option with two lateral extensions to help reduce lateral movement

CHARACTERISTICS

- Helps provide higher degree of spinal support and relief
- Helps provide superior sagittal and coronal control and support from T9-S1

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-1696-2	IsoFORM LSO+	S/M (Hip 71-122cm)
11-1696-4	IsoFORM LSO+	L/XL (Hip 125-152cm)
11-1699	Hot/Cold Pack Accessory	-

DONJOY

ISOFORM TLSO



NEW



INDICATIONS

Suitable for people with compression fractures of the thoracic spine, osteoporosis, thoracic mechanical back pain and kyphosis

CHARACTERISTICS

- Anterior panel - helps restrict flexion and trunk rotation and supports T3-S1
- Posterior panel – rigid paneling helps restrict extension, trunk rotation and supports T3-S1
- Lateral Extensions - Removable semi-rigid supports help reduce lateral movement
- Belt Wings – with removable semi-rigid paneling encompasses the torso for customizable support
- Easy to apply and remove with pivoting quick-release buckles
- Over-the shoulder or underarm strap configuration options

ORDERING INFORMATION

PART NUMBER	HIP CIRCUMFERENCE	SIZE
11-1698-2	71 - 122cm	S/M
11-1698-4	122 - 152cm	L/XL
11-1699	Hot/Cold Pack Accessory	
11-1700	Small Sternal Bar 9"	XS/S
11-1701	Large Sternal Bar 15"	XL/XXL

FITTING GUIDE

DONJOY

ISOFORM POSTURAL EXTENSION



INDICATIONS

Suitable for people with thoracic spine conditions such as kyphosis, osteoporosis, acute pain, mechanical back pain, compression fractures, strain and mild scoliosis.

CHARACTERISTICS

- Adjustable Belt Wings – accommodates different shapes for a custom fit
- Paneling – provides durable rigid support
- Dual Pulley System – provides targeted compression with independent top and bottom compression
- Flexible Pull Tabs – for ease of use with adjustable lace length feature
- Donning Mitt – assists with ease of application
- Semi-Universal – comfortably fits waist sizes from 28-60"
- Hexagonal Spacer Mesh – durable and breathable materials for comfort
- Optional Hot/Cold Pack – for additional relief
- EVA lined Panels – provides comfortable support
- Enclosed Lacing

ORDERING INFORMATION

PART NUMBER	HIP CIRCUMFERENCE	SIZE
11-1697-2	71 - 122cm	S/M
11-1697-4	122 - 152cm	L/XL
11-1699	Hot/Cold Pack Accessory	

FITTING GUIDE



DONJOY
ACTISTRAP



INDICATIONS

Acute, subacute or chronic lower back pain, lumbar sciatica.

CHARACTERISTICS

- Strong compressive elastic material to enhance posture control with double bilateral straps to adjust compression level
- Four rigid dorsal stays aid stabilisation of the sacro-lumbar area and lordosis correction
- Anatomical design for improved comfort and fit
- Two flexible abdominal stays
- Brace height 26cm (10")

ORDERING INFORMATION

PART NUMBER	WAIST CIRC.	SIZE
82-0231-1	56cm - 68cm (22" - 26.5")	XS
82-0231-2	68cm - 82cm (26.5" - 32")	S
82-0231-3	82cm - 98cm (32" - 38.5")	M
82-0231-4	98cm - 116cm (38.5" - 45.5")	L
82-0231-5	116cm - 136cm (45.5" - 53.5")	XL



DONJOY
IMMOSTRAP™



INDICATIONS

Acute or subacute lower back pain, lumbar arthrosis, lumbar sciatica, post disc herniation.

CHARACTERISTICS

- Quadruple bilateral straps to adjust compression level
- Independent superior and inferior compression adjustment for optimal anatomical fit
- Ergonomic front panel and a strong compressive elastic material enhances posture control
- Four semi rigid dorsal stays for stabilisation of the sacro-lumbar area
- Rolling buckles for effortless tightening
- Brace height 26cm (10")

ORDERING INFORMATION

PART NUMBER	WAIST CIRC.	SIZE
82-0229-1	56cm - 68cm (22" - 26.5")	XS
82-0229-2	68cm - 82cm (26.5" - 32")	S
82-0229-3	82cm - 98cm (32" - 38.5")	M
82-0229-4	98cm - 116cm (38.5" - 45.5")	L
82-0229-5	116cm - 136cm (45.5" - 53.5")	XL



DONJOY
SPINOSTRAP



INDICATIONS

DonJoy Spinostrap is designed to provide support and stabilization and partial immobilization (restriction of movement) of the spine.

CHARACTERISTICS

- Semi-rigid exoskeletal dorsal frame to provide support and aid relief of the thoracic spine, aluminum core can be molded to fit the individual curvature of the back
- Adjustable lumbar belt uses accessible strap tabs and finger pockets to provide adjustable support for the spine while leaving the abdomen free. Fixes easily using hook and loop fastening
- Adjustable shoulder straps use fixable hook-and-loop strap tabs to provide adjustable tension, aiding correction of kyphotic posture by pulling the shoulders back and straightening the back
- PVC and Latex free

ORDERING INFORMATION

PART NUMBER	SPINE LENGTH C7-S1	SIZE
82-0057-1	42cm-52cm (16½"-20½")	S
82-0057-2	52cm-62cm (20½"-24½")	M
82-0057-3	62cm-72cm (24½"-28½")	L



PROCARE
MATERNITY BELT



INDICATIONS

Relieve discomfort of pregnancy due to stretching and expansion of the abdomen.

CHARACTERISTICS

- 20cm (8") plush elastic belt helps relieve discomfort of the lower abdomen and back during all stages of pregnancy
- Abdominal strap to support stomach
- Elastic back provides support for the lumbar area


ORDERING INFORMATION

PART NUMBER	HIP CIRCUMFERENCE	SIZE
79-89293	80cm - 105cm (32" - 42")	S
79-89295	105cm - 130cm (42" - 52")	M
79-89297	130cm - 155cm (52" - 62")	L
79-89298	155cm - 180cm (62" - 72")	XL

Measurement taken during pregnancy

SHOULDER & ARM

ULTRASLING PRO



SHOULDER & ARM BRACING

THE BEST FIT GUIDE

SHOULDER		
SMART	STANDARD	VALUE
UltraSling Quadrant	Ultrasling PRO	Sully Shoulder Stabiliser
-	Humeral Fracture Brace/Over Shoulder	Shoulder Immobiliser with Foam Straps
-	Humeral Stabilising System (HSS)	Clinical Shoulder Immobiliser
-	Universal Clavicle	-

ARM		
SMART	STANDARD	VALUE
X-Act ROM Elbow	Surround Tennis Elbow	Pneumatic Armband
Mayo Clinic Elbow Brace	Healthcare Bear Arm Sling	Chieftan Arm Sling
-	CondiLax	Universal Arm Sling
-	-	Blue Vogue Arm Sling

COLD THERAPY SHOULDER OPTIONS

For more details see pages 92-101



DURA*SOFT™ SHOULDER WRAP



CRYO/CUFF™ SHOULDER



ICEMAN® SHOULDER WRAP

COLD THERAPY ARM OPTIONS

For more details see pages 92-101



CRYO/CUFF™ ELBOW



ICEMAN® WRAP ON COLD PAD

DONJOY

SULLY® SHOULDER STABILISER



INDICATIONS

Shoulder immobilisation and controlled range of motion for anterior, multi-directional, inferior, and posterior instabilities, rotator cuff deceleration, shoulder separations, and muscle strains.

CHARACTERISTICS

- Stabilises, assists and restricts movement of the shoulder post-operatively
- Can pull the humerus posteriorly as it restricts external rotation and abduction
- Elastic straps attach with Velcro at any point, in any direction for improved patient fit
- Smooth, controlled assistance or restraint follows the natural movement of the muscles and joints
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	CHEST CIRC.	BICEP CIRC.	SIZE
11-0525-1	86.5cm - 101.5cm	22.5cm - 26.5cm	XS
11-0525-2	91.5cm - 101.5cm	28cm - 33cm	S
11-0525-3	96.5cm - 106.5cm	33cm - 38cm	M
11-0525-4	106.5cm - 116.5cm	35.5cm - 43cm	L
11-0525-5	116.5cm - 129.5cm	38cm - 43cm	XL

Measurement taken at bicep in non-flex position



PROCARE HUMERAL FRACTURE BRACE/OVER THE SHOULDER



DONJOY HUMERAL STABILISING SYSTEM (HSS)



INDICATIONS

Indicated for management of mid-shaft and distal third humeral fractures.

CHARACTERISTICS

- Pre-formed orthopaedic grade thermoplastic cast-like construction for the humerus
- Hook and loop straps for easy adjustment
- Deltoid cap extension, providing a greater alignment control and prevents slippage
- Supplied with two double thickness stockinettes
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	BICEP CIRC.	LENGTH	SIZE
79-97953	20.5cm-28cm	25.5cm	S
79-97955	25.5cm-33cm	26cm	M
79-97957	30.5cm-38cm	28.5cm	L
79-97958	35.5cm-43cm	32.5cm	XL

INDICATIONS

Immobilisation of stable midshaft humeral fractures, stable supracondylar fractures and biceps/triceps repairs.

CHARACTERISTICS

- Lightweight design offers three point stabilisation
- 900g weight for axial traction
- Full circumference shells and padding with the ability to break down the upper arm cuff to suit all patients
- Alternative to serial casting

ORDERING INFORMATION

PART NUMBER	ARM	SIZE
11-0113-9-13066	Left	U
11-0112-9-13066	Right	U

PROCARE

SHOULDER IMMOBILISER WITH FOAM STRAPS



INDICATIONS

Post surgery for rotator cuff injuries, shoulder (sub)luxations and for other shoulder joint pathology.

CHARACTERISTICS

- Made of a soft cotton/poly material for patient comfort
- The waist and shoulder straps are a padded foam laminated design for patient comfort with contact closure
- Exerts upward pressure to immobilise arm and shoulder
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	DEPTH x LENGTH	SIZE
79-84162	15cm x 28cm	XS
79-84163	18cm x 33cm	S
79-84165	20.5cm x 38cm	M
79-84167	23cm x 46cm	L
79-84168	25.5cm x 51cm	XL

PROCARE

CLINICAL SHOULDER IMMOBILISER



INDICATIONS

Preventing shoulder rotation, post trauma after shoulder (sub)luxations and for over use shoulder injuries.

CHARACTERISTICS

- High quality, durable poly cotton sling
- Contact closure adjustment on the shoulder and body straps
- Shoulder pad to provide additional comfort
- Selected sizes available as a bulk pack of 25 slings
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER		DEPTH x LENGTH	SIZE
PACK OF 1	PACK OF 25		
79-84011	-	12cm x 23cm	XXS
79-84012	-	15cm x 28cm	XS
79-84013	24-84013	18cm x 33cm	S
79-84015	24-84015	20.5cm x 38cm	M
79-84017	24-84017	23cm x 46cm	L
79-84018	-	25.5cm x 51cm	XL

DONJOY

ULTRASLING QUADRANT™



INDICATIONS

Post-operative immobilisation in 15° of abduction or greater for shoulder hemiarthroplasty, total shoulder replacement, shoulder reconstruction, superior labral repair and shoulder debridement, protocols involving stabilisation following fractures to the proximal humerus, elbow and forearm, Bankart procedures, rotator cuff repair, biceps tendon surgery, shoulder dislocation and elbow ligament and tendon procedures.

CHARACTERISTICS

- TempGuard™ - Cool to touch, anti-microbial material and optional padded shoulder strap for extended duration of use
- Optimal versatility provides four angles of immobilisation in one brace; internal, external, abduction or adduction
- ProTherapy™ - Precision based system which allows the patient to flex and extend their arm, promoting therapeutic movement throughout recovery
- Elbow opening to eliminate pressure on the olecranon bursa
- Patient friendly, single handed application for patient compliance and optional shoulder strap
- Optional extension belt available separately

ORDERING INFORMATION

PART NUMBER	ARM	WAIST CIRCUMFERENCE
11-0211	Left	61cm - 114cm
11-0210	Right	61cm - 114cm
25-0212	Extension Belt	114cm - 152cm

Measurement taken at waist circumference 1" (2.5cm) above belly button

DONJOY

ULTRASLING PRO™



INDICATIONS

Post-operative rotator cuff repairs. Post-operative Bankart procedures. Capsular shifts. Glenohumeral dislocation/subluxation. Soft tissue strains/repairs

CHARACTERISTICS

- Universally designed to fit the right or left arm
- Low profile, adjustable ring easily slides over the non-injured shoulder, helping to evenly distribute weight and take pressure off neck
- Easy, one-handed application
- Post-op step-down: waist belt straps easily adhere directly to sling when pillow is removed
- UltraSling PRO provide 15° abduction without rotation
- UltraSling PRO ER provides 15° abduction with 15° or 30° internal or external rotation
- UltraSling PRO AB provides 45° or 60° abduction without rotation

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-0447-9	UltraSling PRO	Universal
11-0471-9	UltraSling PRO ER/IR 15°	Universal
11-0472-9	UltraSling PRO ER/IR 30°	Universal
11-0473-9	UltraSling PRO AB	Universal

FITTING GUIDE



DONJOY

CONDILAX™



INDICATIONS

Support and compression following elbow injury and mild elbow osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Prevention of injury. Elbow strains or sprains. Medial or lateral epicondylitis. Tendinopathies



CHARACTERISTICS

- Epicondyle pads help provide support and pressure relief for the elbow
- Anatomical 3D knit design to help provide effective and targeted compression (20-36mmHg, medical standard) and support
- Removable and adjustable strapping helps to control the degree of stabilization and compression
- Lightweight, breathable, elastic, soft knit for enhanced comfort
- Stretch zones in anterior elbow area to aid breathability, comfort, and freedom of movement
- Not made with natural rubber latex
- Easy to clean, machine washable

ORDERING INFORMATION

PART NUMBER	ELBOW CIRC.	SIZE
82-0720-1	18cm-20cm	XS
82-0720-2	21cm-24cm	S
82-0720-3	25cm-28cm	M
82-0720-4	29cm-31cm	L
82-0720-5	32cm-35cm	XL
82-0720-6	36cm-39cm	XXL

Measurement across elbow when arm at 90°

DONJOY

X-ACT ROM ELBOW



INDICATIONS

Fixed or controlled range of motion post injury or post operatively, stable fractures, biceps repair, mild to severe soft tissue injuries.

CHARACTERISTICS

- User-friendly range of motion control with flexion settings from -10° to 120° in 10° increments and extension settings from -10° to 90° in 10° increments
- Quick-lock function from 0° to 90° in 15° increments
- Cuffs telescope independently allowing a precise fit on any patient, adjustability from 36cm to 46cm
- Durable lateral aluminium hinge maximises protection and minimises brace weight
- Unique semi-rigid thermoformed and mouldable liners with anti-migration properties

ORDERING INFORMATION

PART NUMBER	ARM	SIZE
11-9122	Left	U
11-9121	Right	U

AIRCAST

MAYO CLINIC ELBOW BRACE



INDICATIONS

Bicep/tricep tendon rupture, ulnar nerve transposition, total elbow arthroplasty, ligament repair, radial fixation, osteoarthritis and epicondylitis.

CHARACTERISTICS

- Lock/unlock mechanism enables a rapid switch between immobilisation, static stretch and free range of motion
- Calibrated control knob is easy to read and assists in monitoring therapy progress
- Repositionable hinge stops control range of motion
- Removable control knob for day to day patient comfort

ORDERING INFORMATION

PART NUMBER	ARM	SIZE
05EL	Left	U
05ER	Right	U



PROCARE
CHIEFTAIN ARM
SLING



INDICATIONS

Cast support, injuries to the arm, wrist, hand or humeral fractures.

CHARACTERISTICS

- Breathable Tietex mesh material with hook and loop closure
- Thumb loop prevents wrist drop
- Supplied as a bulk pack of 12 slings
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	DEPTH x LENGTH	QUANTITY	SIZE
79-84172	16.5cm x 33cm	Pack of 12	XS
79-84173	17.5cm x 40.5cm	Pack of 12	S
79-84175	19cm x 43cm	Pack of 12	M
79-84177	20.5cm x 49.5cm	Pack of 12	L
79-84178	21.5cm x 54.5cm	Pack of 12	XL



PROCARE
BLUE VOGUE ARM
SLING



INDICATIONS


Cast support, injuries to the arm, wrist or hand.

CHARACTERISTICS

- Generous size envelope constructed of durable poly/cotton material provides the patient with comfortable, custom fit whilst ensuring immobility
- Fully padded 5cm wide foam strap with contact closure and D-rings for easy adjustment and application
- Thumb loop helps prevent sling from migrating
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	DEPTH x LENGTH	SIZE
79-84153	19cm x 32cm	S
79-84155	20cm x 38cm	M
79-84157	21.5cm x 44.5cm	L
79-84158	22.5cm x 50.5cm	XL



PROCARE
HEALTHCARE BEAR
ARM SLING



INDICATIONS

Injuries to the arm, wrist or hand of paediatric patients.

CHARACTERISTICS

- Made from a cotton/poly envelope with a cotton web strap with quick release buckle for easy application
- Paediatric friendly, bear print design on the sling
- Suitable for use in conjunction with a cast
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	DEPTH x LENGTH	SIZE
79-99132	11.5cm x 28cm	XS / Paediatric
79-99133	16.5cm x 33cm	S / Child



AIRCAST
PNEUMATIC
ARMBAND



INDICATIONS

Epicondylitis (tennis elbow), medial epicondylitis (golfer's elbow).

CHARACTERISTICS

- Single pre-inflated aircell concentrates compression on the extensor muscle, providing more support, less constriction
- Latex free
- Fits left or right elbow

ORDERING INFORMATION

PART NUMBER	FOREARM CIRC.	SIZE
05A-B	20cm - 35.5cm	U

WRIST & HAND



MANULAX

WRIST & HAND

THE BEST FIT GUIDE

WRIST & HAND		
SMART	STANDARD	VALUE
-	ManuLax	ComfortForm Wrist
-	ErgoForm	ComfortForm Wrist with Abducted Thumb
-	RhizoForm	Universal Wrist-O-Prene
-	Wrist Wraps	Universal Thumb-O-Prene
-	-	Elastic Wrist Brace
-	-	Curved Aluminium Finger Splints
-	-	Staxx Plastic Finger Splints

WRIST & HAND

COLD THERAPY WRIST & HAND OPTIONS

For more details see pages 92-101



CRYO/CUFF™ HAND

DONJOY
MANULAX



INDICATIONS

Support and compression following mild wrist injury and mild wrist osteoarthritis. Proprioceptive and neuromuscular stimulation through compression. Prevention of injury. Mild wrist sprains. Wrist tendinitis and mild tenosynovitis. Wrist instability. Carpal tunnel syndrome.

CHARACTERISTICS

- Removable and adjustable strapping helps to control the degree of stabilization and compression
- Adjustable hand stay helps provide additional stability and support, and can be removed after the acute phase
- Styloid pad helps provide support and pressure relief
- Anatomical 3D knit design to help provide effective and targeted support
- Lightweight, breathable, elastic, soft knit for enhanced comfort
- Not made with natural rubber latex
- Easy to clean, machine washable

ORDERING INFORMATION

PART NUMBER		WRIST CIRCUMFERENCE	SIZE
LEFT	RIGHT		
82-0719-1	82-0718-1	13cm-14.5cm	XS
82-0719-2	82-0718-2	14.5cm-16cm	S
82-0719-3	82-0718-3	16cm-17.5cm	M
82-0719-4	82-0718-4	17.5cm-19cm	L
82-0719-5	82-0718-5	19cm-20.5cm	XL
82-0719-6	82-0718-6	20.5cm-22cm	XXL

Measurement taken around wrist

DONJOY
ERGOFORM™



INDICATIONS

All indications requiring immobilisation of the thumb joints. Indicated for sprains fo the thumb, subluxations, dislocations, rhizoarthrosis.

CHARACTERISTICS

- Easy application via adjustable strap and "Stop Strap" system
- Malleable aluminium thumb stay
- Fits left or right hand
- Advanced Ozone™ material promotes airflow
- Lightweight, low-profile design

ORDERING INFORMATION

PART NUMBER	WRIST CIRC.	SIZE
82-96190	12cm - 16cm	S
82-96191	14cm - 18cm	M
82-96192	18cm - 22cm	L

DONJOY

COMFORTFORM™ WRIST



INDICATIONS

Treatment of sprains, strains, and control of carpal tunnel syndrome symptoms.

CHARACTERISTICS

- Durable, lightweight Lycra lined foam laminate for breathability and patient comfort
- Dual mouldable aluminium stays support the palmar surface
- Contact closure straps

ORDERING INFORMATION

PART NUMBER		WRIST CIRC.	BRACE LENGTH	SIZE
LEFT	RIGHT			
81-97892	81-97882	11.5cm - 14cm	19cm	XS
81-97893	81-97883	14cm - 16.5cm	19cm	S
81-97895	81-97885	16.5cm - 19cm	20cm	M
81-97897	81-97887	19cm - 21.5cm	21.5cm	L
81-97898	81-97888	>21.5 cm	21.5cm	XL

DONJOY

COMFORTFORM™ WRIST WITH ABDUCTED THUMB



INDICATIONS

Treatment of sprains, strains, scaphoid injuries, carpal tunnel syndrome, Gamekeeper's thumb and De Quervain's syndrome symptoms.

CHARACTERISTICS

- Durable, lightweight Lycra lined foam laminate for breathability and patient comfort
- Dual mouldable aluminium stays support the palmar surface and the extensor surface of the thumb
- Contact closure straps

ORDERING INFORMATION

PART NUMBER		WRIST CIRC.	BRACE LENGTH	SIZE
LEFT	RIGHT			
81-87312	81-87302	11.5cm - 14cm	19cm	XS
81-87313	81-87303	14cm - 16.5cm	19cm	S
81-87315	81-87305	16.5cm - 19cm	20cm	M
81-87317	81-87307	19cm - 21.5cm	21.5cm	L
81-87318	81-87308	>21.5 cm	21.5cm	XL



PROCARE
UNIVERSAL
WRIST-O-PRENE



INDICATIONS

Carpal tunnel syndrome and wrist sprains and strains.

CHARACTERISTICS

- Wraparound neoprene brace provides compression and support
- Removable contoured aluminium palmar stay
- Inside contact closure allows for ease of application and proper fit

ORDERING INFORMATION

PART NUMBER		WRIST CIRC.	SIZE
LEFT	RIGHT		
79-82471	79-82470	<28cm	U



PROCARE
UNIVERSAL
THUMB-O-PRENE™



INDICATIONS

Tendonitis, basel joint arthritis and instability of the carpometacarpal joint.

CHARACTERISTICS

- Neoprene construction provides compressive wrist support with thumb abduction
- Unique design combines two stays, rigid and spiral plus a movable elastic strap to address various pathologies of the thumb
- One size fits most patients
- Fits left or right hand

ORDERING INFORMATION

PART NUMBER	SIZE
79-82700	U

PROCARE
ELASTIC WRIST BRACE



INDICATIONS

Mild sprains or strains, carpal tunnel syndrome, post cast removal.

CHARACTERISTICS

- Soft, durable cotton/elastic design provides compressive support
- Pre-formed aluminium stay allows for an anatomically correct fit
- Length 17.5cm

ORDERING INFORMATION

PART NUMBER			WRIST CIRC.	SIZE
LEFT	RIGHT	QUANTITY		
79-87082	79-87072	Pack of 1	<14cm	XS
79-87083	79-87073	Pack of 1	14cm - 16.5cm	S
79-87085	79-87075	Pack of 1	16.5cm - 19cm	M
79-87087	79-87077	Pack of 1	19cm - 21.5cm	L
79-87088	79-87078	Pack of 1	>21.5cm	XL

PART NUMBER			WRIST CIRC.	SIZE
LEFT	RIGHT	QUANTITY		
24-87083	24-87073	Pack of 25	14cm - 16.5cm	S
24-87085	24-87075	Pack of 25	16.5cm - 19cm	M
24-87087	24-87077	Pack of 25	19cm - 21.5cm	L



DONJOY
WRIST WRAPS



INDICATIONS

Limitation of range in one direction only, such as extension, flexion, ulnar or radial deviation. Ideal for sports e.g. swimming and diving.

CHARACTERISTICS

- Woven nylon outer shell for durability and neoprene inner lining for patient comfort
- Interchangeable foam inserts for variable levels of support
- Opposing straps limit rotation of the brace

ORDERING INFORMATION

PART NUMBER	WRIST CIRC.	SIZE
11-0386-2-06060	12.5cm - 15cm	S
11-0386-3-06060	15cm - 17.5cm	M
11-0386-4-06060	17.5cm - 20cm	L



PROCARE
CURVED ALUMINIUM
FINGER SPLINTS



INDICATIONS

Transverse fractures of the phalanges, dislocations, phalangeal and metacarpal fractures.

CHARACTERISTICS

- Lightweight aluminium design is curved to conform to the finger in a semi-flexed position
- Padded for comfort
- Sold as a pack of 12 splints

ORDERING INFORMATION

PART NUMBER	LENGTH	WIDTH	SIZE
79-71923	3.8cm	2.2cm	S
79-71925	7.6cm	2.2cm	M
79-71924	10.2cm	2.2cm	M (Long)
79-71926	14cm	2.2cm	M (X-Long)
79-71927	15.2cm	2.2cm	L
79-71928	22.9cm	2.2cm	XL

PROCARE
STAXX PLASTIC
FINGER SPLINTS



INDICATIONS

Repair of ruptured tendons, fingertip and nail bed injuries and mallet finger.

CHARACTERISTICS

- Lightweight plastic splints
- Pack of padded finger straps available separately for securing the splints
- Starter set of clear finger splints available (various sizes)

ORDERING INFORMATION

PART NUMBER		LENGTH	DIAMETER	QUANTITY	SIZE
CLEAR	SKIN				
79-72251	79-72241	5.3cm	1.5cm	Pack of 12	1
79-72252	79-72242	5.2cm	1.7cm	Pack of 12	2
79-72253	79-72243	5.8cm	1.9cm	Pack of 12	3
79-72254	79-72244	6.3cm	1.9cm	Pack of 12	4
79-72255	79-72245	5.9cm	1.9cm	Pack of 12	5
79-72256	79-72246	6.1cm	2.3cm	Pack of 12	5.5
79-72257	79-72247	6.4cm	2.3cm	Pack of 12	6
79-72258	79-72248	6.3cm	2.5cm	Pack of 12	7

PART NUMBER	DESCRIPTION	QUANTITY
79-72240	Finger Straps	Pack of 25
79-72260	Clear Splint Starter Set	Pack of 30

HOT/COLD THERAPY



HOT/COLD THERAPY

DONJOY

DURA*SOFT™ SHOULDER WRAP



INDICATIONS

Decrease pain, oedema and secondary hypoxic tissue injury.

CHARACTERISTICS

- Nylon wraps and elastic straps provide a universal fit and effective compression
- Supplied with two ice inserts to provide extended cold therapy
- May be used repeatedly to provide years of long term rehab use
- Dimensions 30.5cm x 33cm
- Fits left or right shoulder
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	ITEM	SIZE
11-0915-0-02000	Wrap	U
11-0917-0-00000	Additional Ice Inserts	Pack of 12

DONJOY

DURA*SOFT™ KNEE WRAP



INDICATIONS

Decrease pain, oedema and secondary hypoxic tissue injury.

CHARACTERISTICS

- Nylon wrap and elastic straps provide a universal fit and effective compression
- Supplied with two ice inserts to provide extended cold therapy
- May be used repeatedly to provide years of long term rehab use
- Size 28cm x 30.5cm
- Fits left or right knee

ORDERING INFORMATION

PART NUMBER	ITEM	SIZE
11-0913-0-02000	Wrap	U
11-0917-0-00000	Additional Ice Inserts	Pack of 12

AIRCAST

CRYO/CUFF™ IC COOLER



INDICATIONS

Integrated cold and compression clinically proven to help reduce post-operative swelling and regain range of motion. Ideal for post-operative recovery, trauma, athletic training rooms and home use.

CHARACTERISTICS

- Provides both gravity and motorised cold and compression therapy for flexibility of care
- Integrated pneumatic pump within the cooler lid provides automated cold and compression therapy, with simplicity of design and ease of use operation
- Compatible with all Cryo/Cuff cuffs

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
52A	Cryo/Cuff IC Cooler
REPLACEMENT PARTS	DESCRIPTION
25-0238	Cryo/Cuff IC Lid
25-2740	Cryo/Cuff IC Power Supply
2130	Cryo/Cuff Cooler Tube
2002	Cryo/Cuff Cooler Replacement Insulation Disc

AIRCAST

CRYO/CUFF™ GRAVITY COOLER



INDICATIONS

Ideal for post-operative recovery, trauma, athletic training rooms and home use. An integrated approach to cold therapy, combining cold and compression therapy.

CHARACTERISTICS

- Cold water delivered using a simple, gravity fed system
- Cooler holds water and ice needed for 6 to 8 hours of cryotherapy
- Compatible with all Cryo/Cuff cuffs and ArcticFlow wraps

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
2125	Cryo/Cuff Gravity Cooler
REPLACEMENT PARTS	DESCRIPTION
2130	Cryo/Cuff Cooler Tube
2002	Cryo/Cuff Cooler Replacement Insulation Disc

AIRCAST

CRYO/CUFF™ WRAPS & CUFFS

INDICATIONS

The anatomically designed cuffs provide the therapeutic benefits of controlled cold compression to minimise haemarthrosis and swelling and reduce pain.

CHARACTERISTICS

- Cuffs are anatomically designed to cover all major joints of the body
- Cuffs provide maximum coverage for optimum cryotherapy treatment
- Soothing cold covers the affected area reducing swelling and providing pain relief
- Cuffs are simple to apply and remove and may be detached from the cooler at any time for mobility without interruption of treatment



Foot



Ankle



Thigh



Elbow



Shoulder



Back/Hip/Rib



Knee



Hand



Calf



Paediatric Ankle

PART NO.	DESCRIPTION	SIZE		
10C01	Foot CC	M	Foot Circ.	23-33 cm
10B01	Foot CC	L		25-43 cm
10A01	Ankle CC	U	Universal	-
11C01	Knee CC	S	Circ. 15 cm	25-50 cm
11A01	Knee CC	M	(6") above patella	45-60 cm
11B01	Knee CC	L		50-80 cm
13C01	Calf CC	U	Calf Circ.	36-51 cm
13A01	Thigh CC	L	Thigh Circ.	43-58 cm
13B01	Thigh CC	XL		51-69 cm

PART NO.	DESCRIPTION	SIZE		
16A01	Hand/Wrist CC	U	-	-
15A01	Elbow CC	U	-	-
12A01	Shoulder CC	U	Torso Circ.	81-122 cm
12B01	Small Shoulder CC	S		< 86 cm
12AXL01	Shoulder CC w/ XL Strap	XL	-	107-137 cm
1220XL	Shoulder Strap	XL	-	-
14A01	Back/Hip/Rib CC	U	-	-
10P01	Paediatric Ankle CC	Ped	-	-

ICEMAN® CLEAR3™



INDICATIONS

The IceMan CLEAR³ cold therapy unit helps reduce pain and swelling, speeding up rehabilitation. It provides extended cold therapy for a variety of indications and protocols. It utilises DonJoy's patented semi-closed loop recirculation system, which maintains more consistent and accurate temperatures than other cold therapy units, in a preset configuration.

CHARACTERISTICS

- Sure fit, stackable, locking lid with watertight seal
- Low profile design, small cooler for easy portability and storage
- Transparent cooler design for easy refill indication
- Preset inlet temperature for simple and convenient, plug and play use
- Practical, integrated hose hanger for ease of use
- Compatible with a variety of pads

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
11-1639	IceMan CLEAR ³



DONJOY
ICEMAN® MCGUIRE
KNEE WRAP



INDICATIONS

Reduce pain, swelling and facilitate rehabilitation.

CHARACTERISTICS

- Anatomical design provides effective cold therapy to the knee
- Wrap dimensions 27cm x 28.5cm
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	SIZE
11-1268-9-00000	U



DONJOY
ICEMAN® FOREFOOT
WRAP



INDICATIONS

Bunionectomies and procedures of the forefoot.

CHARACTERISTICS

- Anatomical design suitable for treatment of the forefoot
- Includes both cold pad and wrap
- Fits left or right foot
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	SIZE
11-1521	U



DONJOY
**ICEMAN® HIP
 WRAP-ON PAD**



INDICATIONS

Procedures of the hip.

CHARACTERISTICS

- Anatomical design suitable for treatment of the hip
- Includes both cold pad and wrap

ORDERING INFORMATION

PART NUMBER		SIZE
LEFT	RIGHT	
11-9112	11-9111	U



DONJOY
**ICEMAN®
 WRAP-ON COLD PAD**



INDICATIONS

Reduce pain, swelling and facilitate rehabilitation.

CHARACTERISTICS

- Universal shape, suitable for use on many parts of the body
- Wrap dimensions 27cm x 28.5cm
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	SIZE
11-0679-9-00000	U



DONJOY
**ICEMAN® ANKLE
 CUFF**



INDICATIONS

Procedures of the ankle, forefoot and hind foot.

CHARACTERISTICS

- Anatomical design suitable for treatment of the ankle
- Includes both cold pad and wrap
- Fits left or right ankle

ORDERING INFORMATION

PART NUMBER	SIZE
11-1522	U



DONJOY
**ICEMAN® SHOULDER
 WRAP-ON PAD**



INDICATIONS

Reduce pain, swelling and facilitate rehabilitation of the shoulder.

CHARACTERISTICS

- Anatomical design provides effective cold therapy to the shoulder
- Includes both cold pad and wrap
- Fits left or right shoulder
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	SIZE
11-1334-9-00000	U



INDICATIONS

These durable stainless steel water tanks house the Hydrocollator Moist Heat HotPacs. The units are thermostatically controlled to ensure a constant temperature of between 71°C and 74°C, an ideal therapeutic temperature for HotPacs. All models include a thermal safety cut-out. Hydrocollator remains the benchmark heating unit against which all other units are judged.

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	TANK CAPACITY	HEATING TIME	HOTPACS WITH UNIT
2107	E-1 Stationary	15L	2 Hours	4 x Standard
2806	E-2 Stationary	43L	4 Hours	3 x Standard, 2 x Oversize, 1 x Neck
2407	M-2 Mobile	69L	8 Hours	12 x Standard
2307	SS-2 Mobile	49 L	4 Hours	8 x Standard

ACCESSORY

PART NUMBER	DESCRIPTION	DIMENSION (W/L/H)	FEATURE
4002	Trolley MB	48cm x 43cm x 51cm	For E-2
4004	Trolley MB-T	48cm x 43cm x 76cm	For E-1/E-2
4018	Trolley MB-T	48cm x 43cm x 76cm	For E-1/E-2 with Drawer
22100	Mini Lotion Warmer	-	Fits front of E-1/E-2



CHATTANOOGA HOTPAC TERRY COVERS

INDICATIONS

Terry covers are suitable for use with Moist Heat HotPacs, these absorb excess water and create a layer between the HotPac and the patient.



CHARACTERISTICS

- All-Terry Covers feature hook and loop closures to prevent unwrapping or slipping
- Foam-Filled Terry Covers have a polyurethane foam layer between two layers of terry
- 1 year warranty

ORDERING INFORMATION

PART NUMBER	MATERIAL	FITS HOTPACS	DIMENSIONS
1118	All-Terry Cover	Standard Size	48cm x 69cm
1126	All-Terry Cover	Oversize	61cm x 91cm
1104	All-Terry Cover	Neck Contour	64cm x 46cm

PART NUMBER	MATERIAL	FITS HOTPACS	DIMENSIONS
1102	Foam Filled Cover	Standard Size	48cm x 69cm
1124	Foam Filled Cover	Oversize	61cm x 91cm
1120	Foam Filled Cover	Neck Contour	64cm x 41cm



CHATTANOOGA DIAL THERMOMETER

INDICATIONS

Easy monitoring of water temperature within Hydrocollator heating units. Can also be used with chilling units.

CHARACTERISTICS

- Stainless steel thermometer with Celsius and Fahrenheit readings
- Temperature ranges: 50°C and 100°C
- Clip for attachment to Hydrocollator wire racks
- 1 year warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
4228	Dial Thermometer

HYDROCOLLATOR HOTPACS



INDICATIONS

Hydrocollator Moist Heat HotPacs have been the first choice for physiotherapists, physicians, chiropractors, athletic trainers and kinesiologists for over 60 years. Chattanooga is the inventor of the original Hydrocollator HotPac and is proud to still be producing this all natural product today. The benefits of moist heat therapy over dry heat are well established for pain relief due to arthritis, myofibrosis, back pain, muscle strains and spasms.

CHARACTERISTICS

- Seven HotPac sizes designed to fit every body contour
- Provides up to 30 minutes of deep, soothing, moist heat
- Coloured tabs on most HotPacs for easy rotation in busy clinics
- Cost effective, provides hundreds of treatments from one HotPac
- All natural pack containing bentonite, a natural hydrocolloid filling, wrapped in a woven cotton blend
- 1 year warranty

ORDERING INFORMATION

PART NUMBER	SIZE
1008	23cm x 39cm



LONG

PART NUMBER	SIZE
1010	25cm x 61cm



HALF SIZE

PART NUMBER	SIZE
1014	13cm x 30cm



OVERSIZE

PART NUMBER	SIZE
1004	38cm x 61cm



STANDARD

PART NUMBER	SIZE
1006	25cm x 30cm



NECK

PART NUMBER	SIZE
1002	61cm



KNEE/SHOULDER

PART NUMBER	SIZE
1012	25cm x 50cm

CHATTANOOGA
COLPAC®

INDICATIONS

Cold therapy, when applied directly to the affected area, is used to control pain due to inflammation, swelling, multiple sclerosis, muscle strains, sprains and spasms.

CHARACTERISTICS

- Available in seven different sizes of vinyl ColPac which conform to body area or shape
- Effective relief of pain
- Filled with a non-toxic silica gel that remains pliable throughout treatment
- Delivers up to 30 minutes of soothing relief from acute pain, swelling and fever
- Latex free
- 1 year warranty



OVERSIZE

PART NUMBER	SIZE
1512	28cm x 53cm



STANDARD

PART NUMBER	SIZE
1500	28cm x 36cm



CHATTANOOGA
FLEXIPAC

INDICATIONS

Hot or cold therapy, when applied directly to the affected area, is used to control pain due to inflammation, swelling, multiple sclerosis, muscle strains, sprains and spasms.

CHARACTERISTICS

- Flexible gel pack conforms to any body contour at any temperature
- Microwaveable to provide soothing heat
- Freezable to provide penetrating cold therapy
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	QUANTITY	SIZE
4026	Pack of 48	3cm x 15cm
4020	Pack of 24	13cm x 25cm
4029	Pack of 12	20cm x 35cm



HALF SIZE

PART NUMBER	SIZE
1506	19cm x 28cm



QUARTER SIZE

PART NUMBER	SIZE
1504	14cm x 19cm



NECK CONTOUR

PART NUMBER	SIZE
1508	58cm



LONG

PART NUMBER	SIZE
1502	8cm x 28cm



CLOSED EYELIDS

PART NUMBER	SIZE
1510	Eye



CHATTANOOGA
COLPAC® CHILLING UNIT

CHARACTERISTICS

Large capacity and faster chilling time mean your patients will never have to wait for a Colpac. The Colpac chilling unit helps you keep pace throughout a busy day. No plumbing required, simply plug into a standard electrical outlet. Colpac chilling unit includes:

- Closed-cell foam insulation for efficiency
- Stainless steel, welded construction for long life
- Heavy-duty compressor for more efficient cooling
- Specially designed interior chilling coils cool packs faster
- 8cm swivel-type rubber casters for silent, friction-free movement of mobile units
- Drain valve for easy cleaning and defrosting
- Supplied with six Standard ColPac and six Half Size ColPac
- 1 year warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
3114	ColPac Chilling Unit



CMF

CMF BONE GROWTH STIMULATION

RECOMMENDED SIZING

FRACTURE SITE	DUAL COIL Part No. 01-203-0001 L9" x W61/2 x h7" 1.85lbs	SIZE 1 Part No. 01-211-0001 L4" x W41/4 x h11/4" 0.9lbs	SIZE 2 Part No. 01-211-0002 L9" x W81/2 x h41/2" 1.9lbs	SIZE 3 Part No. 01-211-0003 L101/2" x W111/2 x h53/4" 2.45lbs	SIZE 4 Part No. 01-211-0004 L13" x W14 x h9" 4.95lbs	SPINE Part No. 01-207-0007 L101/2" x W9 x h1" 1.85lbs
CLAVICLE	Proximal		○			
	Mid		○			
	Distal		○			
HUMERUS	Proximal			○		
	Mid			○		
	Distal	○		○		
RADIUS/ULNA	Proximal	○		○		
	Mid	○		○		
	Distal	○				
HAND/WRIST	Phalanges	○	○			
	Metacarpals	○	○			
FEMUR	Proximal			○	○	○
	Mid			○	○	○
	Distal			○	○	
TIBIA/FIBULA	Proximal			○	○	
	Mid	○		○		
	Distal	○		○		
FOOT/ANKLE	Phalanges	○	○			
	Metacarpals	○	○			
	Ankle	○		○		
SCAPHOID	○	○				
LUMBAR/SPINE						○

CMF

CMF SPINALOGIC®



INDICATIONS

The CMF SpinaLogic is a portable, battery powered, micro controlled, non-invasive bone growth stimulator, indicated as an adjunct electromagnetic treatment to primary lumbar spinal fusion surgery for one or two levels. Anterior or posterior approach in fusion surgery.

CHARACTERISTICS

- Treatment time of 30 minutes per day
- Combined magnetic field technology at a frequency of 76.6Hz promotes bone healing
- Lightweight, portable battery powered medical device with simple one-button operation
- Built-in patient compliance monitor

ORDERING INFORMATION

PART NUMBER	SIZE
01-207-0007	U

CMF

CMF OL1000™ DUAL COIL



INDICATIONS

Non-invasive treatment of an established non-union fracture acquired secondary to trauma, excluding all vertebrae and flat bones.

CHARACTERISTICS

- Treatment time of 30 minutes per day
- Combined magnetic field technology at a frequency of 76.6Hz promotes bone healing
- Lightweight, portable battery powered medical device with simple one-button operation
- Built-in patient compliance monitor
- Can be used with internal or external fixation or over a cast

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
01-203-0001	OL1000 Dual Coil

CMF

CMF OL1000™ SINGLE COIL



INDICATIONS

CMF OL1000 Bone Growth Stimulators are portable, battery-powered medical devices indicated for use in the non-invasive treatment of an established non-union fracture acquired secondary to trauma, excluding all vertebrae and flat bones. This treatment has been shown in pre-clinical studies to help the body's own healing process begin working. Clinical studies have shown an increased chance of healing of 60.7% in patients with non-unions that averaged 29.3 months from injury. Registry data on over 2300 patients reports a heal rate of 75.1%.

CHARACTERISTICS

- Lightweight and comfortable
- Easy to use & non-invasive
- Requires simple, one-button operation
- Device is worn for only 30 minutes per day
- Can be used with internal or external fixation or over a cast
- Maximum 270 day treatment cycle from first use

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	DIAMETER
01-211-0002	OL1000 Single Coil, Small	18cm
01-211-0003	OL1000 Single Coil, Medium	25cm
01-211-0004	OL1000 Single Coil, Large	30cm

CMF

CMF OL1000™ SINGLE COIL



INDICATIONS

Non-invasive treatment of a delayed or established non-union fracture acquired secondary to trauma, excluding all vertebrae and flat bones.

CHARACTERISTICS

- Lightweight, portable battery powered medical device with simple one-button operation
- Lightweight model, specifically designed for scaphoid and clavicle fractures
- Treatment time of 30 minutes per day
- Combined magnetic field technology at a frequency of 76.6Hz promotes bone healing
- Built-in patient compliance monitor
- Can be used with internal or external fixation or over a cast

ORDERING INFORMATION

PART NUMBER	SIZE
01-211-0001	Size 1 (Universal)



VASCULAR

AIRCAST

VENAFLOW® ELITE



INDICATIONS

VenaFlow is the only VTE prophylaxis compression device that is proven to mimic ambulation, thus reducing VTE rates by 50% vs. other devices. This is the best technology to protect patients and healthcare practices from the rising risk of the VTE crisis

- Prophylaxis for deep vein thrombosis
- VenaFlow Elite System combines normal inflation and graduated sequential compression to accelerate venous velocity and enhance fibrinolysis
- For the superior emptying of veins, VenaFlow Elite also incorporates asymmetric compression

CHARACTERISTICS

- Low profile, lightweight design
- Compliance monitor/alarm
- Telescoping bed hanger for custom adjustability
- Automatic Cuff Detection - system identifies and adjusts pressure based on cuffs attached
- Single Leg Operation Button - one push of a button selects single leg operation
- Preset pressures and alarms - no adjustments necessary
- One pump for calf, thigh and foot cuffs provides for ease of use and minimises inventory
- Breathable and comfortable calf, foot and thigh cuffs available

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
30BI	VenaFlow Elite
30BI-B	VenaFlow Elite, with Battery
30BI-S	VenaFlow Elite, S Mode
30BI-SB	VenaFlow Elite, S Mode with Battery



AIRCAST VENAFLOW® ELITE THIGH CUFF

PART NUMBER	QUANTITY	SIZE
3045	Pair	U



AIRCAST VENAFLOW® ELITE FOOT CUFF

PART NUMBER	QUANTITY	SIZE
3046	Pair	U
24-3046	20 Pairs	U



AIRCAST VENAFLOW® ELITE CALF CUFF

PART NUMBER	QUANTITY	MAX CALF CIRC.	SIZE
3040	Pair	48cm	U
24-3040	20 Pairs	48cm	U
3042	Pair	56cm	XL
24-3042	20 Pairs	56cm	XL
3043	Pair	76cm	XXL (Bariatric)



AIRCAST VENAFLOW® TUBING

PART NUMBER	TUBE LENGTH
3008	1.7m (5.5')
3008XL	2.6m (8.5')
3008XXL	3.2m (10.5')
3008XXXL	3.8m (12.5')

AIRCAST
VENAPURE®



VenaPure anti-embolism stockings have been designed with Venous Thromboembolism (VTE) prevention in mind. Enovis has been able to use the years of clinical experience gained in the arena of VTE prophylaxis to combine all of these positive attributes into VenaPure.

INDICATIONS

For use in the prophylaxis of deep vein thrombosis and pulmonary embolism in hospitalised patients and for use in other conditions as directed by the clinician.



CHARACTERISTICS

- Graduated compression according to Sigel profile (18mmHg at the ankle, decreasing to 8mmHg at the thigh) to minimise venous stasis and venous dilation
- Inspection toe opening on top of the foot to minimise risk of pressure point and ensure that skin integrity is maintained
- Popliteal break is easily identified by stitch change to ensure that pressure is relieved around the popliteal vein
- Pressure relieving panel at the inner thigh to ensure reduced compression of the femoral vein
- Continuous silicone band at upper thigh provides excellent fixation to the skin, without compromising compression
- Inlay circumferential knit maintains correct positioning of stocking every time
- Colour coded toe and top ensure that the correct size is selected according to patient measurement
- Tape measure included in the dispenser to make measuring of the patient easier
- Available in two styles with lengths and sizes to fit all patients

ORDERING INFORMATION KNEE LENGTH (PACK OF 12 PAIRS)

CALF CIRC.	LEG LENGTH	CODE	SIZE
<30.5cm	<41cm	400R-2	S
	>41cm	400L-2	
30.5cm - 38cm	<43cm	400R-3	M
	>43cm	400L-3	
38cm - 44.5cm	<46cm	400R-4	L
	>46cm	400L-4	
44.5cm - 51cm	<46cm	400R-5	XL
	>46cm	400L-5	
51cm - 58.4cm	<46cm	400R-6	XXL
	>46cm	400L-6	
58.5cm - 66cm	<46cm	400R-7	XXXL
	>46cm	400L-7	

ORDERING INFORMATION THIGH LENGTH (PACK OF 6 PAIRS)

THIGH CIRC.	CALF CIRC.	LEG	CODE	SIZE
<30.5cm		<74cm	401S-2	S
		74cm - 84cm	401R-2	
		>84cm	401L-2	
<63.5cm	30.5cm - 38cm	<74cm	401S-3	M
		74cm - 84cm	401R-3	
		>84cm	401L-3	
	38cm - 44.5cm	<74cm	401S-4	L
		74cm - 84cm	401R-4	
		>84cm	401L-4	
63.5cm - 81.5cm	44.5cm - 55cm	<74cm	401S-5	XL
		74cm - 84cm	401R-5	
		>84cm	401L-5	
81cm - 91cm	55cm - 66cm	<74cm	401S-6	XXL
		74cm - 84cm	401R-6	
		>84cm	401L-6	



CLINICAL THERAPY






CLINICAL THERAPY

Chattanooga is the world's largest manufacturer of rehabilitation equipment for treating musculoskeletal, neurological and soft tissue disorders. For over 75 years, Chattanooga has set a benchmark for leadership, reliability and excellence. Our products contribute to better treatment outcomes in hospitals, clinics and home settings worldwide. We lead the physical therapy industry by example and, through continuous innovation, provide real world solutions for clinicians and their patients.

As the world's leading brand in rehabilitation, Chattanooga is sold in over 80 countries worldwide. Our Chattanooga customers are supported for order placement, training, product delivery and after sales service by a dedicated team of highly skilled Enovis employees.





INDICATIONS

Shockwaves are used to treat a wide range of pathologies in areas such as bone healing, tendinopathies, muscular trigger points/myofascial pain, dermatology, wound healing, spasticity and impotence (erectile dysfunction). Focus shockwave has high energy (around one hundred times more powerful than a Radial Pressure Wave) making it suitable for a range of conditions in urology, orthopaedic indications, cardiology and neurology. Focus Shockwave’s power (shockwaves can penetrate approximately 12cm with pinpoint focus) allows it to disintegrate calcifications, stones, and growths as well as providing the standard benefits of shockwave therapy; activation of the healing cascade, enhanced blood flow, and tissue regeneration.

CHARACTERISTICS

- 12.5cm penetration with pinpoint focus
- Enhanced energy: 0.01 - 0.55mJ/mm2
- Broad frequency range: 1 - 8 Hz
- 3 standoff variations for tailored treatment
- LCD touch screen
- 2 year warranty, or up to 1 million pulses on the handpiece


ORDERING INFORMATION

PART NUMBER	DESCRIPTION
21090	Intelect Focus Shockwave

PART NUMBER	OPTIONAL ACCESSORIES
4650	Focus Shockwave Cart
19100	30mm Stand-off I
19200	15mm Stand-off II
19300	Transparent Closing Ring for Stand-off I and II



STAND-OFFS AND PENETRATION DEPTHS OF THE F-SW HANDPIECE



TARGETED PAIN RELIEF

Shockwaves are used to treat a wide range of pathologies in areas such as bone healing, tendinopathies, muscular trigger points/myofascial pain, dermatology, wound healing, spasticity and impotence (erectile dysfunction). It can be easier to find deeper, painful, diffuse trigger points with shockwave therapy than with palpation.



EPICONDYLITIS



CALCIFIC TENDINITIS



PATELLAR TENDINITIS



TRAPEZIUS MUSCLE



TIBIAL STRESS SYNDROME



HEEL SPUR



ACHILLODYNIA



PLANTAR FASCIITIS

INDICATIONS

The Intellect® RPW 2 allows for the treatment of indications with radial, pneumatically-generated, low-energy acoustic waves, or ‘pressure pulses’. The applied pressure pulse propagates radially within the tissue, generating a therapeutic effect. The V-ACTOR® HF optional hand piece allows mechanical vibration, for vibration massage.

TREATMENT GUIDELINES

- Localize the painful areas you plan to treat
- Select most appropriate transmitters based on treatment goal and targeted tissue
- Apply ultra sound water based gel on the skin overlying the treatment area
- Always start with lower level of intensity
- Perform a post-treatment test to assess improved function and decrease in pain

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
2174KIT	Intellect RPW 2 230V unit (set includes the below items)
18080	FALCON® handpiece, including 2x projectile; 2x guiding tube; 1 cleaning brush, 1x wrench; 1x allen key
13-28660*	Quick start guide booklet (in multiple languages)
22652	Transmitters tray
22651	Gel bottle support cup
4248	1x Transmission gel, 250 ml bottle
28178	Transmitter Ro40
29724	D20-S Transmitter
82-0274	Intellect RPW 2 pouch
15-1140	USB Key





SUGGESTED PROTOCOLS LIBRARY

Available at a finger click; treatment pictures, anatomical and pathological library, treatments guidelines.

LARGE GRAPHIC COLOR TOUCH SCREEN

High resolution 7" capacitive touch screen for greater customer experience.

NEW FALCON® HANDPIECE

allowing lower Bar pressures for an increased range of indications and applicators.



COMFORT MODE

Unique automatic energy ramp up designed to provide better patient comfort.

V-ACTOR® HF*

Vibration frequency up to 50Hz for Vibration therapy and muscle relaxation.

LARGE CHOICE OF TRANSMITTERS*

17 different transmitters available for tendinopathies, trigger points, spine and fascia treatments.

EASY FIT FOR EVERY TREATMENT ROOM

Slim design 41x42x110cm (LxWxH), not a lot of space needed for this new RPW generation.

ERGONOMIC DESIGN

What you need close in hand; handpiece, transmitter tray, gel bottle, cables holders, soft pouch.

EXTENDED PRESSURE RANGE

The Falcon handpiece provides a greater pressure range from 0.3 to 5 bar, making it more comfortable for the patient and allowing Fascia and Spine treatments. NOTE: PERI-V-ACTOR® and SPINE-V-ACTOR® can only be used under 3 Bar pressure.



ORDERING INFORMATION

PART NUMBER	DESCRIPTION
18080	FALCON handpiece, including 2x projectile; 2x guiding tube; 1 cleaning brush; 1x wrench; 1x allen key
29104	Revision kit short including 2x projectiles; 2x guiding tubes short; 2x o-ring set; 1x o-ring guide
29801	A6 transmitter black
29802	T10 transmitter black
29729	C15 transmitter black
29728	DI15 transmitter black
29726	F15 transmitter black
29724	D20-S transmitter black
29725	D35-S transmitter black
28736	D20-T transmitter black

PART NUMBER	DESCRIPTION
28178	Ro40 transmitter red
18638	R15 transmitter black
29539	Atlas RPW transmitter
19365	V-ACTOR® set: V-ACTOR®, V25 and V40
29742	V-ACTOR V10 transmitter
28740	V-ACTOR V25 transmitter
28741	V-ACTOR V40 transmitter
31800	Transmitter box including 9 transmitters
28945	PERI-ACTOR® fascia transmitter set
28946	SPINE-ACTOR® transmitter set
82-0274	Intellect RPW 2 pouch
22651	Gel bottle support cup
22652	Transmitter tray
4248	24x Transmission gel, 250 ml
14820	Power cord EU

CHATTANOOGA

MOBILE 2 RPW SHOCKWAVE



INDICATIONS

Offers a broader range of shock frequency and greater power (up to 5 Bar of pressure), making it a more flexible tool for helping patients recover. The versatility of the Mobile 2 RPW makes it suitable for a range of environments and applications: orthopaedics, physiotherapy, aesthetic therapy, chiropractors, rheumatology, sports clubs, veterinary, physical medicine, acupuncture and podiatrists.

TREATMENT GUIDELINES

- All the features of the original Mobile RPW but with improved functionality
- Single and Continuous shocks: Frequency D-ACTOR: 1 - 21 Hz, Frequency V-ACTOR: 1 - 31 Hz
- 1-5 Bars of pressure, infinitely variable with precise adjustments of 0.1 increments
- 2 year warranty, 1 year on accessories

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
2910	Mobile 2 RPW
31800	Set of 9 transmitters in carry case



ENQUIRE NOW

CHATTANOOGA

RPW SHOCKWAVE ACCESSORIES



MOBILE RPW CARRY BAG

PART NUMBER	DESCRIPTION
28745	Mobile RPW Carry Bag




REVISION KIT

PART NUMBER	DESCRIPTION
28701	Revision Kit (2x Guide Tube, Projectile, Sealing Set)



V-ACTOR KIT

PART NUMBER	DESCRIPTION
28730	V-ACTOR Kit



D-ACTOR KIT

PART NUMBER	DESCRIPTION
28699	D-ACTOR Kit (no transmitters)

CHATTANOOGA

INTELECT® RPW LITE

INDICATIONS

For today's busy physical therapist, the need to treat patients efficiently is crucial. The Intellect® RPW Lite is a compact and user-friendly device that delivers effective Radial Pressure Wave (RPW) therapy with the minimum of fuss. Its range of essential features have been designed to remove the complexity from RPW treatment, leaving therapists free to concentrate on their patients.

TREATMENT GUIDELINES

- Localize the painful areas you plan to treat
- Select the most appropriate transmitters based on treatment goal and targeted tissue
- Apply water-based ultrasound gel to the skin overlying the treatment area
- Always start with a lower level of intensity
- Perform a post-treatment test for improved function and decrease in pain



PART NUMBER	DESCRIPTION
30680-230V	Intellect® RPW Lite 230V set including standard accessories
30885-1	Intellect RPW Lite SPARROW™ handpiece kit grey, including 1x projectile; 1x guiding tube; 1x wrench; 1x allen key
30873	Transmitter D20-S Black RPW Lite, 20mm



ENQUIRE NOW

PART NUMBER	DESCRIPTION
00032.012	EU Power cord
30891	Operating manual and quick start guide
22601	Transmission gel, 250 ml
28739	Guiding tube cleaning brush
30890	Tube set (includes 3 tubes)



PORTABILITY

Compact and easy to handle.

LARGE LCD TOUCH SCREEN

For quick and easy control of treatment parameters.

PULSE COUNTER

INTENSITY LEVEL

6 preset energy levels reduce set-up time whilst providing a useful range of settings for a variety of possible RPW indications.

EASY CONNECT

SPARROW™ handpiece technology with quick-release connection allows easy interchange of handpieces, ideal for treatments that require multiple different transmitters.

CHOICE OF 5 TRANSMITTERS



RPW TRANSMITTERS

A FOR USE WITH INTELECT RPW 2, MOBILE 2 RPW, RPW LITE

<p>C15</p>  <ul style="list-style-type: none"> • Ø 15 mm, CERAmA-x® transmitter • Penetration depth: 0–30 mm 	<p>Ro40 BEAM</p>  <ul style="list-style-type: none"> • Ø 6 mm, AkuSt transmitter • Penetration depth: 0–20 mm • RPW Lite version is black 	<p>DI15 GOLDEN DEPTH</p>  <ul style="list-style-type: none"> • Ø 10 mm, trigger point transmitter • Penetration depth: 0–40 mm 	<p>D20-S STANDARD OSCILLATOR</p>  <ul style="list-style-type: none"> • Ø 6 mm, AkuSt transmitter • Penetration depth: 0–20 mm 	<p>D20-T GOLDEN OSCILLATOR</p>  <ul style="list-style-type: none"> • Ø 10 mm, trigger point transmitter • Penetration depth: 0–40 mm
--	--	--	---	---

B FOR USE WITH INTELECT RPW 2, MOBILE 2 RPW

<p>A6 POINT</p>  <ul style="list-style-type: none"> • Ø 6 mm, AkuSt transmitter • Penetration depth: 0–20 mm 	<p>T10 FINGER</p>  <ul style="list-style-type: none"> • Ø 10 mm, trigger point transmitter • Penetration depth: 0–40 mm 	<p>F15 WHITE SOFT</p>  <ul style="list-style-type: none"> • Ø 15 mm, CERAmA-x® transmitter • Penetration depth: 0–30 mm 	<p>D35-S LARGE OSCILLATOR</p>  <ul style="list-style-type: none"> • Ø 15 mm, CERAmA-x® transmitter • Penetration depth: 0–30 mm 	<p>V10 VIBRATION THERAPY</p>  <ul style="list-style-type: none"> • Ø 10mm • Treatment of superficial trigger points 	<p>V25 VIBRATION THERAPY</p>  <ul style="list-style-type: none"> • Ø 25mm • Lymphatic drainage • Muscle smoothing • Massage • Increase in local blood circulation • Connective tissue activation 	<p>V40 VIBRATION THERAPY</p>  <ul style="list-style-type: none"> • Ø 40mm • For large Muscle areas • Lymphatic drainage • Muscle smoothing • Massage • Increase in local blood circulation • Connective tissue activation
--	---	---	---	--	--	--

C FOR USE WITH INTELECT RPW 2

<p>ATLAS</p>  <ul style="list-style-type: none"> • Soft flexible tipped transmitter 15mm for highly sensitive regions and triggerpoints in the cervical spine • Penetration depth 0-15mm 	<p>SPINE-ACTORS®</p>  <ul style="list-style-type: none"> • Double Tip shape to allows treatments of the paravertebral region 	<p>FASCIA TREATMENTS</p>  <p>PERI-ACTOR® KNUCKLE / SCOOP / SCRAPER / SPHERE</p> <ul style="list-style-type: none"> • Unique shape especially designed to dissolve clots in the Fascia structures
--	---	---

PART NO.	DESCRIPTION
----------	-------------

A INTELECT RPW 2, MOBILE 2 RPW

29729	C15
28178	Ro40 Beam
29728	DI15 Golden Depth
29724	D20-S Oscillator
28736	D20-T Golden Oscillator

A RPW LITE

30875	C15
30842	Ro40 Beam
30874	DI15 Golden Depth
30873	D20-S Oscillator
30876	D20-T Golden Oscillator

PART NO.	DESCRIPTION
----------	-------------

B INTELECT RPW 2, MOBILE 2 RPW

29801	A6 Point
29802	T10 Finger
29726	F15 White Soft
29725	D35-S Large Oscillator
29742	V-ACTOR V10
28740	V-ACTOR V25
28741	V-ACTOR V40

PART NO.	DESCRIPTION
----------	-------------

C FOR USE WITH INTELECT RPW 2

29539	Atlas
28946	SPINE-ACTOR® set
28945	PERI-ACTOR® set

CHATTANOOGA**RUBBER ELECTRODES****INDICATIONS**

Silicone rubber electrodes for use with Stim devices.

CHARACTERISTICS

- 2mm pin connector
- Sold as a pack of 1
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
79966	4cm x 6cm
79967	6cm x 8cm
79968	8cm x 12cm

**CHATTANOOGA****ELECTRODE SPONGES****INDICATIONS**

Electrode sponges for use with Silicone Rubber Electrodes.

CHARACTERISTICS

- Sponge pockets available in three sizes to match electrode sizing
- Sold as a pack of 1
- 3 month warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
79969	4cm x 6cm
79970	6cm x 8cm
79971	8cm x 12cm

**CHATTANOOGA****NYLATEX INDIVIDUAL****INDICATIONS**

Ideal for holding sponge pockets, electrodes and cold therapy packs in place.

CHARACTERISTICS

- Elastic, reusable design
- Wraps sold in packs of 3
- 3 month warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
1202	6cm x 46cm
1204	6cm x 91cm
1206	6cm x 122cm
1208	10cm x 46cm
1210	10cm x 91cm
1212	10cm x 122cm
1218	15cm x 122cm
1264	6cm x 61cm

CHATTANOOGA

INTELECT® MOBILE 2

INDICATIONS

The Intellect® Mobile 2 is the next generation in physical medicine modalities, its intuitive design is clever and its features, usability and simplicity make Intellect Mobile 2 the choice for today's modern practice.

The Intellect Mobile 2 is simple to use and set up thanks to the built in suggested protocol settings developed around the current clinical practices. The Intellect Mobile 2 has an easily accessible high resolution anatomical library which allows you to explain visually to the patient their condition and the treatment pathway.

TREATMENT GUIDELINES

STIM

TENS

- Symptomatic relief of chronic pain
- Management of post-operative PAIN

NMES


- Muscle re-education
- Increasing local blood supply
- Relaxation of muscle spasms
- Maintaining/increasing range of motion

ULTRASOUND

- Relief of pain from muscle spasm
- Relief of pain from joint contracture
- Relief of pain associated with ligament sprains, tendinitis and muscle spasm

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
15-1205	Intellect® Mobile 2 Combo Set (Combo set includes the below items)
79967	Carbon electrodes
42198	Gel electrodes
70010	Stim lead wires
6522055	Chattanooga® straps
15-0142	5cm2 ultrasound applicator
4248	Transmission gel bottle
15-0146	Power cord
15-1140	USB Key





LARGE GRAPHIC COLOR TOUCH SCREEN

High resolution 7" capacitive touch screen for access to intuitive user interface and extensive anatomical library. Suggestions for the placement of electrodes and areas for ultrasound applicator placement are in easy to use and find guidelines. The touch screen also gives real time information as to the progress of the individual treatment and a summary after treatment ends.



NEW ULTRASOUND APPLICATOR

Ergonomic design fit for the hand. Improved coupling detection and monitoring.

MOBILE

Lightweight device with ergonomically designed handle and optional integral battery, eliminating the need to carry mains wires on home visits. Optional custom carry case available for transporting the device safely and securely.

ERGONOMIC DESIGN

New contemporary modern design to fit in any therapy office, but versatile enough with the integral option of the battery to make the device simple and convenient to use in the patients own home.

HOME USE

Intellect® Mobile 2 was designed to be compliant with IEC 60601-1-11 and can be used for treatment in the patients own home as well as the therapist office, giving flexibility and versatility for both clinician and patient alike.

MEMORY

The Intellect Mobile 2 has built in suggested protocols based on currently available clinical evidence, the user can also write and store custom protocols which can be stored with a specific ID or treatment.



STIM ORDERING INFORMATION

PART NUMBER	DESCRIPTION
15-1203	INTELLECT® MOBILE 2 STIM SET (Stim set includes the below items)
79967	Carbon electrodes
42198	Gel electrodes
70010	Stim lead wires
6522055	Chattanooga® straps
15-0146	Power cord

ULTRASOUND ORDERING INFORMATION

PART NUMBER	DESCRIPTION
15-1201	INTELLECT® MOBILE 2 ULTRASOUND SET (US set includes the below items)
15-0142	5cm ² ultrasound applicator
4248	Transmission gel (1 bottle)
15-0146	Power cord

ACCESSORIES

PART NUMBER	DESCRIPTION
15-1136	Intellect Mobile 2 cart
15-1137	Intellect Mobile 2 cart with vacuum
70010	Stim channel 1/2 lead wire kit
70012	Stim channel 1/2 lead wire kit XL
15-1140	USB drive
79967	6x8cm carbon electrodes (4x)
15-0140	1cm ² ultrasound applicator
15-0141	2cm ² ultrasound applicator
15-0142	5cm ² ultrasound applicator
15-0143	10cm ² ultrasound applicator
4248	Transmission gel (24 bottles)
14-108	Intellect Mobile 2 battery

CHATTANOOGA

LIGHTFORCE® THERAPY LASERS




INDICATIONS

Deep Tissue Laser Therapy™ is used to treat painful muscles and joints associated with acute and chronic soft tissue injuries. It is also indicated for helping relieve pain and stiffness associated with osteoarthritis via its ability to increase localized blood flow. These factors make therapeutic laser a sensible modality choice for post-activity recovery.

POWER

Ranging from 15 all the way to 40 watts of power, LightForce® Therapy Lasers provide the power needed to reach deep tissues. The increased power directly impacts the time needed to apply a therapeutic dose of light energy, allowing clinicians to decrease the time required to treat effectively. The higher power equips clinicians to treat deeper and faster while covering more tissue area in less time.



INFLUENCE™ TECHNOLOGY

Influence™ Technology takes into account the patient's unique absorption characteristics and specific clinical condition to help provide the most effective treatment. Patient specifications such as condition, body type, and skin type are taken into account. The in-depth protocol settings provide recommended parameters on treatment heads and power levels to aid in clinician confidence and reduced learning time.

EMPOWER IQ™ DELIVERY SYSTEM

MODELS: 40W XLI, 25W XPI

The LightForce® Empower IQ™ Delivery System is changing the game in laser therapy emission. Designed to be user friendly and responsive, Empower IQ™ is the new generation of delivery technology.

FIBER

Ultra-premium workmanship and next level durability. Extra-long for user convenience.

QUICK DISCONNECT HEADS

Quick release and change-out treatment heads for better ease of use and versatile delivery options. The smart system senses appropriate applicators, treatment heads, associated power, and treatment hand speed.



SMART HANDPIECE

Built-in dosing feedback technology with embedded finger switch and indicator lights. In order to achieve effective treatments and improve dosing accuracy, the Smart Handpiece assesses the operator's speed and gives real-time visual and sensory feedback. The Smart Handpiece is ergonomically styled for operator comfort.



RECOMMENDED SPEED



SPEED TOO FAST



SPEED TOO SLOW

CHATTANOOGA

LIGHTFORCE® XLI



40W

Treat larger, deeper areas and provide faster pain relief with real time dosing feedback. Chattanooga's most powerful laser, the 40W LightForce® XLI impacts deep tissue faster for pain relief in less time. The LightForce smart hand piece and onscreen recommendations increase dosing accuracy and help clinicians treat with confidence.

SPECIFICATION

Laser Type : Class IV, Solid State
Laser Wavelength: 980/810 nm
Laser Power: 1W – 40W
Operating Modes: CW or Pulsed
Aiming Beam: 650 nm, < 4.5 mW
Dimensions: 43 cm x 28cm x 26cm
Weight: 9kg

ENQUIRE NOW

WARRANTY

Device: 24 month

PART NUMBER	DESCRIPTION
SYS-000656-000-00	LightForce XLI (40W) Therapy Laser System

CHATTANOOGA

LIGHTFORCE® XPi



25W

Smarter than your average laser. The 25W LightForce® XPi gives real time visual and haptic feedback, allowing clinicians to provide consistent, accurate treatments. The XPi's built-in intelligence tells users which treatment heads are acceptable or not recommended for different power levels, which helps to provide for safer treatments and less guesswork from clinicians.

SPECIFICATION

Laser Type : Class IV, Solid State
Laser Wavelength: 980/810 nm
Laser Power: 1W – 25W
Operating Modes: CW or Pulsed
Aiming Beam: 650 nm, < 4.5 mW
Dimensions: 43cm x 28cm x 26cm
Weight: 9kg

ENQUIRE NOW

WARRANTY

Device: 24 month

PART NUMBER	DESCRIPTION
SYS-000837-000-00	LightForce XPi (25W) Therapy Laser System

CHATTANOOGA

LIGHTFORCE® FXi



15W

Treat quickly and consistently on the go. The 15W LightForce® FXi is ideal for those practices or athletic training rooms needing to treat in multiple locations in or out of the office. Easily transportable with an integral battery, the lightweight FXi quickly impacts pain while weighing in at only 3.2kg.

SPECIFICATION

Laser Type : Class IV, Solid State
Laser Wavelength: 980/810 nm
Laser Power: 0.5W – 15W
Operating Modes: CW or Pulsed
Aiming Beam: 650 nm, < 4.5 mW
Dimensions: 438m x 28cm x 20cm
Weight: 3.2kg

ENQUIRE NOW

WARRANTY

Device: 24 month

PART NUMBER	DESCRIPTION
SYS-000829-000-00	LightForce FXi (15W) Therapy Laser System

EXCLUSIVE LIGHTFORCE® TREATMENT HEADS



SMALL CONE

Systems: 25W XPi, 15W FXi
When treating very narrow or precise areas, the Small Cone is ideal due to its smaller spot size.




LARGE CONE

Systems: 40W XLi, 25W XPi, 15W FXi
The Large Cone is designed for off-contact treatments that require higher power settings and is ideal for treating larger areas or when deeper tissues need to be addressed. The Large Cone is not recommended for treatment around the face or over smaller areas due to its larger spot size.



SMALL MASSAGE BALL

Systems: 40W XLi, 25W XPi, 15W FXi
Designed to be used with low to moderate power, the Small Massage Ball is ideal for treating trigger points, around the face, or the digits of the feet and hands.



XL CONE

System: 40W XLi
Coupled with the patented Smart Handpiece, the XL Cone has the unique ability to keep high power treatment comfortable for most patients. With a spot size that is twice the size of our Large Cone, it allows clinicians to safely treat large areas and deeper conditions at over 25 Watts of power.



FLAT WINDOW

Systems: 40W XLi, 25W XPi, 15W FXi
The Flat Window is great for treating superficial tissues with small to moderate surface areas, especially in the presence of boney prominences.



ROLLING PIN

System: 40W XLi
LightForce created the Rolling Pin to allow large areas to be treated with the 40W XLi with a one of a kind, on-contact head that makes treatments simple to administer.



LARGE MASSAGE BALL

Systems: 40W XLi, 25W XPi, 15W FXi
The most commonly used on-contact head, the Massage Ball reduces light scatter on the surface while allowing physical compression of tissue. Paired with higher power settings, it is the ideal tool to address problems in the back, hips, and larger muscle groups.

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
PRD-000232-000-L	Small Cone (XPi & FXi)
PRD-000231-000-L	Small Massage Ball (XLi, XPi & FXi)
PRD-000233-000	Flat Window (XLi, XPi & FXi)
PRD-000229-010-00	Large Massage Ball (XLi, XPi & FXi)
PRD-000230-010-L	Large Cone (XLi, XPi & FXi)
PRD-000644-000	XL Cone (XLi)
PRD-000721-000	Rolling Pin (XLi)

ACCESSORIES

PART NUMBER	DESCRIPTION
PRD-000807-000-01	Carry Bag (FXi)
PRD-000683-000-01	Carry Bag (XLi, XPi)
PRT000473	Safety Goggles
PRD-000424-010-00-L	Empower IQ Delivery System Handpiece (FXi)
PRD-000758-010-00-L	Empower IQ Delivery System Handpiece (XLi, XPi)
PRD-000204-010-L	Hard Rolling Case (FXi)
PRD-000215-010-L	Hard Rolling Case (XLi, XPi)
PRD-000207-010-L	Laser Cart and Bin (XLi, XPi, FXi)



OPERATION

2:29
Time (Min:Sec)

8.0
Power (W)

ON
Treatment

Ready

2:29

Time Remaining (Min:Sec)

Lightforce



CHATTANOOGA

VITALSTIM® PLUS

VitalStim Plus neuromuscular electrical stimulation is an interactive therapy that aids muscle strengthening to rehabilitate swallowing. Pre-set or customised programmes allow the provider to truly customise treatments while biofeedback and visuals create a patient-focused experience that encourages engagement and achievement of goals.

INDICATIONS

To recruit and re-educate muscles used in the process of swallowing.



CHARACTERISTICS

- Four separate channels deliver stimulation concurrently and pause at the same time
- Adjustable treatment protocol, such as phase duration, frequency, work and rest time
- Modern touch screen design with stylus for ease of use.
- 256Mb removable micro SD card for storing patient protocols
- Live treatment progress can be mirrored onto a computer screen, allowing clinicians to guide patients through personalised, targeted swallowing exercises

- sEMG biofeedback helps to increase effort and duration of swallowing attempts and to improve coordination
- Supplied with coloured snap lead wires, EMG lead wire, operator remote hand switch, stand and protective rubber sleeve
- Powered by four AA type alkaline batteries (included)
- 2 year warranty, 6 month warranty on accessories

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
5923-3	VitalStim Plus Electrotherapy System

PART NUMBER	VITALSTIM ADULT ELECTRODES
59000	Pack of 12
59042	Pack of 30
59043	Pack of 50
59044	Pack of 100





CHATTANOOGA TENS/NMES

PHYSIO INDICATIONS

For all indications of Electrotherapy: Standard (TENS/NMES), Full (Neurology, ACL) and Denervated muscles, Iontophoresis, Hyperhidrosis, Urology or oedema prevention.

THETA INDICATIONS

For both standard (Pain, NMES) and specific protocols in Rehabilitation especially when triggered or co-contractions are requested. Neurology (Spasticity, Hemiplegic foot, Agonist-Antagonist) or specific pathologies (eg Post ACL Ligamentoplasty) are protocols provided by the unit.

REHAB INDICATIONS

For pain management, treatment of vascular diseases and muscle stimulation with basic programmes.



CHARACTERISTICS

- Four channel stimulator with Muscle Intelligence™ Technology on all channels (Mi-Scan ready, Mi-Tens ready, Mi-Range ready and Mi-Action ready, Physio Includes 4 x MI Sensors, Theta includes 2 MI Sensors)
- Four programme categories: Rehabilitation I, Pain Relief I, Vascular and Conditioning I
- 2+2 feature allows one programme to run on channels 1 and 2, and another programme on channels 3 and 4


- Trigger function for working phase voluntary when requested in protocols
- Motor Point Pen is included for optimal electrode placement
- Create and personalise programmes in order to customise the patient specific parameters
- Guidelines and treatment recommendations as well as technical parameters for all programmes
- Kit includes rechargeable battery (charging time will take a maximum of 2.5 hours)
- 3 year warranty, 6 months for cables and accessories

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
2535115	Physio
2534815	Theta
2533115	Rehab

CHATTANOOGA

WIRELESS PRO 4 CHANNEL (FULL)



INDICATIONS

A 4 channel stimulator developed to manage the main indications of electrotherapy (TENS/NMES) and full (Neurology, ACL) with the benefits of wireless freedom. The Trigger function will allow management of new types of sessions where voluntary contractions are required, or when the physiotherapist wants to control stimulation manually.

CHARACTERISTICS

- Eight programme categories: Rehabilitation I&II, Pain Relief I&II, Vascular, Neurologic Rehabilitation and Conditioning I&II.
- Wireless technology - removes constraints for you and your patient. Gives your patient total freedom of movement facilitating patient manipulations and allowing you to combine active exercise with electro stimulation for a better result
- Software allows you to create patient history list and to create customised programmes
- Includes four modules, two hand controllers, smart docking station with storage tray to manage several sessions simultaneously. Also includes a carry case for on the go treatment
- Motor Point Pen is included for optimal electrode placement
- 3 year warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
2532981	Wireless Pro 4 Channel



BATTERY PACK

PART NUMBER	SIZE
941213	2.0Ah Battery Pack



CHARGER

PART NUMBER	SIZE
683012	UK 3 Pin 9v Charger - 1.4Ah




LEAD WIRE WITH SNAP CONNECTOR

PART NUMBER	SIZE
601131L	Set of Snap Connector Lead Wires



MI-SENSOR LEAD WIRE

PART NUMBER	SIZE
6601160	Mi-Sensor Lead Wire




LEAD WIRE WITH PIN CONNECTOR

PART NUMBER	SIZE
601132L	Set of Pin Connector Lead Wires



TRANSMISSION GEL

PART NUMBER	SIZE
602048	Transmission Gel (250g)



SNAP ADAPTERS

PART NUMBER	SIZE
6111944	Snap Adapters (Pack of 12)



MOTOR POINT PEN

PART NUMBER	SIZE
980026	Motor Point Pen





CHATTANOOGA
CEFAR TENS

INDICATIONS

Combines TENS with NMES to bring pain relief and muscle rehabilitation from the clinic to the home.

CHARACTERISTICS

- 10 preset treatment programmes, plus 5 open slots
- 1+1 Function: Perform a treatment using two different programmes simultaneously
- Programme lock, compliance, pain score record
- Large easy to read backlit screen
- Mains rechargeable battery. Fully charges in 90 minutes and charge will last 6-10 hours
- 3 year warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
2106	Cefar TENS



CHATTANOOGA
REHAB X2

INDICATIONS

TENS programmes to treat pain and NMES rehab programmes for neurological and orthopaedic rehabilitation where alternating stimulation, strength programmes with active rest, warm up, recovery and massage programmes is required.

CHARACTERISTICS

- Cefar Easy Touch™ automatic intensity adjustment - at the press of a button the intensity is automatically set and adjusted according to the results of the auto test analysis, making treatment easier and ensuring stimulation takes place at an optimal and effective level
- Cefar Flow Tens: a preset programme featuring a comfortable stimulation method. Use both channels at the same time to combine pain relief with a pleasant massaging sensation
- Three personalised programmes can be added
- Pause function, programme lock, key lock, compliance, self check
- Large easy to read backlit screen
- Powered by two AA alkaline batteries (included)
- Includes Hand Control Trigger
- 2 year warranty, 6 months for accessories

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
2942	Cefar Rehab X2

CHATTANOOGA
DURA-STICK

CHARACTERISTICS

- Cost effective electrodes in which comfort and performance are not sacrificed
- White cloth backing aids in the conformity of the electrode to the skin



5cm x 9cm - PACK OF 2

PART	SHAPE	TYPE	PIN
42189	Rectangle	Cloth	2mm Pin



5cm x 5cm - PACK OF 2

PART	SHAPE	TYPE	PIN
42188	Square	Cloth	2mm Pin



5cm x 9cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42191	Rectangle	Cloth	2mm Pin



5cm x 5cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42190	Square	Cloth	2mm Pin

DURA-STICK® PLUS

CHARACTERISTICS

- Precise control and optimisation of electrical current
- Comfortable and effective stimulation is achieved through the electrode
- Available in cloth or foam backing
- Available with lead wires, proprietary MI snap connection or standard clip



5cm - PACK OF 40

PART	SHAPE	TYPE	PIN
42192	Round	Foam	2mm Pin



5cm x 5cm - PACK OF 40

PART	SHAPE	TYPE	PIN
42193	Square	Foam	2mm Pin



5cm x 9cm - PACK OF 40

PART	SHAPE	TYPE	PIN
42194	Rectangle	Foam	2mm Pin



5cm x 10cm - PACK OF 40

PART	SHAPE	TYPE	PIN
42196	Oval	Foam	2mm Pin



5cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42197	Round	Cloth	2mm Pin



5cm x 5cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42198	Square	Cloth	2mm Pin



5cm x 9cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42199	Rectangle	Cloth	2mm Pin



5cm x 10cm - PACK OF 2

PART	SHAPE	TYPE	PIN
42200	Rectangle	Cloth	2mm 2x Pin



5cm x 5cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42201	Square	Cloth	Clip



5cm x 10cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42202	Rectangle	Cloth	Clip



5cm x 10cm - PACK OF 2

PART	SHAPE	TYPE	PIN
42203	Rectangle	Cloth	2x Snap



5cm x 5cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42204	Square	Cloth	Snap



5cm x 10cm - PACK OF 2

PART	SHAPE	TYPE	PIN
42223	Rectangle	Cloth	Snap



3.2cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42219	Round	Cloth	2mm Pin

CHATTANOOGA

DURA-STICK® PREMIUM


CHARACTERISTICS

- Stainless steel electrodes
- Superior conductivity
- Unmatched electrode conformity – pliable material conforms to irregular body contours yet still maintains superior conductivity properties
- Blue gel electrodes - multi-layer adhesive gel that optimally bonds to the electrode layer and mesh (scrim) layer
- Gel is formulated to allow multiple applications to sensitive skin



3.2cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42205	Round	Mesh	2mm Pin




4cm x 6cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42207	Oval	Mesh	2mm Pin



8cm x 13cm - PACK OF 2

PART	SHAPE	TYPE	PIN
42208	Oval	Mesh	2mm Pin



5cm x 5cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42209	Square	Mesh	2mm Pin



5cm x 9cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42210	Rectangle	Mesh	2mm Pin



5cm x 5cm - PACK OF 4

PART	SHAPE	TYPE	PIN
42211	Square	Mesh and Blue Gel	2mm Pin



SAUNDERS

LUMBAR HOME TRACTION DEVICE

INDICATIONS

Lumbar stretch, traction of the lower back in the home, workplace or whilst travelling.

CHARACTERISTICS

- Friction free, actively moving treatment surface ensures efficient application of traction and a split table design replicates smooth clinical traction and promotes consistency of treatment
- Easy to read pressure gauge provides total control of the force at all times and a pressure release valve limits the amount of force that can be applied
- Carry case with detachable wheels offers enhanced portability
- Non-slip harnesses incorporate belt organisers for quick set up
- 1 year warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
199603	Saunders Lumbar Home Traction Device



SAUNDERS

CERVICAL HOME TRACTION DEVICE

INDICATIONS

Cervical traction in the home, workplace or whilst travelling.

CHARACTERISTICS

- Traction forces are directed to the back of the head to prevent compression of the temporomandibular joint
- Easy to read pressure gauge provides total control of the force at all times and a pressure release valve limits the amount of force that can be applied
- Carry case offers enhanced portability
- Self adjusting neck wedges ensure a perfect fit for the patient
- 1 year warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
199594	Saunders Cervical Home Traction Device

ARTROMOT®-K1



INDICATIONS

- Arthrotomy and arthroscopy procedures in combination with synovectomy, arthrolysis
- Treatment following mobilisation of joints in anaesthetised patients
- Surgical treatment of fractures and pseudoarthrosis
- Exercise-stable osteosynthesis
- Operations on soft tissue in the joint area
- Patellectomy
- Corrective osteotomy
- Meniscectomy
- Total knee/hip replacement
- ACL/PCL reconstruction
- Complex knee surgery
- Articular cartilage surgery, e.g. autologous chondrocyte transplantation (ACT)

CHARACTERISTICS

- ARTROMOT-K1 combines the best features of the time-tested ARTROMOT knee machine with state of the art technology
- ARTROMOT-K1 was developed for the demanding use in clinics and home care/ rental programmes
- Using the ARTROMOT-K1 is made easy by intuitive operation of the hand-held programming unit with self-explanatory symbols
- 2 year warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
80.00.045	ARTROMOT-K1 Classic
80.00.040	ARTROMOT-K1 Standard
Please Enquire for Rentals	

PRODUCT WARRANTY

Our products have a limited warranty from the date of purchase (see table below) against defects in materials and workmanship. Normal wear and tear during the usage of a product is not considered a defect.

PRODUCT FAMILY	DONJOY®	PROCARE®	AIRCAST®
CUSTOM DEFIANCE® STEEL HINGE AND FRAME	Lifetime	N/A	N/A
CUSTOM DEFIANCE® STANDARD HINGE AND FRAME	5 years	N/A	N/A
PREFABRICATED RIGID KNEE BRACE FRAME AND HINGE	1 year	6 months	N/A
ASSOCIATED SOFT GOODS*	6 months	6 months	N/A
ELASTIC	6 months	6 months	6 months
WALKER BOOTS	6 months	6 months	6 months
NEOPRENE/DRYTEX	6 months	6 months	6 months
POST-OP BRACING	6 months	6 months	6 months
COLD THERAPY	6 months	N/A	6 months
ASSOCIATED SOFT GOODS*	6 months	N/A	N/A

* Softgoods accessories for rigid and functional bracing requested beyond the 6 month warranty will be sold separately

TERMS & CONDITIONS

1. PARTIES

- a. The "Company" means DJO UK Ltd. (Trading as Enovis)
- b. The "Customer" means the person, firm or company buying goods from the Company to whom the goods are to be invoiced (unless the context otherwise requires) and any person to whom or for the benefit of whom the goods may be delivered.
- c. The "Goods" means the goods described on our invoice.

2. ORDERS

- a. In the UK, credit account customer orders with a value of less than £50.00 (excluding VAT, net discounts) will incur a £7.50 delivery charge. In the UK, credit account customer orders with a value of more than £50.00 (excluding VAT, net discounts) will incur a £4.99 delivery charge. In the Republic of Ireland, credit account customer orders with a value of less than €50.00 (net discounts) will incur a €10.00 delivery charge. Any NHS Supply Chain (NHS SC)/Supply Chain Coordination Limited (SCCL) credit account customer orders must have a Minimum Order Value (MOV) of €115.00 (excluding VAT, net discounts), otherwise they will be rejected.
- b. Any orders submitted or variations made orally must be confirmed in writing within 5 days.
- c. All orders are to be forwarded to the address below.
- d. V.A.T. Exemption – Buyers who qualify for exemption from paying V.A.T. should include a valid V.A.T. exemption certificate (or other documentation approved by the Customs and Excise for the purpose) with each order.
- e. Orders to non credit account customers require payment to accompany the order before the Goods can be dispatched.

3. QUOTATIONS

Any estimate or quotation shall comprise an invitation to treat and will remain open for a period of 30 days provided that the Company has not previously withdrawn it.

4. SPECIFICATION

Every effort is made to ensure our catalogue descriptions and illustrations are clear and accurate. All colours and sizes are approximate. The reproduction of colours is as accurate as photographic and printing processes allow. Due to our policy of continual improvement, some products may vary from those described. The company reserves the right to amend the specifications, construction and/or design without prior notice.

5. PRICES

- a. All prices quoted are exclusive of Value Added Tax.
- b. The Company reserves the right to vary the price of the Goods to take account of any increases in the cost of raw materials, manufacture, packing, transport, wages or otherwise.
- c. Prices charged will be those effective at the date of receipt of Customer order.

6. PAYMENT

- a. Unless the Company otherwise agrees, settlement is due within 30 days of date of invoice to DJO UK Ltd at the address below.
- b. The Company reserves the right to charge interest on all overdue payments and outstanding balances at the rate of 3% per month.
- c. The Customer shall not be entitled to withhold payment either wholly or in part once the same has become due except where expressly agreed with the Company.
- d. Non credit account customers and sales to members of the public require either payment of a pro-forma invoice or payment to accompany the order before the Goods will be dispatched.
- e. Any invoice dispute needs to be raised within 20 days from receipt of invoice

7. DELIVERY / CARRIAGE CHARGES

- a. Unless the Company otherwise agrees, delivery shall be made to the delivery point specified in the order.
- b. Stated delivery times are an estimate on the part of the Company and shall not be binding upon the Company which shall not be liable for any costs or damages caused by reason of any delay in delivery.
- c. The Company reserve the right to make delivery by instalments and to tender a separate invoice in respect of each instalment.
- d. Orders received before 3.00pm Monday to Friday will only be dispatched the same day for UK customers with credit accounts that do not have overdue payments. This is subject to products being in stock and does not apply to custom manufactured products and non stock items. Delivery will be on a next day basis with some exceptions. Deliveries to remote areas of Scotland and the Republic of Ireland may take up to 72 hours.
- e. Standard Delivery Charges
 - (i) All deliveries to UK residential addresses including webshop orders and drop ship addresses will incur a delivery charge of £7.50
 - (ii) Timed and Saturday deliveries can be arranged at a cost. This will be detailed at the time of ordering.
 - (iii) Oversized items are subject to additional delivery and handling charges, regardless of the value of the order.
- f. Oversized products are delivered to a ground floor location and are not installed in-situ. Should additional delivery services be required, they must be specified at the time of order and may be subject to additional costs.
- g. Delivery for custom manufactured products may vary from 1 week to 6 weeks from date of order.
- h. Spare parts, unusual sized products and special order items can take between 10 and 15 working days for delivery.
 - i) DJO UK Shared Services (SS) opening times are Monday to Friday (excluding bank holidays) from 9:00am to 4:00pm.

8. CANCELLATIONS

1. No notice of cancellation of an order shall be valid unless given in writing and agreed to by an authorised representative of the Company.
2. If written notice of cancellation is received by the Company, the Company reserves the right to levy a cancellation charge up to 50% of the amount of the order.
3. If the contract provides for delivery by instalments, each instalment shall be deemed to be the subject of a separate contract, and non-delivery or delay in the delivery of any instalment shall not entitle the Customer to repudiate the contract as regards the balance.

9. PASSING OF TITLE

The legal title to the Goods shall not pass to the Customer until all sums due or payable by the Customer to the Company whether in respect of the Goods or otherwise have been received by the Company.

10. PASSING OF RISK

The risk in the Goods shall pass to the Customer when the Goods are tendered for delivery at the delivery point specified in the order. In the event that the Goods are collected by the Customer from the Company, the risk of the Goods shall pass to the Customer when they are loaded onto the vehicle of the Customer or of his carrier or other agent.

11. WARRANTY

The Company warrants that the Goods shall at the time of delivery correspond to the specifications published by the Company when used for purposes for which goods of that type are ordinarily used. In the absence of any such published specifications applicable to the Goods, the Company warrants that the Goods shall at the time of delivery be of merchantable quality. If the Goods do not conform to that warranty, the Company will, at its option:

- a. Replace the Goods found not to conform to the warranty;
- b. Bring the Goods into conformity with the published specifications of the Company or (as the case may be) merchantable quality; or
- c. Take back the Goods found not to conform to the warranty and refund the appropriate part of the purchase price.

12. FORCE MAJEURE

The Company shall not be liable for any loss or damage caused by circumstances beyond the Company's control, including but not limited to, an act of God, war, civil disturbance, governmental restrictions, import or export regulations, industrial disputes and difficulties in obtaining labour or materials. Should any such event occur, the Company might cancel or suspend this contract without incurring any liability for any loss or damage caused.

13. NOTICE OF DAMAGE, NON-DELIVERY OR WRONG DELIVERY

- a. The Company gives no guarantee of the exact quantity of goods delivered, and in the event that quantity of Goods delivered does not correspond with the quantity stated on the order the customer shall only be liable to pay for the quantity delivered.
- b. If upon receipt the Goods appear damaged or incomplete, please indicate on the carrier's documentation and notify our Customer Services Department of the situation no later than 3 days after delivery. If in doubt, please sign the carrier's documentation "Unchecked". Should an expected consignment not arrive, please notify our Customer Services Department within 21 days from date of invoice, failing which the company shall be under no liability for non-delivery.
- c. Any surplus goods delivered shall remain the property of the Company, and the Customer shall take all reasonable precautions for the safe custody and protection of such surplus goods until the time of their removal by the Company.
- d. In the event of goods being wrongly delivered, under no circumstances shall any goods be returned to the Company without its prior consent, otherwise liability is not accepted. The Company will issue a Returns Authorisation Number (RAN) and this must be marked on the exterior of every package returned. Returned Goods must be suitably packaged for transportation by the carrier. Goods received without suitable packaging will not be credited. Where the Company has made the error, full credit will be given on uplifted goods in good condition.
- e. The Company is under no obligation to accept goods for return, e.g. in the case of Customer ordering error, surplus stock, etc., but where it agrees to do so the following handling charges will apply:-
 - 30 days of delivery, no restocking fee
 - 31 Days – 365 Days of delivery, a 15% restocking/packaging fee
 - Over 12 months of delivery, no credit will be given.The Customer will arrange carriage and pay all freight charges. The Customer should obtain a Returns Authorisation Number, which should be marked clearly on each package being returned. Only complete units of sale may be returned and all goods are to be in good condition. The Customer will be liable for a 15% restocking fee for any goods returned in a non saleable condition or with missing accessories, regardless of the number of days since delivery, with the balance of the goods being paid for in full by the Customer.
- f. Refused Deliveries – We reserve the right to make a carriage charge plus 15% handling charge on all consignments which are tendered for delivery and refused without justification. Goods returned in a non-saleable condition will be paid for in full by the Customer. Goods should not be returned without the prior agreement of the Company.
- g. Only unused and complete goods may be returned, and they must be undamaged, complete with all accessories and less than 12 months old. Goods not conforming to these criteria will not be accepted for credit. The Company reserves the right to dispose of any non saleable and non creditable returned goods.
- h. Custom manufactured products and non stock products are non returnable.

14. DISPUTES AND SET-OFF

Any liability of the Company under the contract shall be subject to and conditional upon the due performance and observance by the Customer of all its obligations under these conditions. Subject to these conditions, the Customer shall not be entitled to withhold or delay payment or exercise any right or set-off whatsoever and howsoever arising or arisen which might otherwise be available to it.

15. COPYRIGHT

All copyrights shall remain the property of the Company alone unless otherwise expressly agreed in writing with the Customer. All product images remain the property of the Company and may not be reproduced in print or electronic format without the express written permission of the Company.

16. RETAIL/RESELLING

Reselling of DJO branded products is forbidden without written consent from the Company.

17. PATENT RIGHTS

The Customer warrants that any instructions furnished or given by the Customer shall not be such as will cause the Company to infringe any letters patent or trademark or other intellectual property right in execution of the Customer's order and shall indemnify the Company against all actions, costs, claims and such demands arising therefrom to the extent that such instructions lead to such infringements.

18. GOVERNMENT CHARGES

The Customer will reimburse the Company for any increase in tax or governmental charge, or for any new tax or governmental charge hereafter becoming effective which the Company may be required to pay to any government upon the sale, production, or transportation of the Goods and which has the effect of increasing the cost of the Goods sold hereunder.

19. ASSIGNMENT

This contract is personal to the Customer and may only be assigned by it with the prior written consent of the Company.

20. BANKRUPTCY / LIQUIDATION

If the Customer being an individual (or, when the Customer is a firm, any partner in that firm) shall at any time become bankrupt, or having a receiving order made against them or shall make any composition or arrangement with or for the benefit of the creditors or if the Customer, being a company, an order for winding up is made or a resolution for winding up is passed or if a receiver or manager is appointed on behalf of a creditor or otherwise, the Company shall have the right to stop Goods in transit, to suspend further deliveries and to determine any subsisting contract with the Customer and (in any case where title to the Goods has to be passed to the Customer but payment is due and has not been made) the Company shall have a lien over the Goods and right of resale, provided that any such action shall be without prejudice to any right or claim the Company may have against the Customer.

21. CUSTOMER CARE

In our quest to "Never Stop Getting Better", we may from time to time monitor and record customer communications. This is for internal training purposes only and the confidentiality of our Customers and their transactions will be fully protected at all times.

22. PERSONAL CALLERS

We regret that we are unable to deal directly with personal callers to our office.

23. TERMS AND CONDITIONS

Terms and conditions quoted in this catalogue supersede all previous terms and conditions unless otherwise stated.

24. COMPANY DETAILS

DJO UK Ltd
1a Guildford Business Park Guildford, Surrey GU2 8XG Tel:
+44 (0) 1483 459659 Fax: +44 (0) 1483 459470
Company Registration Number: 4276829 VAT Registration Number: 785 0868 83
Eire VAT Registration Number: IE 998 7359 J



enovis™

1a Guildford Business Park, Guildford, Surrey, GU2 8XG, UNITED KINGDOM

[enovis.com/uk](https://www.enovis.com/uk) | 01483 459659

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

MKT-000-00011973-EN - Rev D

Copyright © 2024 DJO, LLC, a subsidiary of Enovis Corporation