



Playmaker® II

Optimal anvendelse

Selv lette ligamentskader kan sætte en stopper for helt almindelige aktiviteter. Playmaker® II fra DonJoy giver den stabilitet og komfort som kræves for at styrke bevægeevnen.

Indikationer:

Egnet til mild instabilitet i: • Forreste korsbånd (ACL) • Bagerste korsbånd (PCL) • Mediale kollaterale ligament (MCL) • Laterale kollaterale ligament (LCL) • Menisk

Fordele og egenskaber

- Hi10™ (bagside i slidstærkt nylon/elastanstof) giver enestående komfort og pasform
- DonJoys effektive DIII-hængselled tilbyder et design der er tilpasset til at håndtere ROM-indstillinger
- Fremstillet i et materiale (spacer fabric) som er åndbart og ikke bliver for varmt
- Formstøbte spænder muliggør en hurtig og nem justering af remmene til en sikker og behagelig pasform
- Sleeve eller wraparound model
- Forcepointled eller standardled

Bestillingsinformationer



Playmaker® II	
Varenummer	Beskrivelse
11-3495-x	Playmaker II, Standard Hinge, Sleeve
11-3497-x	Playmaker II, Standard Hinge, Wrap
11-3499-x	Playmaker II, Standard Hinge, Sleeve with Donut
11-3500-x	Playmaker II, Neoprene Sleeve
11-3502-x	Playmaker II, Neoprene Wrap
11-3503-x	Playmaker II, Neoprene Sleeve with patella donut
11-3504-x	Playmaker II, Forcepoint, Sleeve
11-1482-x	Playmaker I, DropLock, Wraparound w/ Popliteal Cut-Out, Drytex®