


enovis™

DONJOY®
OA ADJUSTER™ 3

Mighty. Tough.



DONJOY®

OA ADJUSTER™ 3

The improved OA Adjuster™ 3 is reinforced with a tough aluminum frame, making it an extremely durable and functional DonJoy OA off-loader knee brace.

Mighty. Tough.



CLINICALLY-PROVEN*

DonJoy's Adjuster™ technology offloads a little—or a lot

FLOATING HINGES

allow a contoured, personalised fit for all levels of activity, with double upright design

FOUR POINTS OF LEVERAGE™ SYSTEM PIONEERED BY DONJOY PROVIDES CLINICALLY PROVEN LIGAMENT STABILITY**

- ① Anterior thigh anchor
- ② Posterior tibia anchor
- ③ Lower thigh strap pushes femur anteriorly
- ④ Anterior tibia strap applies an active constant load to prevent anterior tibial translation



CONTOURED DESIGN

conforms specifically to the calf, allowing minimal migration and maximum comfort for even the toughest patients



PART NUMBER				THIGH CIRCUMFERENCE 15 cm ABOVE MID-PATELLA	SIZE
LEFT LEG		RIGHT LEG			
MEDIAL	LATERAL	MEDIAL	LATERAL		
11-1591-1	11-1593-1	11-1590-1	11-1592-1	33 – 39 cm	XS
11-1591-2	11-1593-2	11-1590-2	11-1592-2	39 – 47 cm	S
11-1591-3	11-1593-3	11-1590-3	11-1592-3	47 – 53 cm	M
11-1591-4	11-1593-4	11-1590-4	11-1592-4	53 – 60 cm	L
11-1591-5	11-1593-5	11-1590-5	11-1592-5	60 – 67 cm	XL
11-1591-6	11-1593-6	11-1590-6	11-1592-6	67 – 75 cm	XXL
11-1591-7	11-1593-7	11-1590-7	11-1592-7	75 – 81 cm	XXXL



PART NUMBER	DESCRIPTION	SIZE
25-1376-0-00000	Flexion/Extension Stop Kit	Universal
11-0075-x-00000	Cotton Lycra Suspension Sleeve	S - XXL

*Matthew Nadaud, Richard Komistek et al, 2006
 (J Bone Joint Surg AM, 87:114-119, 2005)
 **Fleming, Beynnon, et al 2000
 (Am J Sports Med. 2000;28:815-824.)

DONJOY®

Asia Pacific

1905, Tower 2, Grand Central Plaza
Shatin
HONG KONG
T: 852 31051415
E: marketing.ap@enovis.com

Australia – New Zealand

Suite 1, 25 Frenchs Forest Rd E
Frenchs Forest NSW 2086
AUSTRALIA
T: 1300 66 77 30
E: customerservice.au@enovis.com

Benelux

Kleinhoefstraat 5, bus 39
Geel 2440
BELGIUM
T: +32 (0)14 248090
E: orders.bx@enovis.com

Canada

6485 Kennedy Rd
Mississauga
Ontario, L5T2W4
CANADA
T: 1-866-866-5031
E: canada.orders@enovis.com

China

A312 SOHO ZhongShan Plaza
1055W, West Zhongshan Rd
Shanghai, 200051
CHINA
T: (8621) 6031 9989
E: po.meng@enovis.com

France

Centre Européen de Fret
3 rue de Bethar
64990 Mouguerre
FRANCE
T: +33 (0)5 59 52 86 90
E: sce.cial@enovis.com

Germany

Bötzingen Straße 90
79111 Freiburg
GERMANY
T: +49 (0) 180 1 676 333
E: kundenservice@enovis.com

India

Plot No. 5277, 8th Street, Ram Nagar
North Extension (Ground Floor)
Madipakkam, Chennai – 600 091
INDIA
T: (+91) 44 4978 7760
E: sunil.mathew@enovis.com

Italy

Centro Direzionale Milanofiori
Strada 1, Palazzo F9
20090 Assago (MI)
ITALY
T: +39 02 835 98001
E: vendite@enovis.com

Nordic

Murmansgatan 126
212 25 Malmö
SWEDEN
T: +46 40 39 40 00
E: info.nordic@enovis.com

South Africa

Unit 6, 2 on London, 2 London Circle
BrackenGate Business Park
Brackenfell, 7560
SOUTH AFRICA
T: +27 (0) 21 276 2968
E: info.southafrica@enovis.com

Spain

Ctra. D'Esplugues, 225 2ª B
08940 Cornellà de Llobregat
Barcelona
SPAIN
T: +34 93 480 32 02
E: ventas@enovis.com

UK & Ireland

2000 Cathedral Square Cathedral Hill
Guildford, Surrey, GU2 7YL
UK
T: +44 (0)1483 459 659
E: ukorders@enovis.com

USA

2900 Lake Vista Drive
Lewisville
TX 75067
U.S.A.
T: 760 727 1280
E: customercare@enovis.com

Export

Centre Européen de Fret
3 rue de Bethar
64990 Mouguerre
FRANCE
E: orders.exp@enovis.com

enovis[™]

Creating Better Together[™]

enovis-medtech.eu

2000 Cathedral Square, Cathedral Hill, Guildford, Surrey, GU2 7YL, UNITED KINGDOM

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.