

enovis™

**DONJOY®**  
KIT ICE WRIST

The power of cold



# DONJOY®

## KIT ICE WRIST

Discover the power of combined cryotherapy and compression in a single product.

Consisting of a wrist wrap with an integrated gel cold pack and inflatable aircell, the Kit Ice Wrist is a convenient way to provide reusable post-operative or post-trauma care for the wrist.

### INTENDED USE/INDICATIONS

Post-operative care of wrist surgery, non-surgical trauma, wrist sprain, fractures of the hand and wrist, and rehabilitation. The combined cold therapy and compression may help provide analgesic and anti-haemorrhagic effects and reduce edema.

### FEATURES AND BENEFITS

- **Combined cryotherapy and compression**
  - **Integrated aircell** inflated using manual hand pump for easily controlled compression, as a built-in compression feature may help enhance therapeutic effects of cryotherapy
  - Even when frozen, **integrated gel cold pack** remains moldable to adapt to the wrist
- **Hook-and-loop closure** for easy, secure fitting
- **Bilateral** – fits either hand
- **Universal sizing** – one size fits all



- **Included in kit:**
  - 1 wrist cold pack wrap
  - 1 skin protective fabric
  - 1 hand pump
  - 1 freezer bag

### ORDERING INFORMATION

PART NUMBER	DESCRIPTION
82-0171	Kit Ice Wrist

**enovis**<sup>™</sup>  
Creating Better Together<sup>™</sup>

[enovis-medtech.eu](http://enovis-medtech.eu)

2000 Cathedral Square, Cathedral Hill, Guildford, Surrey, GU2 7YL, UNITED KINGDOM

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice.

The contents of this document do not constitute medical, legal, or any other type of professional advice.

Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

MKT-10702885262-EN-12 - Rev A

Copyright © 2026 DJO, LLC, a subsidiary of Enovis Corporation